



Aug 17, 2020

Cell harvest and RNA prep using Qiagen RNAeasy Kit

Katrin Blondrath¹¹RIKEN**1** Works for me This protocol is published without a DOI.

Katrin Blondrath

PROTOCOL CITATION

Katrin Blondrath 2020. Cell harvest and RNA prep using Qiagen RNAeasy Kit. **protocols.io**
<https://protocols.io/view/cell-harvest-and-rna-prep-using-qiagen-rnaeasy-kit-bf4sjqwe>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

May 07, 2020

LAST MODIFIED

Aug 17, 2020

PROTOCOL INTEGER ID

36722

MATERIALS

NAME	CATALOG #	VENDOR
RNAeasy mini kit	74106	Qiagen
Parafilm		
beta mercaptoethanol		
low bind nuclease free eppendorf tubes		
RNAse zap		
Ice		
Ethanol		
Nano drop or similar		
384well plate		
seal for 384well plate		
One-step PrimeScript RT-PCR kit for real-time RT-PCR		
One Step TB Green™ PrimeScript™ RT-PCR Kit II	RR086A	Takara

BEFORE STARTING

have ice ready
wipe surfaces and pipettes with RNAse zap
prepare 70% ethanol

Cell harvesting and RNA extraction

- 1 Aspirate medium from 6 well plates, wash once with PBS
- 2 Add **500 µl** lysis buffer containig freshly added beta-mercaptoethanol to each well (lysis buffer provided in RNAeasy kit)



Beta-mercaptoethanol should be added freshly every time performing the experiment

- 3 (OPTIONAL) Wrap plates in parafilm, store at -80°C until extraction
- 4 (Defrost plates and) Scrap with 1000 ul pipette tip and pipett up and down the well, transfer into low bind nuclease-free tube
- 5 Extract according to manufacturer's protocol using Qiagen RNAeasy kit (including DNase treatment step) 1h

Storage

- 6 Measure RNA concentration using Nanodrop and store at -80°C or continue with qPCR

qPCR

- 7 Dilute samples to have equal concentrations of 2500 ng of each sample.
Use $2\text{ }\mu\text{l}$ of this diluted RNA directly in 10 ul reaction to achieve 50nM per well of 384 well reaction plate
- 8 Mix contents of One Step TB Green™ PrimeScript™ RT-PCR Kit II (in case running a whole 384well plate)

Pre Master: Mix contents of entire kit:

- $2500\text{ }\mu\text{l}$ TB Green buffer
- $200\text{ }\mu\text{l}$ RT/Polymerase enzyme
- $100\text{ }\mu\text{l}$ ROX
- $800\text{ }\mu\text{l}$ H₂O
- $3600\text{ }\mu\text{l}$ total

Adjust accordingly to the total number of required reactions

1x Mastermix:

	1x in micro-liters
ENZYME	0.4
TB green Buffer 4	5
Primer 30uM F/R mix	0.8
RNA	2
ROX	0.4
H2O	1.8

9 Final master mix:

Per reaction/ per well of 384well plate = 10 µl

Multiply by amount of samples and add ~20% dead volume.