



Mar 13, 2022

Rapid Run v2 Primer Rehyb

Amanda Chilaka¹, Sumeet Gupta¹¹Whitehead Institute Genome Technology Core

1



dx.doi.org/10.17504/protocols.io.b46nqzde



Amanda Chilaka

Switch the flowcell to another sequencer if first base fails - HiSeq

DOI

dx.doi.org/10.17504/protocols.io.b46nqzde

Amanda Chilaka, Sumeet Gupta 2022. Rapid Run v2 Primer Rehyb. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b46nqzde>



_____ protocol ,

Feb 16, 2022

Mar 13, 2022

58286

If a rapid run has failed at first base and you want to switch the flowcell to another side or sequencer, follow these steps (**ONLY applicable to single read runs, changing sequencer for paired end run is not possible**)

These steps are to be performed on the HiSeq where the flow cell is being moved.

ONLY applicable to single read runs, changing sequencer for paired end run is not possible

Rapid Run v2 Primer Rehyb 2h 16m 30s

1

1m

Wash the sequencer that will be used (MWB and then water). [Link to wash protocol]

Load the SBS rack (with reagents) to the sequencer.

1m

2

3 Load **water** in the Paired End rack and the template tubes. 1m

4 Place a dummy flowcell on the deck. Engage vacuum. 30s

5 Setup the run as normally you would, but select “cluster generation on cbot.” even if it was 1m
onboard clustering. This will save an additional hour of time.

6 Enter the other parameters used for the flowcell, such as read lengths, flowcell ID, etc...This 1m
will create the files needed to allow you to do a rehyb later.

7 Start the run and wait for the first base report to come up. End the run. 1h 30m

8 Thaw a Rehybridization kit and replace or top-off the reagents according to the “HiSeq Rapid 30m
Primer Rehybridization Reference Guide.” A general summary is below.

8.1 For **Single-Read Runs**, replace/top-off the following reagents with these 5m
from the Rehyb kit:

Position	Reagent	Description
5	USB	Add to current USB bottle
18	HP10	Replaces Read1 primer
17	HP12	Replaces Index1 primer (i7)
16	HP9	Replaces Index2 primer (i5)

9 Leave the dummy flowcell on the sequencer. You will prime following reagent twice: 5m

1. FDR (position 15) - aspiration rate of 1500, dispense rate of 2000, and using 250uL
2. Flow USB (5) - aspiration rate of 1500, dispense rate of 2000, and using 250uL

10 Return to the Main screen and select "Sequence," then "Rehyb Run."

1m

Replace the dummy flowcell with sequencing flowcell and proceed with rehybridization.