



Oct 24, 2020

Agarose Gel Electrophoresis

Zhujun Wei¹¹Northeast Forestry University

1

Works for me

dx.doi.org/10.17504/protocols.io.bnwsmfee

2020 iGEM NEFU China



Zhujun Wei

ABSTRACT

The concentration of the agarose gel is 1%.

DOI

dx.doi.org/10.17504/protocols.io.bnwsmfee

PROTOCOL CITATION

Zhujun Wei 2020. Agarose Gel Electrophoresis. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bnwsmfee>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Oct 24, 2020

LAST MODIFIED

Oct 24, 2020

PROTOCOL INTEGER ID

43698

MATERIALS TEXT

Agarose
50X TAE
Gelred

SAFETY WARNINGS

Please wear gloves for the experiment, Be careful when heating, protect your hands.

ABSTRACT

The concentration of the agarose gel is 1%.

- 1 Weigh 0.20 g of agarose.
- 2 Add 20 ml of TAE 1X.
Need to prepare 1X TAE with 50X TAE.

- 3 Heat up until the solution is homogeneous.
- 4 Add 1 μ l Gelred into solution and mix well.
- 5 Put the solution into bed for polymerize, make sure "comb" is well placed and the solution is balanced.
- 6 20 min later, use the gel.