



JUL 28, 2023

🌐 Cellular lipid uptake with flow cytometry readout

Igor Beletchi¹, rosanne.wouters¹, Peter Vangheluwe¹

¹KU Leuven



rosanne.wouters

ABSTRACT

Cellular lipid uptake with flow cytometry readout

OPEN  ACCESS



DOI:

dx.doi.org/10.17504/protocols.io.14egn2pnpg5d/v1

Protocol Citation: Igor Beletchi, rosanne.wouters, Peter Vangheluwe 2023. Cellular lipid uptake with flow cytometry readout.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.14egn2pnpg5d/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: May 31, 2023

harvest cells

10m

1 harvest cells by detachment with



TrypLE Express or preferred cell dissociation reagent Thermo Fisher Scientific

1.1 cells should be grown to 70% confluency, not higher



1.2 wash cells with



TrypLE Express or preferred cell dissociation reagent Thermo Fisher Scientific

1.3 incubate cells with



TrypLE Express or preferred cell dissociation reagent Thermo Fisher Scientific

(3ml/T175 flask) at room temperature

1.4 when cells detach resuspend with 7ml



PBS without Ca²⁺ or Mg²⁺ Gibco, ThermoFisher Catalog #10010-031 and keep on ice in 15 mL falcon

2




collect cells by centrifugation



400 x g, 4°C,
00:05:00

5m

3

wash pellet with with  PBS without Ca²⁺ or Mg²⁺ Gibco, ThermoFisher Catalog #10010-031

4

collect cells by centrifugation

400 x g, 4°C,
00:05:00

5m

5

resuspend pellet in



HBSS Gibco - Thermo Fischer Catalog #14170-112

6

count cells

7

prepare cell suspension with 1×10^6 cells in 500µl

HBSS Gibco - Thermo Fischer Catalog #14170-112

prepare lipid suspension

8

work under fume hood equipped with N2 flow

9

wash glass syringe 3x with 100% ethanol and 3x with chloroform, let dry

10


with glass syringe pipette needed volume of NBD-lipid into glass vial





NBD-PC (18:1-06:0 NBD PC) Avanti Polar Lipids, Inc. Catalog #810132C




NBD-lyso-PC (12:0 lyso NBD PC) Avanti Polar Lipids, Inc. Catalog #810128C

 NBD-PS (18:1-06:0 NBD PS) Avanti Polar Lipids, Inc. Catalog #810194C

 NBD-PE (18:1-06:0 NBD PE) Avanti Polar Lipids, Inc. Catalog #8105155C


 NBD-SM (C6 NBD sphingomyelin) Avanti Polar Lipids, Inc. Catalog #810218C

 NBD-GluCer (C6-NBD Glucosyl Ceramide) Avanti Polar Lipids, Inc. Catalog #810222C


10.1 keep glass bottle with NBD-labeled lipids in chloroform solution cold, by working on ice or in freezer block


10.2 for each sample 500µl lipid suspension with 2µM final concentration is needed









11 dry lipids by evaporating chloroform under N₂ flow


12 resuspend lipids in  HBSS Gibco - Thermo Fischer Catalog #14170-112 to final concentration of 2µM with magnetic stirbar until fully dissolved, keep lipids at 4°C while dissolving

lipid uptake

13 equilibrate cell suspension  0 rpm, 37°C, 00:15:00

14 equilibrate lipid suspension  0 rpm, 37°C, 00:15:00

- 15 prepare collection tubes with 200µl  HBSS Gibco - Thermo Fischer Catalog #14170-112
containing 5%  Albumin Bovine Serum Fraction V Fatty Acid-Free Merck Catalog #126575
(per sample 1 collection tube is needed for each timepoint) and keep on ice
- 16 for timepoint 0, take 200µl of equilibrated cell suspension and add to collection tube T=0
(containing 200µl of ice-cold  HBSS Gibco - Thermo Fischer Catalog #14170-112
containing 5%  Albumin Bovine Serum Fraction V Fatty Acid-Free Merck Catalog #126575)
on ice
- 17 add 500µl of lipid suspension to remaining 300µl of cell suspension
- 18 incubate at 37°C with thermoshaker with interval 1 min shaking every 5 min (700 rpm).
- 18.1 keep samples away from light during incubation
- 19 after 30min, take 200µl of lipid+cell suspension and add to collection tubes T=30 (containing 200µl of ice-cold  HBSS Gibco - Thermo Fischer Catalog #14170-112 containing 5%  Albumin Bovine Serum Fraction V Fatty Acid-Free Merck Catalog #126575), keep on ice
- 20 after 60min, take 200µl of lipid+cell suspension and add to collection tubes T=60 (containing 200µl of ice-cold  HBSS Gibco - Thermo Fischer Catalog #14170-112 containing 5%  Albumin Bovine Serum Fraction V Fatty Acid-Free Merck Catalog #126575), keep on ice

- 
- 21** add sodium dithionite to all samples (1/100 dilution of freshly prepared 1M stock in tris-buffer, pH10) and vortex samples (make sure to mix well with quencher)
 - 22** measure total internal fluorescence with flow cytometer