



1

Mar 01, 2022

Open Field V.1

Haley Geertsma¹

¹University of Ottawa

1



dx.doi.org/10.17504/protocols.io.b5qzq5x6

Haley Geertsma
University of Ottawa

This protocol is used to test mice in the Open Field behavioural test.

DOI

dx.doi.org/10.17504/protocols.io.b5qzq5x6

Haley Geertsma 2022. Open Field. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5qzq5x6>



protocol ,

Mar 01, 2022

Mar 01, 2022

58873

45cm³ open field box

- 1

Habituate mice in a separate room for 60 minutes prior to testing.

1h
- 2

Confirm testing room lighting is set to 100 lux and EthoVision is set up properly.

5m
- 3

Place mice into open field box and record their movement for 10 minutes.

10m

- 4 When testing is complete, place mice back into their home cage and clean the open field boxes well before testing the next mice. 5m