

JUN 13, 2023

OPEN ACCESS

dx.doi.org/10.17504/protocol s.io.3byl4jyr8lo5/v1

Protocol Citation: michela. deleidi, Bianca Marchetti, Federico Bertoli, Carmela Giachino 2023. Gene Expression. protocols.io https://dx.doi.org/10.17504/p rotocols.io.3byl4jyr8lo5/v1

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: Mar 31, 2023

Last Modified: Jun 13, 2023

PROTOCOL integer ID:

79817

Keywords: RT-PCR, quantify

gene expression

Gene Expression

Federico

michela.deleidi^{1,2}, Bianca Marchetti^{3,4}, Bertoli^{1,2}, Carmela Giachino³

¹Mitochondria and Inflammation in Neurodegenerative Diseases, DZNE, Tübingen-Germany;

²Hertie Institute for Clinical Brain Research, University of Tübingen; ³Neuropharmacology Laboratory, Oasi Research Institute-IRCCS, Troina,

⁴Biomedical and Biotechnological Sciences, Pharmacology Section, University of Catania-Italy



carmela.giachino

ABSTRACT

RT-PCR is a technique to quantify gene expression in samples.

ATTACHMENTS

k4drb48if.docx

Materials

QIAzol Lysis Reagent (200ml) **Qiagen Catalog** #79306

RNeasy Lipid Tissue Mini Kit (50) Qiagen Catalog

- Nanodrop-ND 1000
- Retroscript Kit (Ambion #1710, Austin, Texas)
- 2X TaqMan Universal PCR Mastrer Mix (Applied Biosystems)
- Step One Detection System (Applied Biosystems)
- lysis buffer (Qiagen)
- DNase I, RNase-free (Qiagen)
- NanoDrop ND-100 (Nano Drop Technologies)
- the Retroscript Kit (Ambion)
- QIAquick PCR Purification kit (Qiagen)

Gene Expression

1 Perform gene expression on tissue samples of the ventral midbrain (containing the substantia nigra pars compacta, SNpc).

RNA Extraction

- 2 At due time points, homogenize tissue samples in 1 mL of QIAzol Lysis Reagent (Qiagen, #79306) using a rotor-stator homogenizer.
- 3 Isolate the total RNA from homogenized tissue samples using RNeasy Lipid Tissue Kit (Qiagen, #74804) including Dnase digestion.

Reverse Transcription

- 5 Perform cDNA synthesis using Retroscript Kit (Ambion #1710, Austin, Texas) according the manufacturer's instructions.
- 6 Add Δ 50 μL of water solution containing Δ 0.5 μg of each pool to an equal volume of 2X TaqMan Universal PCR Mastrer Mix (Applied Biosystems).

Real Time PCR

- Perform real-time quantitative PCR using TaqmanTM Assay Reagents on an Step One Detection System (Applied Biosystems) according to manufactures protocol.
- **8** Process tissue samples as above.
- Lyse the cell samples in lysis buffer (Qiagen) and store at -80 °C until the RNA extraction following manufacture's instructions.
- Remove residual genomic DNA by incubating with DNase I, RNase-free (Qiagen) and elute from the RNeasy mini columns with RNase-free water.
- Quantify the amount of total RNA was using a NanoDrop ND-100 (Nano Drop Technologies) and synthesize the cDNA from \mathbb{Z} 2 μ g of total RNA using the Retroscript Kit (Ambion).

- After purification using QIAquick PCR Purification kit (Qiagen), use 250 ng of cDNA for Real-time PCR using pre-developed Taqman Assay Reagents (Applied Biosystems).
- Perform real-time quantitative PCR with Step One Detection System (Applied Biosystems) according to manufacturer protocol.



Note

We used the housekeeping gene, β -actin, as normalizer and embryonic mouse brain as calibrator.

- Using the delta delta C_t ($2^{-\Delta\Delta Ct}$) comparative method, quantification of the abundance of target gene expression was determined relative to β -actin with respect to the control group, express the results as arbitrary units (AU).
- 15 Indicate relative fold changes over WT in the respective treatment groups.