

OCT 03, 2023

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.e6nvwdz6zlmk/v1

Protocol Citation: Lucian Dipeso, Emily Hatch 2023. m6A visualization/immunofluorescence of DamID. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.e6nvwdz6zlmk/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


Protocol status: Working
 We use this protocol and it's working

Created: Oct 03, 2023

🌐 m6A visualization/immunofluorescence of DamID

Lucian Dipeso^{1,2}, Emily Hatch^{1,2}

¹University of Washington; ²Fred Hutchinson Cancer Center



Lucian Dipeso

ABSTRACT

Used for visualizing Dam activity by immunofluorescence in mammalian cells. This protocol stains N-6 methyladenosine (m6A) modified DNA by immunofluorescence while preserving other epitopes, avoiding harsh denaturation steps like heat or chemical treatments.

MATERIALS

100% MeOH
 PBS
 BSA
 TritonX-100
 Sodium azide
 DpnI (NEB #r0176)
 CutSmart buffer (NEB #b7204)
 Rabbit anti-m6A (Synaptic #202 003)
 Secondary antibody

Last Modified: Oct 03, 2023

PROTOCOL integer ID:
88738

Keywords: DamID, m6A,
Immunofluorescence

Protocol

- 1 Fix cells in 100% MeOH for 10 min at -20 °C
- 2 Block, and RNase treat in 2 ug/mL RNase A in blocking buffer (3% BSA, 0.04% TritonX-100, 0.02% sodium azide in PBS) at 37 °C for 30 min
- 3 Wash cells 2x in PBS for 5 min at room temperature
- 4 Permeabilize cells in 2% TritonX-100 in PBS for 10 min at room temperature
- 5 Digest in 50 U/mL DpnI (NEB #r0176) in 1x CutSmart buffer for 30 min at 37 °C
- 6 Wash 3x in PBS for 5 min at room temperature
- 7 Incubate cells in rabbit-anti-m6A (Synaptic, #202 003) diluted 1:500 in PBS for 30 min at room temperature in a humidified chamber

- 8 Wash 3x in PBS for 5 min at room temperature
- 9 Incubate cells in anti-rabbit, AlexaFluor conjugated antibody in blocking buffer, diluted according to manufacturer's suggestion, for 30 min at room temperature in a humidified chamber in the dark
- 10 Proceed with further immunofluorescence steps and/or DNA staining and mount cells