



Nov 18, 2020

Curva de Padrão Fosfato 16/11/20

Mariana Jacques¹¹Universidade Federal do Rio de Janeiro*In Development*

This protocol is published without a DOI.

Phosphate**Mariana Jacques**

PROTOCOL CITATION

Mariana Jacques 2020. Curva de Padrão Fosfato 16/11/20. **protocols.io**
<https://protocols.io/view/curva-de-padr-o-fosfato-16-11-20-bpummnu6>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Nov 18, 2020

LAST MODIFIED

Nov 18, 2020

PROTOCOL INTEGER ID

44653

Procedimento:

- 1 Em uma placa de ___ poços colocar :
50 µL de H₂O MiliQ
- 2 A1: 50 µL Pi 40 µM
(Realizar diluições seriadas)
- 3 Adicionar 25 µL MR em cada poço da placa
- 4 Preparar Revelador
4:1 Verde de Malaqueta: Molibdato de Amônio
esperar 30' e filtrar
- 5 Adicionar 200 µL de Revelador por poço
- 6 Esperar 5-10 min
ler no DO 660 nm

