



Jul 07, 2022

TX-100 FRACTIONATION PROTOCOL

Svermily ¹¹University of Minnesota1 *Works for me* Sharedx.doi.org/10.17504/protocols.io.5qpvob21zl4o/v1

Team Lee



jbalster

ABSTRACT

To isolate insoluble and soluble proteins from dissected brain regions frozen for biochemical analysis.

ATTACHMENTS

[452-952.docx](#)

DOI

dx.doi.org/10.17504/protocols.io.5qpvob21zl4o/v1

PROTOCOL CITATION

Svermily 2022. TX-100 FRACTIONATION PROTOCOL. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvob21zl4o/v1>



KEYWORDS

insoluble proteins, soluble proteins, TX-100 FRACTIONATION, brain

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited


CREATED

Jun 03, 2022

LAST MODIFIED

Jul 07, 2022

OWNERSHIP HISTORY

Jun 03, 2022  renuka.sJun 15, 2022  jbalster

PROTOCOL INTEGER ID

63788

MATERIALS TEXT

Materials

- phosphatase
- protease inhibitors

A	B
1x TNE (TNE) + Inhibitors	1 mL
2x TNE	500 µL
100x HALT inhibitor	10 µL
100x HALT EDTA	10 µL
H2O to final	480 µL





A	B
Complete 2x TNE (2x CTNE) + Inhibitors	1 mL
2x TNE	500 µL
1% SDS	100 µL
0.5% NP-40	100 µL
0.5% DOC	100 µL
100x HALT inhibitor	10 µL
100x HALT EDTA	10 µL
H2O to final	180 µL

Thaw Tissue

8h

1 






8h

If tissue stored at  **-80 °C**, place in  **-20 °C** for at least  **04:00:00** or  **Overnight** to thaw prior to homogenization.

Prepare 1X TNE




- 2 Prepare 1X TNE with phosphatase and protease inhibitors.


Homogenize Tissue

- 3 Weigh tissue out in mg.
- 4 Add in 10 volumes of 1X TNE.
 1.  **10 µL** of TNE per  **1 mg** of tissue.
 2. Ex:  **50 mg** tissue =  **500 µL** of 1X TNE.
- 5 Either by mechanical (Dounce) or homogenizer machine, homogenize tissue gently and  **On ice**.
- 6 This is TNE crude lysate (no detergents).

Prepare for Soluble v Insoluble

3h 20m 36s

- 7 Take specific volume of tissue in TNE and add in equal volume of 1X TNE w/ 2% Triton X-100 (Tx100).
 - Ex.  **150 µL** of TNE tissue +  **150 µL** 1X TNE+2% Tx100.
- 8 Sonicate @  **4 °C**.

8.1 3 pulses:  **00:00:10** ON /  **00:00:02** OFF.

12s

9



Spin down.

9.1 Option 1:  **16000 x g, 4°C, 00:15:00** .







15m

9.2 Option 2:  **20000 x g, 4°C, 01:00:00** .

1h

10 After spin:

10.1 **Supernatant = soluble**. Save supernatant and add equal volume of ^{25m 12s}complete TNE (cTNE).

1. Sonicate @  **4 °C** : 3 pluses:  **00:00:10** ON /  **00:00:02** OFF.
2. Boil:  **00:10:00** at  **95 °C** .
3. Spin down:  **16000 x g, 4°C, 00:15:00** .
4. Supernatant from this is the **soluble fraction**.

10.2 **Pellet = insoluble**.

11

Wash pellet in  **150 µL** of 1X TNE+1% Tx100.

Note: Same volume that was used above from TNE.



11.1

Resuspend pellet via pipette.

11.2

1h 15m

Spin down.

1. Option 1:  **16000 x g, 4°C, 00:15:00** .
2. Option 2:  **20000 x g, 4°C, 01:00:00** ( **25000 x g**) .

12 Resuspend pellet in  **75 µL** to  **100 µL** of cTNE.

13 Sonicate @  **4 °C** : 3 pluses:  **00:00:10** ON /  **00:00:02** OFF. 12s

14 Boil:  **00:10:00** at  **95 °C** . 10m

15  15m

Spin down:  **16000 x g, 4°C, 00:15:00** .

16 Supernatant from this is the **insoluble fraction**.