



VERSION 2

APR 12, 2023

OPEN  ACCESS

Protocol Citation: Haim Barr, Noa Lahav 2023. SARS-COV-2 Main Protease (Mpro) Fluorescence Dose Response. **protocols.io** <https://protocols.io/view/sars-cov-2-main-protease-mpro-fluorescence-dose-response-csbwwape> Version created by ASAP Discovery

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Apr 03, 2023**Last Modified:** Apr 12, 2023**PROTOCOL integer ID:**
79958

SARS-COV-2 Main Protease (Mpro) Fluorescence Dose Response V.2

Haim Barr^{1,2}, Noa Lahav^{1,2}¹The Weizmann Institute of Science; ²ASAP Discovery

ASAP Discovery

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

This is a functional, biochemical assay used to identify treatments for viral infectious disease in SARS-COV-2 Main Protease.

Utilizing a direct enzyme activity measurement method, the experiment was performed in a 384-well plate reading the fluorescence intensity. This assay tested the mode of action of inhibition.

GUIDELINES

Note: Inhibitor compounds stock concentration is **20 mM**. Compounds are pre-dispensed into 384 plates and stored at -20 C until use.

Plate Information:

Total Assay Volume: 20 μ L

Compounds Top Assay Concentration: 100 μ M

Dilution Factor: 2

Dose Response Points: 12

Number of Replicates: 2

Backfill with DMSO: Yes

MATERIALS

Assay Buffer Reagents (Concentration listed is the final concentration within the plate)

- [M] 20 millimolar (mM)
⊗ HEPES Buffer (pH 7.3) Fisher Scientific Catalog #BP299-1 (or similar)
- [M] 50 millimolar (mM)
⊗ Sodium Chloride Fisher Scientific Catalog #S271 (or similar)
- [M] 10 % volume ⊗ Glycerol Contributed by users Catalog #G5516 (or similar)
- [M] 0.01 % volume
⊗ Tween 20 Bio-Rad Laboratories Catalog #170-6606-MSDS (or similar)
- [M] 1 millimolar (mM) ⊗ TCEP HCl P212121 Catalog #SV-TCEP (or similar)

(all components are added fresh to the assay buffer before each experiment)

Additional Reagents

- [M] 5 nanomolar (nM) SARS Mpro Enzyme*

***Note:** Enzyme original stock was originally [M] 710 micromolar (μM) and was diluted to create smaller aliquots of [M] 20000 nanomolar (nM) (Storage Buffer was : 50 mM Tris pH 7.5, 1 mM DTT, 50 mM NaCl, 1 mM EDTA, 50% Glycerol). The 20000 nM aliquots [M] 20000 micromolar (μM) SARS Substrate were then diluted with **assay buffer** to [M] 10 nanomolar (nM), before each experiment. It was then diluted to be [M] 100 micromolar (μM) with **DMSO**.

- [M] 375 nanomolar (nM) SARS Substrate*

***Note:** SARS Substrate ([5-FAM]-AVLQSGFR-[Lys(Dabcyl)-K-amide] was dissolved in DMSO with an original concentration of [M] 20000 micromolar (μM). Following an intermediate dilution with DMSO, the SARS Substrate was diluted to a final concentration of 750 nM in assay buffer when it was first loaded into the Multi-Drop Combi Reagent Dispenser. When the entire plate was fully dispensed, the final concentration for SARS Substrate was then 375 nM





SAFETY WARNINGS





Please be sure to wear proper Personal Protective Equipment (PPE) while performing this experiment.

Prepare 384 Well Plate

15m

- 1 **PRIME** with **Assay Buffer** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes are filled completely.
 - 1.1 **DISPENSE**  10 µL Assay Buffer to Columns **1** and **23** of assay plate
 - **Note:** These will represent the ***inhibitor control columns*** (Contain: Substrate, Assay Buffer, DMSO, **no experimental compounds**)
 - 1.2 **EMPTY** Multi-Drop Combi Tube Dispensing Cassette
- 2 **PRIME** with  10 nanomolar (nM) SARS MPro by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes are filled completely.
 - **Note:** Be sure to cycle dispensing several times on a clean plate lid
- 2.1 **DISPENSE**  10 µL  10 nanomolar (nM) SARS MPro to Columns **2 through 22** and Column **24**

Note:

 -  10 nanomolar (nM) SARS Mpro is two times the final concentration for the assay. It is diluted to be a final concentration of  5 nanomolar (nM) SARS Mpro
 - Column 2 and Column 24 are ***neutral control columns*** (Contain: Enzyme, Substrate, DMSO, no compounds)

2.2 **EMPTY** Multi-Drop Combi Tube Dispensing Cassette (by selecting the **EMPTY** button on the Combi Dispenser until the tubes of the cassette are emptied.)

- **Note:** Discard the 100 nM SARS Mpro discharged from the cassette.

3 **CENTRIFUGE** 15000 rpm, Room temperature, 00:01:00 the plate to remove bubbles

1m

4 **INCUBATE** plate for 00:15:00 at Room temperature

15m

5 **PRIME** with **Assay Buffer** by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes are filled completely. Then, immediately **EMPTY** the Multi-Drop Combi Tube Dispensing Cassette.

6 **PRIME** with 750 nanomolar (nM) SARS Substrate by Multi-Drop Combi Tube Dispensing Cassette by selecting the **PRIME** button on the Combi Dispenser until the tubes are filled completely.

- **Note:** Be sure to cycle dispensing several times on a clean plate lid

7 **DISPENSE** 10 µL 750 nanomolar (nM) SARS Substrate to Columns 1 through 23 (the full plate)

Note:

- 750 nanomolar (nM) SARS Substrate is two times the final concentration for the assay. It is diluted to be a final concentration of 375 nanomolar (nM) SARS Substrate

8 **CENTRIFUGE** 15000 rpm, Room temperature, 00:01:00 the plate to remove bubbles

1m

9 **INCUBATE** plate for 00:30:00 at Room temperature

30m

⚠ Make sure the plate is protected from light

- **Recommended:** Clean the Multi-Drop Combi Reagent Dispenser during this incubation step

Read Plate Fluorescence

- 10 **READ** and **RECORD** the plate Relative fluorescence units (RFU) via the "**SARS Endpoint protocol**" on the **PHERASTAR FS Control Software**.

Equipment

PHERASTAR FS	NAME
Microplate reader	TYPE
BMG LABTECH	BRAND
0471B0001A	SKU
https://www.bmglabtech.com/en/pherastar-fsx/?utm_term=pherastar%20plate%20reader&utm_campaign=usa.roi.products&utm_source=adwords&utm_medium&gclid=Cj0KCQjw8qmhBhCIARIsANAtbodGRjigZtEYwcoMXUtxsLn25xp4gjKra3ZNt9jLh9-FwOoFR_5EUHUaAlkREALw_wcB	LINK

Expected result

Gain 300 should yield ~10,000 RFU in full reaction and ~6,000 RFU in Buffer Control