

🌐 Qiagen AllPrep 96 DNA/RNA protocol for bee abdomens (FFAR) V.1

Jocelynz ¹

¹University of Oregon

1 ▼

Oct 18, 2021

1



protocol .

Ponisiolab

Jocelynz

Qiagen AllPrep 96 DNA/RNA protocol for bee abdomens (FFAR)

Jocelynz 2021. Qiagen AllPrep 96 DNA/RNA protocol for bee abdomens (FFAR).
protocols.io
<https://protocols.io/view/qiagen-allprep-96-dna-rna-protocol-for-bee-abdomen-bxxippke>



protocol ,

Sep 02, 2021

Oct 18, 2021

52938

Before beginning protocol

1

⊗ DTT (2M) Fisher

Scientific Catalog #50-198-452

must be added to

⊗ Buffer

RLT Qiagen Catalog #79216

before use. Add **20 μL** [M] **20 Molarity (M)**

⊗ DTT (2M) Fisher

Scientific Catalog #50-198-452

per **1 mL**

⊗ Buffer

RLT Qiagen Catalog #79216

. Stock solution of

⊗ DTT (2M) Fisher

Scientific Catalog #50-198-452

should be prepared fresh or frozen

⊗ Buffer

in single use aliquots. RLT Qiagen Catalog #79216

containing

⊗ DTT (2M) Fisher

Scientific Catalog #50-198-452

can be stored at room

temperature for up to 1 month.



⊗ DTT (2M) Fisher

Scientific Catalog #50-198-452

is highly toxic, review

SDS thoroughly and dispose of as hazardous waste



⊗ Buffer

RLT Qiagen Catalog #79216

⊗ Buffer

RW1 Qiagen Catalog #1053394

⊗ Buffer

AW1 Qiagen Catalog #19081

contain guanidine salt - do not mix with

bleach



⊗ Buffer

RPE Qiagen Catalog #1018013

concentrate can be disposed in the

sink in small volumes – no more than 100 g solute per day, flush with 10-20 fold water thoroughly after sink disposal.



⊗ Buffer

Manage and dispose of AW2 Qiagen Catalog #19072

as hazardous

waste

Prepare samples

18m

- 2 UV sterilize supplies for 2 96 well plates: stainless steel beads, 8 reagent reservoirs, 2 collection microtube plate (light blue) and strip caps (clear), 2 brand new s-blocks, 4 more s blocks, 4 elution microtube plates (dark blue) and strip caps (white), and 96 well microplates (for aliquots). Label foil lids for aliquots. Turn on water bath and set to **70 °C**
- 3 Add 1 stainless steel bead to each well of collection microtube plate
- 4 Working on dry ice to prevent samples from thawing, add 1 bee abdomen to each collection microtube. Enter each sample ID into your plate map as you work and keep used vials stored in order.
- 5 [Buffer](#)
Remove from dry ice and immediately add [RLT Qiagen Catalog #79216](#)

Check reagent bottle to ensure

[DTT \(2M\) Fisher](#)

[Scientific Catalog #50-198-452](#)

has been added before use

- 5.1 For larger bee abdomens such as melissodes, add **500 µL**

[Buffer](#)

[RLT Qiagen Catalog #79216](#)

- 5.2 For smaller bee abdomens such as lasioglossum, add **350 µL**

[Buffer](#)

[RLT Qiagen Catalog #79216](#)

Seal with strip caps and place in tissue lyser, placing a folded kim wipe between strip cap tops

6 and adapter to ensure a tight seal

TissueLyser II

Bead Mill

QIAGEN

85300



TissueLyser Adapter Set 2 x 96

Adapter set

Qiagen

69984



7 Lyse at **25 Hz** for **00:02:00**

2m

Depending on size of sample, you may need to adjust speed and time

8 Remove plates from lyser and adapters, rotate plates 180 degrees, and place back in tissue lyser

9 Lyse at **25 Hz** for **00:02:00**

2m

Depending on size of sample, you may need to adjust speed and time

10 Centrifuge **3486 x g, 20-25°C, 00:07:00**

7m


Eppendorf™ 5810R Centrifuge
Centrifuge

Eppendorf 02-262-8187 [↗](#)


- 11 Open a brand new s-block and place DNA filter plate on top

Do not reuse an old s-block for this step!

- 12 Transfer  **350 µL** supernatant to DNA filter plate and seal with AirPore tape

- 13 Centrifuge  **3478 x g, 20-25°C, 00:07:00** . Check that all liquid has flowed through columns, repeating centrifugation if needed 7m

**** DO NOT DISCARD FLOWTHROUGH ****



- 14 Place DNA filter plate on new s-block and put aside at room temperature to  **4 °C** for later use

You may reuse and old (sterilized) s-block for this step and all future steps

RNA purification 38m

- 15 Place RNA filter plate on s-block













- 16 Add  **350 µL**  **70% Ethanol** **Contributed by users** to s-block containing flowthrough from step 13. Pipette up and down 3 times to mix.
- 17 Transfer samples to RNA filter plate. Seal with AirPore tape
- 18 Centrifuge  **3486 x g, 20-25°C, 00:07:00** , discard flowthrough 7m
- 19  **Buffer**
Add  **800 µL**  **RW1** **Qiagen Catalog #1053394** and seal with AirPore tape
- 20 Centrifuge  **3486 x g, 20-25°C, 00:07:00** , discard flowthrough 7m
- 21  **Buffer**
Add  **800 µL**  **RPE** **Qiagen Catalog #1018013** and seal
- 22 Centrifuge  **3486 x g, 20-25°C, 00:16:00** , discard flowthrough 16m
- 23 Place RNA filter plate on top of elution microtube plate
- 24 Add  **60 µL** RNase-free water directly to center of each filter membrane and seal

Use  **45 µL** -  **70 µL** RNase-free water depending on desired DNA yield/concentration

- 25 Incubate at room temperature ⌚ 00:01:00 1m
- 26 Centrifuge 🌀 3486 x g, 20-25°C, 00:07:00 7m
- 27 Repeat steps 23-26, eluting into same microtube plate
- 28 Distribute product into 📄 10 µL aliquots and seal with labeled foil lids. Seal remaining product in elution microtube plates with white strip caps
- 29 Store at 🌡 -20 °C - 🌡 -80 °C

DNA purification 35m

- 30 🔗 Buffer
Place EB Qiagen Catalog #19086 in 🌡 70 °C water bath
- Place reagent bottle in water bath float - do not fully submerge
- 31 Get out your DNA filter plate and s-block from step 14
- 32 🔗 Buffer
Add 📄 800 µL AW1 Qiagen Catalog #19081 and seal with AirPore tape
- 33 Centrifuge 🌀 3486 x g, 00:07:00 , discard flowthrough 7m

- 34 Add  **800 µL**  **Buffer**  **AW2 Qiagen Catalog #19072** and seal
- 35 Centrifuge  **3486 x g, 00:16:00** , discard flowthrough 16m
- 36 Place DNA filter plate on top of new elution microtube plate
- 37 Add  **50-100 µL** preheated  **Buffer**  **EB Qiagen Catalog #19086** directly to the center of each filter membrane and seal
- 38 Incubate at room temperature for  **00:05:00** 5m
- 39 Centrifuge  **3486 x g, 00:07:00** 7m
- 40 Repeat steps 36-39, eluting into same plate
- 41 Distribute product into  **10 µL** aliquots and seal with labeled foil lids. Seal remaining product in elution microtube plates with white strip caps
- 42 Store at  **-20 °C** -  **-80 °C**