



JUL 31, 2023

OPEN ACCESS



Protocol Citation: Nimalka M Weerasuriya 2023. Pythium Zoospore Production Soaking Solution. [protocols.io](https://protocols.io/view/pythium-zoospore-production-soaking-solution-cxv6xn9e) <https://protocols.io/view/pythium-zoospore-production-soaking-solution-cxv6xn9e>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development
We are still developing and optimizing this protocol

Created: Jul 28, 2023

Last Modified: Jul 31, 2023

PROTOCOL integer ID:
85662

Pythium Zoospore Production Soaking Solution

Nimalka M
Weerasuriya¹

¹Oklahoma State University



Nimalka M Weerasuriya
Oklahoma State University, USDA Agricultural Research Servic...

DISCLAIMER

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to [protocols.io](#) is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with [protocols.io](#), can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

ABSTRACT

Creation of soaking solutions for *Pythium myriotylum* to be used for large-scale zoospore production.

MATERIALS

Soaking solutions:

- 1 L beaker
- 1 g CaCO₃
- Whatman #1 filter
- 1 N KOH – pH adjustment
- 1 N HCl – pH adjustment
- Sucrose

Testing Zoospore Culture:

- Haemocytometer
- Microscope and slides
- Counter
- 0.08% Methylene blue

Preparation

- 1 Have mature colonies of verified *Pythium myriotylum* growing on CMA or 1.5-2% WA. Colony maturity ~7 days, with visible oospores.

CITATION

Jones, B. L., & Woodard, K. E (1986). A Technique for Evaluating Peanut Germ Plasm for Resistance to *Pythium myriotylum*. *Plant Disease*, 70(11), 1038–1043.

LINK

<https://doi.org/10.1094/PD-70-1038>

Soaking Solutions

- 2 Make soaking solutions 1, 2 and 3:

CITATION

Nyochembeng, L. M., Pacumbaba, R. P., & Beyl, C. A (2002). Calcium Enhanced Zoospore Production of *Pythium myriotylum* in vitro. *Journal of Phytopathology*, 150(7), 396–398.

LINK

<https://doi.org/10.1046/J.1439-0434.2002.00759.X>

2.1

Soaking Solution #1, [M] 0.01 Molarity (m) Ca++ (1 L):
1 L RO water at pH 7 in 1 L beaker
 Add 1 g CaCO₃
 Filter through Whatman #1

pH is adjusted from pH 8.5 to pH 10.5 with [M] 1 Mass Percent HCl to duplicate soaking solution in original methods (Pacumbaba, cite)

2.2

Soaking Solution #2, [M] 0.01 Molarity (m) Ca++ + [M] 0.001 Molarity (m) sucrose (500 mL):
 <missing info>

pH is adjusted to pH 7 with [M] 1 Mass Percent HCl

Note

(1 M (mol/l) = 1 N for an acid that releases 1 proton* when dissolved in water)

2.3

Soaking Solution #3, [M] 0.001 Molarity (m) sucrose (500 mL):
 <missing info>

pH is adjusted to pH 7 with [M] 1 Mass Percent HCl

3 Autoclave for 20 minutes liquid cycle

Protocol



NAME

Sterilizer (Consolidated)



CREATED BY



Nimalka M Weerasuriya

PREVIEW

Testing Zoospore Production

1d



- 4 Take  30 mL of each Soaking Solution  Soaking Solution Preparation onto surface of 7-day cultures in plate.

- 5 Incubate under light at  Room temperature for  24:00:00 .
Check for abundant sporangia that will appear after immersion.

1d

- 6 At 1.5 up to 4 h every 30 minutes:



- Take  80 μ L soaking liquid from each immersion to a microcentrifuge tube with  20 μ L 0.08% methylene blue . Gently invert to stain (or gently vortex) and immobilize zoospores.

- 7 Take  10 μ L of liquid into haemocytometer and examine under 40x to quantify.