



Sep 16, 2024

Ethanol precipitation of small DNA fragments

DOI

dx.doi.org/10.17504/protocols.io.n2bvj7jnvk5w/v1

Eva Petrova¹, Roey Angel¹

¹Soil and Water Research Infrastructure

SoWa RI Anaerobic and Molecular Microbiology (public)

Tech. support email: eva.petrova@bc.cas.cz



Eva Petrova

Soil and Water Research Infrastructure

OPEN  ACCESS



DOI: **dx.doi.org/10.17504/protocols.io.n2bvj7jnvk5w/v1**

Protocol Citation: Eva Petrova, Roey Angel 2024. Ethanol precipitation of small DNA fragments. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n2bvj7jnvk5w/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 24, 2018

Last Modified: September 16, 2024

Protocol Integer ID: 12439



Abstract

This protocol is for a simple ethanol precipitation of small fragments and/or small amounts of DNA. Note that this protocol simply concentrates your sample and removes enough salts/enzymes for ligation to be successful. All DNA fragments from your digest will still be present in your pellet.

See [here](#) for a discussion of theory behind ethanol precipitation of DNA fragments.

Guidelines

3M Sodium Acetate, 0.11M MgCl_2 pH 5.5

24.6 g sodium acetate

0.952g MgCl_2 anhydrous

Dissolve in about 60 ml DDW Adjust pH with glacial acetic acid (>8 ml), Dilute to 100 ml with DDW and filter sterilize.

Materials

MATERIALS

 Pellet Paint® NF Co-Precipitant **Sigma Aldrich Catalog #70748-3**

Protocol materials

 Pellet Paint® NF Co-Precipitant **Merck MilliporeSigma (Sigma-Aldrich) Catalog #70748-3** Materials, Step 3

Before start


- Absolute Ethanol (100% = 200 proof) at -80°C
- 70% ethanol at -80°C (don't cool ethanol for too long or water will freeze)
- Tabletop centrifuge
- -80°C freezer





- 1 Transfer the sample to 1.5 mL eppendorf tube.
- 2 Adjust the salt conc. of the sample to 0.3M sodium acetate and 0.01M MgCl_2 .


[M] 0.3 Molarity (M) sodium acetate


[M] 0.01 Molarity (M) MgCl_2
- 3 Add 1 μL of Pellet Paint NF.

 1 μL Pellet Paint NF

 Pellet Paint® NF Co-Precipitant **Sigma Aldrich Catalog #70748-3**
- 4 Add 1 μL of non labeled genomic DNA e.g. from E.coli (this acts as a carrier and increases precipitation yield).


 1 μL genomic DNA (E.coli)
- 5 Add 2 volumes cold absolute ethanol to sample.
- 6 Incubate 2-3 hr at -80°C . The long incubation time is critical for small fragments.

 -80°C

 02:00:00
- 7 Centrifuge for 30 minutes at 0°C at maximum speed (generally $>10000\text{ g}$ at least).

Command**max. speed**


CENTRIFUGE

 0°C  00:30:00

- 8 Remove supernatant.



- 9 Wash with 750-1000 μL cold 70% ethanol. Another critical step for small fragments under 200 base pairs. Washing involves adding the ethanol and inverting several times.

 750 μL cold 70% ethanol


- 10 Centrifuge for 30 minutes at 0°C at maximum speed (generally >10000 g at least).

Command

max. speed

CENTRIFUGE

 0 °C

 00:30:00

- 11 Repeat washing procedure and centrifuge again.

 go to step #9

- 12 Let air dry on benchtop.

- 13 Resuspend in an appropriate volume of H_2O .