



May 11, 2021

Extraction of flagellum basal body complexes from *Buchnera aphidicola*, an endosymbiont of aphids

mschepers¹¹The University of Texas at Austin Institute for Cell and Molecular Biology**1** Works for me dx.doi.org/10.17504/protocols.io.bs5wng7e mschepers

ABSTRACT

Buchnera aphidicola is an intracellular bacterial symbiont of aphids and maintains a small genome of only 600 kbps. *Buchnera* is thought to maintain only genes relevant to the symbiosis with its aphid host. Curiously, the *Buchnera* genome contains gene clusters coding for flagellum basal body structural proteins and for flagellum type III export machinery. These structures have been shown to be highly expressed and present in large numbers on *Buchnera* cells. No recognizable pathogenicity factors or secreted proteins have been identified in the *Buchnera* genome, and the relevance of this protein complex to the symbiosis is unknown. Here, provide a procedure for isolating the *Buchnera* flagellum basal body proteins from the cellular membrane of *Buchnera*. This will facilitate studies of the structure and function of the *Buchnera* flagellum, and its role in this model symbiosis.

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0245710>

DOI

dx.doi.org/10.17504/protocols.io.bs5wng7e

EXTERNAL LINK

<https://doi.org/10.1371/journal.pone.0245710>

PROTOCOL CITATION

mschepers 2021. Extraction of flagellum basal body complexes from *Buchnera aphidicola*, an endosymbiont of aphids. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bs5wng7e>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Mar 09, 2021

LAST MODIFIED

May 11, 2021

PROTOCOL INTEGER ID

48022

Required Reagents

- NaClO solution:** 0.5%; 83.5 mL 6% NaClO in 916.5 mL sterile H₂O
Buffer A: 25mM KCl, 35mM Tris base, 10mM MgCl₂, 250mM anhydrous EDTA in 950 mL sterile H₂O. Add pellets of

KOH until anhydrous EDTA is dissolved/solution becomes clear. Adjust pH to 7.5. Autoclave, then add 50 mL filter-sterilized 2.5M sucrose.

Sucrose solution: 0.5M sucrose, 100mM Tris-HCl. Raise pH to 8 and filter sterilize.

Lysozyme solution: Prepare 2 mg/mL egg white lysozyme in sterile H₂O.

EDTA solution: 100mM anhydrous EDTA in sterile H₂O. Add KOH until EDTA is dissolved/solution becomes clear. Adjust pH to 7.5.

Triton X-100 solution: 10% w/v Triton X-100 detergent in sterile H₂O.

MgSO₄ solution: 100mM MgSO₄ in sterile H₂O.

RNase-free DNase: 1 mg/mL in sterile H₂O.

NaOH, 1N: Dissolve NaOH pellets in sterile H₂O.

TET Buffer: 10mM Tris-HCl, 5mM EDTA, 0.1% w/v Triton X-100. Add KOH or NaOH to dissolve anhydrous EDTA, then adjust pH to 8.0.

Liberating Buchnera from Aphids

- 2 Grow aphids from birth to fourth instar (10 days) and harvest from plants. For this protocol, we used 3-5g live aphids.
- 3 Load aphids into sterile fine mesh tea infuser or other sterile straining device.
- 4 Submerge into 0.5% NaClO solution for 1 min.
- 5 Transfer and submerge into sterile MQ H₂O at 4°C for 1 minute.
- 6 Repeat step 4 (Transfer and submerge into sterile MQ H₂O for 1 minute).
- 7 Remove aphids from strainer, placing in a mortar kept at 4°C. Add 20mL buffer A + EDTA (pH 7.5) and homogenize aphids w/ pestle.
- 8 Pour homogenized aphids into sterile vacuum filter equipped with a 100uM nylon filter. Collect homogenate on ice.
- 9 Remove nylon filter and place aphid filter cake back into mortar. Repeat homogenization with an additional 20mL buffer A.
- 10 Repeat filtration with new 100mM nylon filter, combining filtrate.
- 11 Centrifuge 100uM filtrate for 10 min at 1500g at 4°C in a pre-cooled centrifuge.

- 12 Discard supernatant. Gently resuspend pellet in 40mL cold buffer A with sterile spatula or wide-bore pipette tip. Do not vortex.
- 13 Vacuum filter resuspension through a 20uM nylon filter. Collect filtrate on ice.
- 14 Immediately continue filtration through a 11uM nylon filter, collecting filtrate on ice.
- 15 Filter through a 5uM nylon filter, collecting filtrate on ice.
- 16 Centrifuge 5 uM filtrate for 30 minutes at 1500g at 4C in a pre-cooled centrifuge.
- 17 Discard supernatant. Gently resuspend pellet in 100 mL sterile sucrose solution. Do not vortex.

Extracting Flagellum complexes from Buchnera

- 18 Optional: Use light microscope to check integrity of Buchnera cells. Cells are spherical, averaging 3uM in size.
- 19 Begin to gently stir on ice.
- 20 Add 5mL lysozyme (2 mg/mL, freshly prepared in H₂O at 4°C).
- 21 Add 10mL 100mM EDTA solution pH 7.5 and incubate on ice gently spinning for 40 min.
- 22 Add 10mL 10% Triton X-100 (in H₂O) to lyse cells. Solution should turn from turbid to clear.
- 23 Once solution is transparent, add 100mM MgSO₄ and 1 mg/ml RNase-free DNase. Incubate at room temperature, gently stirring, allowing DNases to degrade cellular DNA. Solution should become less viscous after 30 min, this step can be extended to an overnight incubation at 4°C if DNases require more time.
- 24 Raise pH of the solution to pH 10 with 1N NaOH.

- 25 Spin solution at 5000g in a pre-cooled 4C centrifuge three times, saving the supernatant and decanting into a new tube each time.
- 26 After three spins, transfer supernatant into a polyallomer centrifuge tube and spin at 30,000g for 1h at 4C in a pre-cooled centrifuge.
- 27 Carefully decant or draw off supernatant and discarded. Cover pellet with TET (10mM Tris-HCl, 5mM EDTA, 0.1% Triton X-100, pH 8.0) and leave overnight at 4C to soften protein pellet. Buchnera protein is ready for downstream analyses.