



Sep 06, 2022

Preparing plasmids for nucleofection of hPSCs

In 1 collection

Hanqin Li¹, Oriol Busquets², Steven Poser², Dirk Hockemeyer¹, Frank Soldner²

¹University of California, Berkeley; ²Albert Einstein College of Medicine

1 Works for me

Share

dx.doi.org/10.17504/protocols.io.b4pdqvi6



Devin E Snyder

ABSTRACT

This protocol describes the standard procedure for the preparation of plasmids to be delivered into human pluripotent stem cells (hPSCs) using nucleofection.

General notes

1. Throughout this protocol, the term hPSC is used to collectively refer to both hiPSCs and hESCs. All described procedures have been tested and work equally well for hiPSCs and hESCs.

DOI

dx.doi.org/10.17504/protocols.io.b4pdqvi6

PROTOCOL CITATION

Hanqin Li, Oriol Busquets, Steven Poser, Dirk Hockemeyer, Frank Soldner 2022. Preparing plasmids for nucleofection of hPSCs. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b4pdqvi6>



FUNDERS ACKNOWLEDGEMENT

Aligning Science Across Parkinson's
Grant ID: ASAP-000486

COLLECTIONS ⓘ



Nucleofection (Amaxa) and electroporation (Biorad) of hPSCs

KEYWORDS

ASAPCRN

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Feb 03, 2022

LAST MODIFIED

Sep 06, 2022

PROTOCOL INTEGER ID

57797


PARENT PROTOCOLS

Part of collection

[Nucleofection \(Amaxa\) and electroporation \(Biorad\) of hPSCs](#)

MATERIALS TEXT

Item	Vendor	Catalog #
Synthetic pegRNAs	IDT or Synthego	
Synthetic sgRNAs	Synthego	
pCMV-PE2	Addgene	132775

- 1 Thaw plasmids  **On ice**
- 2 In each nucleofection, use 1 µg total of plasmid
- 3 For prime editing **PE2 strategy**, use:
 - 500 ng pCMV-PE2
 - 500 ng pU6-pegRNA

- 4 For prime editing **PE3 strategy**, use:
 - 500 ng pCMV-PE2
 - 330 ng pU6-pegRNA
 - 170 ng pBPK1520-ngRNA
- 5 Pipet the proper amount of each plasmid into a microcentrifuge tube. If the volume is too small, dilute in autoclaved water first.