



Mar 14, 2022

Workflow for bulk RNAseq of human placenta

Scott Lindsay-Hewett¹, Valentina Stanley¹, Mana Parast¹, Louise Laurent¹¹University of California, San Diego

1

dx.doi.org/10.17504/protocols.io.b59sq96e

Human BioMolecular Atlas Program (HuBMAP) Method Development Community
Tech. support email: Jeff.spraggins@vanderbilt.edu



Scott Lindsay-Hewett

Described here is the workflow used by the Female Reproductive Tissue Mapping Center at UCSD to generate RNAseq data from human placenta.

DOI

dx.doi.org/10.17504/protocols.io.b59sq96e

Scott Lindsay-Hewett, Valentina Stanley, Mana Parast, Louise Laurent 2022.
Workflow for bulk RNAseq of human placenta. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b59sq96e>



protocol ,

Mar 14, 2022

Mar 14, 2022

59410

See accompanying protocols.

Tissue preparation

- 1 As soon as possible after Cesarean section or vaginal delivery, prepare tissue according to the following protocol:

[Human Placenta Tissue Collection and Preservation Methods - UCSD Female Reproductive TMC](#)

For this protocol, use tissue that has been stored in RNAlater.

Total RNA isolation

- 2 Isolate total RNA using a bead beater to disrupt the tissue, followed by organic extraction and ethanol precipitation. Use the following protocol:

[Total RNA extraction from frozen placenta tissue](#)

After passing quality control, proceed to library construction.

Library construction

- 3 Construct libraries using the KAPA RNA HyperPrep Kit with RiboErase (HMR), according to the following protocol:

[Library construction for human placenta bulk RNAseq](#)

After passing quality control, proceed to sequencing.

Sequencing

- 4 For HuBMAP bulk RNAseq samples, the multiplexed pool was sequenced on a NovaSeq 6000 S4 lane using a 100bp paired-end run configuration. Reads were aligned using STAR, and transcript abundances were quantified using RSEM.