



MAY 30, 2023

SJB Artificial Seawater Medium Protocol

Michael
Shelby J Barnes¹, Henson¹, Cameron
Thrash¹

¹University of Southern California

 Shelby J Barnes

ABSTRACT

SJB Artificial Seawater Medium Protocol

SJB Artificial Seawater Medium is an amended version of JW and MWH generations of Thrash Lab artificial seawater medium.

MATERIALS

Sterilized Pyrex Bottle

0.2um PES filters

Base Salts

A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
Salinity			1		2		3		4		4.25		4.5		5	
Base Salts	Chemical Formula	Formula Weight	g/L	Final M	g/L	Final M	g/L	Final M	g/L	Final M	g/L	Final M	g/L	Final M	g/L	Final M
sodium chloride	NaCl	58.443	23.5	0.402101	15.7	0.268067	7.83333	0.134033667	3.91666667	0.067016833	3.012820513	0.05155141	2.008547009	0.034367607	1.004273504	0.017183803
potassium chloride	KCl	74.55	0.746	0.010007	0.497	0.06671333	0.24866667	0.03335667	0.124333333	0.016678333	0.095641026	0.01282949	0.063760684	0.00855299	0.031880342	0.0042765

OPEN ACCESS

DOI:
dx.doi.org/10.17504/protocols.io.rm7vzy615lx1/v1

Protocol Citation: Shelby J Barnes, Michael Henson, Cameron Thrash 2023. SJB Artificial Seawater Medium Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.rm7vzy615lx1/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Apr 14, 2022

Last Modified: May 30, 2023

PROTOCOL integer ID:
60802

Keywords: Artificial Seawater

A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q
sodium bicarbonate	NaHCO ₃	84.01	0.84	0.00999	0.84	0.00666	0.42	0.00333	0.21	0.001665	0.161538462	0.001281923	0.107692308	0.000854615	0.0053846154	0.000427308
sodium sulfate	Na ₂ SO ₄	142.04	4.27	0.030062	2.85	0.02041333	1.42333333	0.01020667	0.71166667	0.05010333	0.547435897	0.03854103	0.364957265	0.02569402	0.182478632	0.01284701
sodium bromide	NaBr	102.89	0.83	0.008067	0.55	0.005378	0.27666667	0.02689	0.13833333	0.013445	0.106410256	0.01034231	0.070940171	0.00689487	0.035470085	0.00344744
boric acid	H ₃ BO ₃	61.83	0.026	0.000421	0.017	0.00280667	0.00866667	0.00140333	0.00433333	7.01667E-05	0.00333333	5.39744E-05	0.00222222	3.59829E-05	0.00111111	1.79915E-05
strontium chloride	SrCl ₂	158.53	0.0142	0.00009	0.009467	0.00006	0.00473333	0.00003	0.00236667	0.00015	0.001820513	1.15385E-05	0.001213675	7.69231E-06	0.000606838	3.84615E-06
sodium fluoride	NaF	41.99	0.0023	0.000055	0.0015	3.6667E-05	0.00076667	1.8333E-05	0.00038333	9.16667E-06	0.000294872	7.05128E-06	0.000196581	4.70085E-06	9.82906E-05	2.35043E-06
magnesium chloride heptahydrate	MgCl ₂ x 6H ₂ O	203.3	10.6	0.05214	7.1	0.03476	3.53333333	0.01738	1.76666667	0.00869	1.358974359	0.00684615	0.905982906	0.0445641	0.452991453	0.02228205
calcium chloride dihydrate	CaCl ₂ x 2H ₂ O	147.01	1.52	0.010339	1.013	0.006892667	0.50666667	0.03446333	0.25333333	0.01723167	0.194871795	0.01325513	0.12991453	0.00883675	0.064957265	0.00441838

Table 1. Base Salts

The range of salinities are shown in order of highest (#1) to lowest (#5).

Trace Metal Mix (100,000x Stock)

A	B	C	D	E
Compound	Chemical Formula	FW	g/100mL	Final M (nM)
manganese dichloride tetrahydrate	MnCl ₂ x 4H ₂ O	197.91	0.018	9
zinc sulfate monohydrate	ZnSO ₄ X H ₂ O	179.47	0.002	1
cobalt(II) chloride	CoCl ₂	129.839	0.001	0.5
sodium molybdate	Na ₂ MoO ₄	205.92	0.001	0.3
sodium selenite	Na ₂ SeO ₃	172.94	0.002	1
nickel(II) chloride	NiCl ₂	129.5994	0.001	1

Vitamin Mix (100,000x)

A	B	C	D	E
Compound	Chemical Formula	FW	g/100mL	Final M (nM)
B1/thiamine	C ₁₂ H ₁₇ ClN ₄ O ₄ S · HCl	337.27	1.69	500 nM
B2/riboflavin	C ₁₇ H ₂₀ N ₄ O ₆	376.36	0.0026	0.7 nM
B3/niacin	C ₆ H ₅ NO ₂	123.12	0.985	800 nM
B5/pantothenate	HOCH ₂ C(CH ₃) 2CH(OH)CON HCH ₂ CH ₂ CO ₂ · 1/2Ca	238.27	1.013	425 nM
B6/pyridoxine	C ₈ H ₁₁ NO ₃ · HCl	205.64	1.028	500 nM
B7/biotin	C ₁₀ H ₁₆ N ₂ O ₃ S	244.31	0.0098	4 nM
B9/folic Acid	C ₁₉ H ₁₉ N ₇ O ₆	441.4	0.0177	4 nM
B12/cyanocobalamin	C ₆₃ H ₈₈ CoN ₁₄ O ₁₄ P	1355.37	0.0095	0.7 nM
myo-inositol	C ₆ H ₁₂ O ₆	180.16	0.901	500 nM
4-aminobenzoic Acid	C ₇ H ₇ NO ₂	137.14	0.0823	60 nM

Amino Acid Mix (5,000x)



A	B	C	D	E
Compound	Chemical Formula	FW	g/100mL	Final M
L-arginine x HCl	$\text{HN}=\text{C}(\text{NH}_2)\text{NH}(\text{C}_6\text{H}_5)\text{CH}_2\text{CH}(\text{NH}_2)\text{COOH} \cdot \text{HCl}$	210.7	0.632	0.0000005999050783
L-cystine x 2HCl	$\text{C}_6\text{H}_{12}\text{N}_2\text{O}_4\text{S}_2 \cdot 2\text{HCl}$	313.22	0.1564	0.00000009986590895
L-histidine x HCl x H ₂ O	$\text{C}_6\text{H}_9\text{N}_3\text{O}_2 \cdot \text{HCl} \cdot \text{H}_2\text{O}$	209.65	0.21	0.0000002003338898
L-Isoleucine	$\text{C}_2\text{H}_5\text{CH}(\text{CH}_3)\text{CH}(\text{NH}_2)\text{CO}_2\text{H}$	131.17	0.2625	0.0000004002439582
L-leucine	$\text{CH}_3\text{CH}_2\text{CH}(\text{CH}_3)\text{CH}(\text{NH}_2)\text{CO}_2\text{H}$	131.17	0.262	0.0000003994815888
L-lysine x HCl	$\text{H}_2\text{N}(\text{CH}_2)_4\text{CH}(\text{NH}_2)\text{CO}_2\text{H} \cdot \text{HCl}$	182.69	0.3625	0.0000003968471181
L-methionine	$\text{CH}_3\text{SCH}_2\text{CH}_2\text{CH}(\text{NH}_2)\text{CO}_2\text{H}$	149.21	0.0755	0.0000001011996515
L-phenylalanine	$\text{C}_6\text{H}_5\text{CH}_2\text{CH}(\text{NH}_2)\text{CO}_2\text{H}$	165.19	0.165	0.0000001997699619
L-threonine	$\text{CH}_3\text{CH}(\text{OH})\text{CH}(\text{NH}_2)\text{CO}_2\text{H}$	119.12	0.238	0.000000399597045
L-tryptophan	$\text{C}_{11}\text{H}_{12}\text{N}_2\text{O}_2$	204.23	0.051	0.00000004994369094
L-tyrosine	$4\text{-(HO)C}_6\text{H}_4\text{CH}_2\text{CH}(\text{NH}_2)\text{CO}_2\text{H}$	181.19	0.18	0.0000001986864617
L-valine	$(\text{CH}_3)_2\text{CHCH}(\text{NH}_2)\text{CO}_2\text{H}$	117.15	0.234	0.0000003994878361

Miscellaneous Mix (20,000x)

A	B	C	D	E	F
Compound	Chemical Formula	FW	g/100mL	mL/100mL	Final M
L-glutamine	$\text{C}_5\text{H}_{10}\text{N}_2\text{O}_2$	146.14	0.14614		0.0000005
dextrose	$\text{C}_6\text{H}_{12}\text{O}_6$	180.16	0.18016		0.0000005
D-ribose	$\text{C}_5\text{H}_{10}\text{O}_5$	150.13	0.15013		0.0000005
sodium pyruvate	$\text{C}_3\text{H}_3\text{NaO}_3$	110.04	0.11004		0.0000005
sodium citrate	$\text{C}_6\text{H}_5\text{Na}_3\text{O}_7$	294.1	0.2941		0.0000005
oxaloacetic acid	$\text{C}_4\text{H}_4\text{O}_5$	132.07	0.13207		0.0000005

A	B	C	D	E	F
sodium acetate	C2H3NaO2	82.03	0.08203		0.0000005
sodium succinate dibasic hexahydrate	NaOOCCH2CH2COONa · 6H2O	270.14	0.27014		0.0000005
α-ketoglutaric acid	C5H6O5	168.08	0.16808		0.0000005
urea	CH4N2O	60.6	0.606		0.000005
glycerol	C3H8O3	92.09	0.09209		0.0000005
glycine betaine	C5H11NO2	153.61	0.15361		0.0000005
choline	(CH3)3N(Cl)CH2CH2OH	139.62	0.13962		0.0000005
sodium thiosulfate	Na2S2O3	158.11	0.15811		0.0000005
potassium cyanate	KOCN	81.11	0.0032444		2.00E-08
dmsO	C2H6OS	78.13	0.0062504	0.05682181818	4.00E-08
dmsp	C5H10O2S	134.967	0.01079736		4.00E-08
L-glycine	C2H5NO2	75.07	0.07507		0.0000005
nucleotides (dNTPs)	-	-	-	0.01	0.00000001
galactose	C6H12O6	180.16	0.18016		0.0000005
lactose monohydrate	C12H24O12	360.31	0.36031		0.0000005
glyoxylic acid sodium monohydrate	C2H3NaO4	114.03	0.11403		0.0000005
tmao	C3H9NO	75.11	0.07511		0.0000005
L-alanine	C3H7NO2	89.09	0.08909		0.0000005
L-asparagine	C4H8N2O3	132.12	0.13212		0.0000005
L-aspartic acid	C4H7NO4	133.1	0.1331		0.0000005
L-cysteine hydrochloride	C3H8ClNO2S	157.62	0.17563		0.0000005
L-glutamic acid	C5H9NO4	147.13	0.14713		0.0000005
L-proline	C5H9NO2	115.13	0.11513		0.0000005

A	B	C	D	E	F
L-serine	C3H7NO3	105.09	0.10509		0.0000005
sodium alginate	C6H9NaO7	216.12	0.21612		0.0000005
catalase	C9H10O3	2,000-5,000 U/mg	0.000542857		10 U
glucosamine	C6H13NO5	221.21	0.22121		0.0000005

Phosphate Mix (1,000x)

A	B	C	D	E	F
Compound	Chemical Formula	FW	g/100mL	mL/100mL	Final M
ortho phosphonate	C3H9O3P	124.076		0.002167	0.0000002
potassium phosphate monobasic	KH2PO4	136.09	0.068045		0.000005

Fatty Acid Mix (2,000,000x)

A	B	C	D	E	F
Compound	Chemical Formula	FW	g/100mL	mL/100mL	Final M
octanoic acid	CH3(CH2)6COOH	144.21		15.84725275	0.0000005
decanoic acid	CH3(CH2)8COOH	172.26	17.226		0.0000005
isobutyric acid	(CH3)2CHCO2H	88.11		9.274736842	0.0000005
butyric acid	CH3CH2CH2COOH	88.11		9.140041494	0.0000005
valeric acid	C5H10O2	102.13		10.87646432	0.0000005
ethanol	CH3CH2OH	46.068		54.86150459	0.000004698025433

Inorganic Nitrogen (2,000x)

A	B	C	D	E
Compound	Chemical Formula	FW	g/100mL	Final M
sodium nitrate	NaNO3	84.99	0.645924	0.000038
sodium nitrite	NaNO2	69	0.0276	0.000002

A	B	C	D	E
ammonium chloride	NH ₄ Cl	53.49	0.05349	0.000005

Iron and NTA Mix (1,000x)

A	B	C	D	E
Compound	Chemical Formula	FW	g/100mL	Final M (nM)
iron (II) sulfate heptahydrate	FeSO ₄ x 7H ₂ O	278.01	0.0028	101
nitrilotriacetic acid disodium salt	NTA NA ₂ salt	235.1	0.0081	345

Stock Preparation

1 Prepare the following stocks prior to medium assembly.

- 1.1** Iron mix (1,000x)
Trace Metals (100,000x)
Vitamins (100,000x)
AA mix (500x)
Misc. mix (20,000x)
Phosphate mix (1000x)
Inorganic Nitrogen mix (2,000x)
Fatty Acid mix (2,000,000x)

1.2 Filter sterilize stocks using 0.2 µm PES filter.

Base Salts

- 2** Using an acid-washed and autoclaved one Liter screw-top Pyrex bottle, combine basic salts in 997.4095 mL MilliQ-filtered water. (stock volumes subtracted)

2.1 Mix continuously with stir bar or invert the tightly capped bottle.

Each chemical should be fully dissolved before next addition.

Stock Addition

3 Stocks should be added in a Biosafety cabinet or hood with appropriate ventilation.

Swirl or invert media between stock additions.

3.1 Add (per Liter):
10 μ L Trace Metal Mix
10 μ L Vitamin Mix
20 μ L Amino Acid Mix
50 μ L Misc. Mix
500 μ L Inorganic Nitrogen Mix
0.5 μ L Fatty Acid Mix
1mL Iron Mix
1mL Phosphate Mix

4 When in the biosafety cabinet or laminar fume hood filter-sterilize medium using a **0.2 μ m filter** into another sterilized one Liter Pyrex bottle.

5 pH and record.

Media can be stored at room temp (~25 C°) in the tightly closed Pyrex bottle.