



JAN 29, 2024

🌐 Assessment of autophagic flux assay in the adult *Drosophila* brain

Mel Feany^{1,2}

¹Department of Pathology, Brigham and Women's Hospital, Harvard Medical School, Boston, Massachusetts 02115;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD 20815

ASAP Collaborative Research Network

OPEN  ACCESS



Beatrice Weykopf

ABSTRACT

Transgenic flies expressing GFP-mCherry-Atg8a (obtained from the Bloomington *Drosophila* Stock Center) under the control of the neuronal nSyb-GAL4 driver are used for assaying autophagic flux.



DOI:

dx.doi.org/10.17504/protocols.io.bp2l6x661lqe/v1

Protocol Citation: Mel Feany 2024. Assessment of autophagic flux assay in the adult *Drosophila* brain. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bp2l6x661lqe/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Jan 29, 2024

Last Modified: Jan 29, 2024

PROTOCOL integer ID: 94348

Keywords: ASAPCRN

- 1 Transgenic flies expressing GFP-mCherry-Atg8a (obtained from the Bloomington *Drosophila* Stock Center) under the control of the neuronal nSyb-GAL4 driver are used for assaying autophagic flux
- 2 Brains from 10-day-old animals are dissected in ice-cold 1X PBS
- 3 The dissected brain are transferred to a glass slide
- 4 The remaining 1X PBS is removed from the slide
- 5 A drop of DAPI Fluoromount is added and the brain mounted with a coverslip
- 6 The dissected brain is imaged immediately within 5 minutes using a Zeiss confocal microscope under a 63X objective lens, and a line scan performed
- 7 For quantification, the intensity of GFP and mCherry is measured per puncta using ImageJ, and the ratio between GFP and mCherry calculated

