



Dec 16, 2020

HuBMAP | Formalin Fixation and Paraffin Embedding of Tissue Samples

Kelley Knizner¹, Christopher Simmons¹¹University of Pittsburgh Medical Center**2** Works for me dx.doi.org/10.17504/protocols.io.bqp6mvre

Human BioMolecular Atlas Program (HuBMAP) Method Development Community GE Research

Liz McDonough
GE Research

ABSTRACT

This method details formalin fixation and paraffin embedding of the HuBMAP tissue specimens.

Protocol adopted from *Marda Jorgensen, Jerelyn Nick (02/24/2020). HuBMAP: Paraffin Embedding Tissue Samples*. <https://dx.doi.org/10.17504/protocols.io.bam9ic96>

DOI

dx.doi.org/10.17504/protocols.io.bqp6mvre

PROTOCOL CITATION

Kelley Knizner, Christopher Simmons 2020. HuBMAP | Formalin Fixation and Paraffin Embedding of Tissue Samples . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bqp6mvre>

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Dec 14, 2020

LAST MODIFIED

Dec 16, 2020

PROTOCOL INTEGER ID

45534

MATERIALS TEXT

Instrumentation

A	B
Tissue Processor	Leica Peloris II
Embedding Center	Leica Embedding Center with Cryo-console (Model EG1160)

Materials/Reagents

A	B
Materials	Catalog Number/Model
Formalin	Fisher Scientific
Ethyl Alcohol	State of Pennsylvania Alcohol Control Board
Xylene	Fisher Scientific
Embedding media/paraffin	Surgipath Formula R (3801470)
4x4 Gauze	Fisher Scientific (1400249)
Embedding molds	Leica Stainless Steel embedding molds (39LC-700-2)
Curved forceps	Fisher Scientific (16-100-110)

ABSTRACT

This method details formalin fixation and paraffin embedding of the HuBMAP tissue specimens.

Protocol adopted from *Marda Jorgensen, Jerelyn Nick (02/24/2020). HuBMAP: Paraffin Embedding Tissue Samples*. <https://dx.doi.org/10.17504/protocols.io.bam9ic96>

BEFORE STARTING

- It is essential to read the complete instruction booklet before starting work.
- Users must be properly trained.
- Equipment and facility must be maintained and in good working order.

1

- During fixation, formalin will penetrate the tissue causing chemical and physical changes that harden and preserve the tissue's morphology and protect it against subsequent processing steps.
- Because melted paraffin wax is hydrophobic, most of the water in a specimen must be removed before it can be infiltrated with wax, so the tissue must go through dehydration in a series of ethanol solutions.
- Xylene is then used to clear the tissue, displacing the ethanol which will be in turn displaced by molten paraffin wax.

TISSUE PROCESSING

- 2 Tissue processing is done using the Leica Peloris II using the following program.

A	B	C
Formalin	30 minutes	Ambient temperature
75% Ethanol	30 minutes	37 degrees C
80% Ethanol	30 minutes	37 degrees C
95% Ethanol	30 minutes	37 degrees C
95% Ethanol	30 minutes	37 degrees C
100% Ethanol	45 minutes	37 degrees C
100% Ethanol	50 minutes	37 degrees C
100% Ethanol	60 minutes	37 degrees C
Xylene	30 minutes	37 degrees C
Xylene	30 minutes	37 degrees C
Xylene	40 minutes	37 degrees C
Xylene	60 minutes	37 degrees C
Paraffin	30 minutes	60 degrees C
Paraffin	60 minutes	60 degrees C
Paraffin	90 minutes	60 degrees C

Tissue processing program the Leica Peloris II

POST-PROCESSING

- 3 Processed cassette should be placed in wax holding chamber of embedding center.
- 4 Pre-heat embedding utensils in the forceps warming station; warming chambers are located on the side of the paraffin dispensing nozzle.
- 5 Turn on the station cold plate module.
- 6 Open the cassette. Be sure to carefully maintain the tissue orientation. Dispose of cassette lid.
- 7 Select an appropriate sized embedding mold.
- 8 Partially fill a mold with molten paraffin by using forceps to tap the dispenser plate on the embedding center's dispensing nozzle.
- 9 Carefully transfer tissue into the center of the mold using heated forceps. The tissue should be placed cut side down.
- 10 Transfer the mold to the cold pad and gently press tissue flat with forceps, allowing the paraffin to solidify in a thin layer.

- 11 Place the lower basket of the emptied cassette on top of the mold and add additional molten paraffin to fill 2/3 of the cassette basket. If there are bubbles, slowly lift and reposition the cassette to remove.
- 12 Transfer the paraffin filled mold to the cold plate module and allow to solidify.