

May 11, 2025

Assay for the Enzymatic Degradation of PET Beads

DOI

dx.doi.org/10.17504/protocols.io.bp2l6xp8klqe/v1

Joao Vitor Molino¹

¹University of California, San Diego



Barbara Saucedo

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bp2l6xp8klqe/v1

Protocol Citation: Joao Vitor Molino 2025. Assay for the Enzymatic Degradation of PET Beads. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bp2l6xp8klqe/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 20, 2023

Last Modified: May 11, 2025

Protocol Integer ID: 91224

Funders Acknowledgements:

US Department of Energy

Grant ID: DE-EE0009671



Disclaimer

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Abstract

This assay focuses on assessing the enzymatic degradation of PET beads through absorbance measurements at 240nm. The primary goal is to quantify the enzymatic activity on PET beads by monitoring the release of BHET. This is achieved using a UV-transparent microplate, with absorbance readings specifically taken at 240nm.

Guidelines

Ensure the supernatant transfer is consistent in volume across all samples.

The enzyme and PET bead concentrations might need optimization based on the sensitivity of the assay and the activity of the enzyme.

A negative control with an inactivated enzyme can help determine the baseline of non-enzymatic degradation.

Adjustments to the protocol may be necessary based on specific experimental parameters, enzyme characteristics, or the sensitivity required for the assay.



Materials

PET plastic beads

1M Potassium phosphate buffer (pH 8)

Plastic-degrading enzyme solution

Thermocycler with a lid temperature set to 105°C

TECAN plate reader

Greiner Bio-One UV-Star 96-well Microplates (UV-transparent)

PCR tubes






Pipettes and pipette tips

Safety warnings



- ⚠ Follow standard laboratory safety protocols, especially when working with enzymes and potentially sharp plastic materials.
Wear appropriate personal protective equipment throughout the experiment.




Reaction Setup

- 1 In each PCR tube, add 30mg of PET plastic beads.
- 2 Add  100 μL of  1 Molarity (M) potassium phosphate buffer to the tubes
- 3 Add  100 μL of the enzyme solution to the tubes with PET beads to initiate the reaction. These constitute the test samples.
- 4 Prepare a blank control by adding  100 μL of buffer and  100 μL of enzyme solution into a PCR tube without PET beads.

Incubation

- 5 Place the PCR tubes in the thermocycler.
- 6 Set the thermocycler to  68 $^{\circ}\text{C}$ with the lid heated to  105 $^{\circ}\text{C}$ to prevent condensation. Incubate for 12 hours.

Sampling

- 7 Following the incubation period, briefly cool the PCR tubes to room temperature.
- 8 Centrifuge the tubes if necessary to ensure that the PET beads settle to the bottom.
- 9 Carefully transfer  100 μL of the supernatant from each tube into the corresponding wells of the UV-Star™ 96-well Microplate.

Absorbance measurement



- 10 Measure the absorbance of each well at 240 nm using the TECAN plate reader. Ensure to measure both the test samples and the blank control.

Data Analysis

- 11 Calculate the enzymatic activity by comparing the absorbance readings of the test samples with that of the blank control. The difference in absorbance is indicative of BHET release and, consequently, of PET degradation.