



Oct 17, 2021

Hybridization and Immobilization

Chia-Hsien Shih¹¹Chung Shan Medical University

1

dx.doi.org/10.17504/protocols.io.by5gpy3w

CSMU_Taiwan

Chia-Hsien Shih
Chung Shan Medical University

This protocol is to assemble GotCha with circular probe, immobilization probe and magnetic beads.

DOI

dx.doi.org/10.17504/protocols.io.by5gpy3w

Chia-Hsien Shih 2021. Hybridization and Immobilization. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.by5gpy3w>



_____ protocol ,

Oct 17, 2021

Oct 17, 2021









54152

Preparation

- 1 Dilute immobilization probe into 10 μ M





Protocol of Hybridization






- 2 Add  4.5 μ L of 1 μ M circular probe

- 3 Add  **9 µL** of 10X  **7.4** PBS buffer
- 4 Add  **0.9 µL** of 10µM immobilization probe
- 5 Add  **75.6 µL** RNase-free water
- 6 Vortex slightly and spin down
- 7 Heated to  **90 °C** for  **00:00:30** then cool down to  **Room temperature** and keep for ^{30m 30s}
 **00:30:00**

Protocol of Immobilization

36m

- 8 Take  **22.5 µL** of 10mg/ml magnetic beads into a eppendorf
- 9 Vortex then spin down in Centrifuge  **15000 rpm, 00:02:00** 2m
- 10 Carefully take the eppendorf onto the DynaMag and wait for  **00:01:00** to collect the beads 1m
- 11 Wash beads with 1X PBS for three times. Make sure that eppendorf should put on DynaMag when removing supernatant.
- 12 Add  **90 µL** of hybridization solution into the eppendorf with dry beads

- 13 Incubate at  **Room temperature** for  **00:30:00** , then spindown in Centrifuge 32m
 **15000 rpm, 00:02:00**
- 14 Carefully take the eppendorf onto the DynaMag and wait for  **00:01:00** to collect the beads 1m
- 15 Wash beads with 1X PBS for three times. Make sure that eppendorf should put on DynaMag when removing supernatant.
- 16 Preserve the functional beads in  **90 µL** 1X PBS