



Mar 01, 2022

Mouse Perfusion

Haley Geertsma¹¹University of Ottawa

1

dx.doi.org/10.17504/protocols.io.b5swq6fe **Haley Geertsma**
University of Ottawa

This protocol is used to perfuse mice and isolate their brain for cryosectioning.

DOI

dx.doi.org/10.17504/protocols.io.b5swq6fe

Haley Geertsma 2022. Mouse Perfusion. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b5swq6fe>



_____ protocol ,

Mar 01, 2022

Mar 01, 2022

58934

- 1 Sedate mice with 120mg/kg Euthanyl then open their chest cavity to expose their heart. ^{2m}
- 2 Insert the needle tip into the left ventricle and cut the right atria. Begin slowly injecting 10mL ^{5m} 1X phosphate buffered saline (PBS) over 5 minutes.
- 3 Once all of the 1X PBS has been perfused, begin injecting 10mL 4% paraformaldehyde or 10% ^{5m} formalin over 5 minutes.
- 4 Extract brain and incubate them in either 4% paraformaldehyde or 10% formalin for 72 hours ^{3d} at

4°C.

- | | | |
|---|---|-----|
| 5 | Incubate brain in 10% sucrose for 24 hours at 4°C. | 1d |
| 6 | Incubate brain in 20% sucrose for 24 hours at 4°C. | 1d |
| 7 | Incubate brain in 30% sucrose for 24 hours at 4°C. | 1d |
| 8 | Decant brain from sucrose solution and flash freeze in isopentane at -35 to -45°C for 1 minute. | 1m |
| 9 | Cryo-section brain at 40µm. | 30m |