



May 09, 2022

NPC to Astrocyte Differentiation

Celeste M M. Karch¹, Jacob Marsh², Rj Martinez²

¹Washington University in St Louis; ²Washington University in Saint Louis - WUSTL (MO)

1



dx.doi.org/10.17504/protocols.io.q26g781qqlwz/v1

Neurodegeneration Method Development Community
Tech. support email: ndcn-help@chanzuckerberg.com



Celeste M. Karch
Washington University in St Louis

Here, we provide a detailed protocol for differentiation of human induced pluripotent stem cell derived neural progenitor cells into astrocytes.

DOI

dx.doi.org/10.17504/protocols.io.q26g781qqlwz/v1

Celeste M M. Karch, Jacob Marsh, Rj Martinez 2022. NPC to Astrocyte Differentiation . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.q26g781qqlwz/v1>



protocol ,

Jan 07, 2021

May 09, 2022

Jan 07, 2021



Jacob Marsh

Washington University in Saint Louis - WUSTL (MO)

Jan 07, 2021



Celeste M. Karch

Washington University in St Louis

46079

Preparation

- 1 Prepare astrocyte medium (ScienCell Astrocyte Medium #1801) by adding FBS, antibiotic, and astrocyte growth supplement to basal media. All components are provided in ScienCell Astrocyte Medium #1801 kit.

Note: Fully prepared media should be used within two weeks when stored at 4°C , but can be aliquoted and stored at -20°C for several months. Alternatively, aliquot the supplements and prepare smaller volumes of complete media, if desired.

- 2 Prior to starting, NPCs are cultured in one well of a six well plate until they reach 70% confluency

NPC to Astrocyte Differentiation - Day 1

3m

- 3 Prepare wells/plates/dishes to be used for the differentiation process with Matrigel

For reference, we coat a single well of a six well plate with 1 mL of Matrigel

- 4 Remove media from NPCs and replace with 3 mL Accutase

- 5 Incubate at 37°C for $00:03:00$

3m

- 6 Tap the plate to dissociate the cells and add 6 mL of room temperature DMEM/F12 and collect cells into a 15 mL tube

- 7 Count the number of cells and move the required volume to a new 15 mL tube. Use the table below as a guide for the number of cells to plate into specific vessels.

Plate Size	# of Cells	Volume of Media
6-well plate	60,000/well	2mL/well
24-well plate	15,000/well	500uL/well
10cm dish	360,000	10mL

8 Centrifuge the required volume at 750-800 rpm for 5 minutes

9 Aspirate the supernatant and resuspend the cell pellet in NPC medium

NPC to Astrocyte Differentiation - Day 2

3m

10 Remove NPC media and replace with equivalent volume of complete astrocyte medium

NPC to Astrocyte Differentiation - Day 4 and Onward

11 Change media every 2-3 days for 30 days and passage with Accutase, centrifuge at 750-800 rpm for 10 minutes and re-plate on a Matrigel coated plate whenever culture becomes 100% confluent

12 Within the first 30 days, keep the number of cells per well low to aid in differentiation

13 After 30 days, cells should be >98% positive for S100b. Some lines may express high levels of GFAP, others will not.

After full differentiation, cell confluency/number at passage is no longer important.

Astrocyte Freezing

14 Dissociate cells with Accutase and collect in complete astrocyte medium in a 15mL conical tube

15 Spin cells at 750-800 rpm for 10 minutes

The slower spinning speed helps to reduce the percentage of cell death

16 Aspirate media from pellet and resuspend in freezing media (50% complete astrocyte media,

40% FBS, 10% DMSO)

- 17 Transfer cell suspension to cryovials and place in -80C for 24-48 hours before transferring to liquid nitrogen for long-term storage

Astrocyte Thawing

- 18 Prepare wells/plates/dishes to be used for the differentiation process with Matrigel

For reference, we coat a single well of a six well plate with  1 mL of Matrigel

- 19 Add 9mL of complete astrocyte medium to a 15mL conical tube
- 20 Thaw frozen vial by swirling in water-bath at 37C for 60-90 seconds
- 21 Add thawed cells to 15mL conical tube containing astrocyte medium
- 22 Spin conical tube at 750-800 rpm for 10 minutes to pellet cells
- 23 Aspirate media from cell pellet and replace with 2mL of complete astrocyte medium
- 24 Transfer cell suspension to Matrigel coated plate and incubate at 37C