



2 ▼

Mar 17, 2022

# Catalase test for bacterial identification

## V.2

lydiariver<sup>1</sup><sup>1</sup>Universidad Autónoma de Nuevo León

1

[dx.doi.org/10.17504/protocols.io.b4yjqxun](https://dx.doi.org/10.17504/protocols.io.b4yjqxun)

lydiariver

The catalase enzyme is found in most aerobic and facultative anaerobic bacteria that contain cytochrome, being the main exception *Streptococcus*. Organisms that do not have the cytochrome system also lack the enzyme catalase and therefore cannot break down hydrogen peroxide.

DOI

[dx.doi.org/10.17504/protocols.io.b4yjqxun](https://dx.doi.org/10.17504/protocols.io.b4yjqxun)<https://doi.org/10.1371/journal.pone.0264273>

lydiariver 2022. Catalase test for bacterial identification. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.b4yjqxun>  
lydiariver



protocol

Palacios-Saucedo GdC, Rivera-Morales LG, Vázquez-Guillén JM, Caballero-Trejo A, Mellado-García MC, Flores-Flores AS, González-Navarro JA, Herrera-Rivera CG, Osuna-Rosales LE, Hernández-González JA, Vázquez-Juárez R, Barrón-Enríquez C, Valladares-Trujillo R, Treviño-Baez JD, Alonso-Téllez CA, Ramírez-Calvillo LD, Cerda-Flores RM, Ortiz-López R, Rivera-Alvarado M

protocol ,

Feb 11, 2022

Mar 17, 2022

58091

- Hydrogen peroxide (Degasa, Cd. Mexico, Mexico).
- Glass slides.
- Plastic loop
- Group B *Streptococcus agalactiae* culture on blood agar plate.

- 1 Transfer 2-4 colonies of the bacterial growth in blood agar plates to a glass slide using a sterile plastic loop, make circles, and let it dry.

NOTE: Culture should be 18 to 24 hours old. Metal loops could cause false positives.

- 2 Place a drop of hydrogen peroxide ( $H_2O_2$ ) on the glass slide with an eyedropper.

NOTE: Store the hydrogen peroxide in a dark bottle and avoid exposure to light. Keep refrigerated when not in use.

- 3 Observe immediate results.  
Positive: The oxygen released will be observed as a formation of bubbles.  
Negative: No or very few bubbles produced.

NOTE: Group B *Streptococcus agalactiae* is catalase-negative.




Figure 2. *Staphylococcus aureus* catalase positive (left). *Streptococcus agalactiae* catalase negative (right).