



Aug 23, 2020

© ELISA for measurement of monocyte chemoattractant protein-1 (MCP-1/CCL2) in human serum.

Angel A Justiz-Vaillant¹, Belkis Ferrer-Cosme²

¹University of the West Indies St. Augustine; ²"Saturnino Lora Torres" Provincial Teaching Clinical Surgical Hospital. Cuba

In Development dx.doi.org/10.17504/protocols.io.bj7fkrjn

University of the West Indies angel.vaillant@sta.uwi.edu

Angel Justiz-Vaillant University of the West Indies St. Augustine

ABSTRACT

MCP-1 is believed to play an important role in monocyte infiltration into tumor tissues. [1] MCP-1 plays an important role in the pathogenesis of atherosclerosis. There is considerable evidence that supports that monocytes containing MCPs and macrophages influence the growth of other cell types within the atherosclerotic lesion. [2]

References

- 1. Deshmane SL, Kremlev S, Amini S, Sawaya BE. Monocyte chemoattractant protein-1 (MCP-1): an overview. *J Interferon Cytokine Res.* 2009;29(6):313-326. doi:10.1089/jir.2008.0027
- Lin J, Kakkar V, Lu X. Impact of MCP-1 in atherosclerosis. Curr Pharm Des. 2014;20(28):4580-4588. doi:10.2174/1381612820666140522115801

DOI

dx.doi.org/10.17504/protocols.io.bj7fkrjn

PROTOCOL CITATION

Angel A Justiz-Vaillant, Belkis Ferrer-Cosme 2020. ELISA for measurement of monocyte chemoattractant protein-1 (MCP-1/CCL2) in human serum.. **protocols.io** https://dx.doi.org/10.17504/protocols.io.bj7fkrjn

LICENSE

G

This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 23, 2020

LAST MODIFIED

Aug 23, 2020

PROTOCOL INTEGER ID

40903

- 1 An anti-human MCP-1/CCL2 coating antibody is adsorbed onto the microwells by incubation overnight at 4°C with carbonate-bicarbonate buffer.
- 2 Add 50 μl of human serum or plasma. Human MCP-1/CCL2 present in the serum or plasma binds to antibodies

Citation: Angel A Justiz-Vaillant, Belkis Ferrer-Cosme (08/23/2020). ELISA for measurement of monocyte chemoattractant protein-1 (MCP-1/CCL2) in human serum.. https://dx.doi.org/10.17504/protocols.io.bj7fkrjn

	adsorbed into the microwells.
3	The microplate is blocked with 3% non-fat milk-PBS buffer and later wash to remove unbound proteins.
4	Fifty (50) μ l of biotin-conjugated anti-MCP-1/CCL2 antibody is added. The optimal dilution must be investigated.
5	The microplate is rewashed with PBS-Tween 20 buffer, pH 7.4.
6	One hundred μ I of streptavidin-HRP conjugate is added and it binds to the biotin-conjugated anti-MCP-1/CCL2 antibody. The optimal dilution of this conjugate must be investigated.
7	The plate is washed following incubation to remove the unbound Streptavidin-HRP.
8	Add 100 µl of 3',3',5',5'- tetramethylbenzidine (TMB; Sigma-Aldrich) into each well.
9	Incubate the microwells in the dark for 20 min.
10	A colored product is formed in proportion to the quantity of MCP-1/CCL2 present in the sample or standard.
11	The reaction is terminated by addition of 100 μ l 3M H2SO4 and the absorbance is measured at 450 nm.
12	A standard curve is made from 7 human MCP-1/CCL2 standard dilutions and the human MCP-1/CCL2 sample concentration is determined.
13	For better results place the microplate on a microplate shaker in every incubation.