

Aug 21, 2020

3: 30mer branch melting temperatures (SABER-FISH)

In 1 collection

Jocelyn Y. Kishi^{1,2,3}, Sylvain W. Lapan^{4,3}, Brian J Beliveau^{1,2,5,3,6}, Emma R. West^{4,3}, Allen Zhu^{1,2}, Hiroshi M. Sasaki^{1,2}, Sinem Saka^{1,2}, Yu Wang^{1,2}, Constance L Cepko^{4,7,6}, Peng Yin^{1,2,6}

¹Wyss Institute for Biologically Inspired Engineering, Harvard University, Boston, MA, USA;

²Department of Systems Biology, Harvard Medical School, Boston, MA, USA; ³These authors contributed equally;

⁴Department of Genetics, Blavatnik Institute, Harvard Medical School, Boston, MA, USA;

⁵Present address: Department of Genome Sciences, University of Washington, Seattle, WA, USA;

⁶Correspondence: py@hms.harvard.edu (P.Y.), cepko@genetics.med.harvard.edu (C.L.C.), beliveau@uw.edu (B. J. B.);

⁷Howard Hughes Medical Institute, Chevy Chase, MD, USA

1 Works for me This protocol is published without a DOI.

Human Cell Atlas Method Development Community

Jocelynkishi

ABSTRACT

This protocol describes the design of 30mer branch sequences.

This protocol is part of the [SABER-FISH collection](#).

EXTERNAL LINK

<http://saber.fish/>

THIS PROTOCOL ACCOMPANIES THE FOLLOWING PUBLICATION

Kishi, J.Y., Lapan, S.W., Beliveau, B.J. et al. SABER amplifies FISH: enhanced multiplexed imaging of RNA and DNA in cells and tissues. Nat Methods 16, 533–544 (2019). <https://doi.org/10.1038/s41592-019-0404-0>

ATTACHMENTS

SABER amplifies
FISH_enhanced
multiplexed imaging of
RNA and DNA in cells and
tissues.pdf

PROTOCOL CITATION

Jocelyn Y. Kishi, Sylvain W. Lapan, Brian J Beliveau, Emma R. West, Allen Zhu, Hiroshi M. Sasaki, Sinem Saka, Yu Wang, Constance L Cepko, Peng Yin 2020. 3: 30mer branch melting temperatures (SABER-FISH). **protocols.io**
<https://protocols.io/view/3-30mer-branch-melting-temperatures-saber-fish-bh9hj936>

MANUSCRIPT CITATION please remember to cite the following publication along with this protocol

Kishi, J.Y., Lapan, S.W., Beliveau, B.J. et al. SABER amplifies FISH: enhanced multiplexed imaging of RNA and DNA in cells and tissues. Nat Methods 16, 533–544 (2019). <https://doi.org/10.1038/s41592-019-0404-0>

EXTERNAL LINK

<http://saber.fish/>

COLLECTIONS ⓘ




SABER-FISH – Signal amplification for multiplexed fluorescence in situ hybridization assays

KEYWORDS

branch melting temperature, temperature

LICENSE

 This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 06, 2020

LAST MODIFIED

Aug 21, 2020

OWNERSHIP HISTORY

Jul 06, 2020



Julia Rossmannith protocols.io

Aug 19, 2020



Jocelynkishi

PROTOCOL INTEGER ID

38921

PARENT PROTOCOLS

Part of collection

[SABER-FISH – Signal amplification for multiplexed fluorescence in situ hybridization assays](#)

SAFETY WARNINGS

For hazard information and safety warnings, please refer to the SDS (Safety Data Sheet).

- 1 It's recommended to use a temperature **at least 1 degree lower than the lowest melting temperature** of all branch sequences you plan to use (see plot of melting temperatures Fig. S2).

You can find these melting temperature curves reported for each sequence, as well as those computed for 20mer imagers, 42mer barcode sequences, and an example set of FISH probes, reported in

 **Supplementary Table 2.xlsx** (available on the [Nature Methods website](#)).




Figure S2: Melting temperatures of branches under different formamide conditions. Melting temperatures of 30mer branch binding sequences are shown for the 50 designed PER primers in 1xPBS with different concentrations of formamide. Modeled with Biopython.⁶⁹