



Sep 08, 2020

SPOT1 assay

Guanhua Xun¹, Huimin Zhao¹, stlane2¹¹University of Illinois at Urbana-Champaign
1 Works for me dx.doi.org/10.17504/protocols.io.bk25kyg6

SPOT

stlane2

DOI

dx.doi.org/10.17504/protocols.io.bk25kyg6

PROTOCOL CITATION

Guanhua Xun, Huimin Zhao, stlane2 2020. SPOT1 assay. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bk25kyg6>

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Sep 08, 2020

LAST MODIFIED

Sep 08, 2020

PROTOCOL INTEGER ID

41789

STEPS MATERIALS


NAME	CATALOG #	VENDOR
Nuclease-free Water - 25 ml	B1500S	New England Biolabs
Bst 2.0 WarmStart DNA Polymerase - 8,000 units	M0538L	New England Biolabs
Deoxynucleotide Solution Mix - 8 umol of each	N0447S	New England Biolabs
Non-CRISPR nuclease	N/A	
Non-primer oligos	N/A	
E gene primer mix	N/A	
Isothermal Amplification Buffer - 6.0 ml	B0537S	New England Biolabs
N gene primer mix	N/A	
Reporter Probe 2	N/A	
WarmStart RTx Reverse Transcriptase - 250 rxns	M0380L	New England Biolabs
Manganese(II) chloride tetrahydrate	M3634	Sigma Aldrich
Magnesium Sulfate (MgSO4) Solution - 6.0 ml	B1003S	New England Biolabs
Reporter probe 1	N/A	
Saliva sample	N/A	
QuickExtract DNA Extraction Solution	QE09050	Lucigen

EQUIPMENT

NAME	CATALOG #	VENDOR
------	-----------	--------

NAME	CATALOG #	VENDOR
SPOT1 device	N/A	
SPOT1 Device	N/A	


- Using the first provided microcap, a sample is collected into capillary A, containing QuickExtract DNA Extraction Solution (Lucigen). Insert the capillary into the SPOT1 device and press the "Start" button to run the 5-minute pretreatment.



QuickExtract DNA Extraction Solution

by Lucigen

Catalog #: QE09050



SPOT1 Device


Incubating fluorometer

University of Illinois N/A

🔥 95 °C

🕒 00:05:00

- After pretreatment, remove capillary A from the SPOT1 device and use the second provided microcap to transfer a small volume of pretreated sample to capillary B, which contains the SPOT assay mastermix. Insert capillary B into the SPOT1 device and press the "Start" button to initiate the 35-minute detection reaction.



SPOT1 device

Incubating fluorometer

University of Illinois N/A

🔥 63 °C 🕒 00:30:00

🔥 98 °C 🕒 00:05:00

SPOT assay mastermix:

	Initial concentration	Final concentration	Amount (μL)
Upper compartment			
WarmStart® Bst 2.0	8000 units/mL	160 units/mL	2
WarmStart® RTx	15,000 units/mL	150 units/mL	1
Isothermal amplification buffer	10X	0.5X	8
dNTPs	10 mM	0.7 mM	5.6
MgSO4	100 mM	4 mM	3.2
N gene primer mix	10X	0.25X	2

E gene primer mix	10X	0.25X	2
Saliva samples			5
Non-CRISPR nuclease	5 mg/mL or 55 μ M	1.375 μ M	2
MnCl ₂	50 mM	0.5 mM	0.8
Non-primer oligos (total 6 oligos)	100 μ M	625 nM	3
Reporter probe 1	100 μ M	156.25 nM	0.125
Reporter probe 2	100 μ M	312.5 μ M	0.25
Nuclease-free water			44.025
Total			80



Bst 2.0 WarmStart DNA Polymerase - 8,000 units
by New England Biolabs
Catalog #: [M0538L](#)

2 μ l



WarmStart RTx Reverse Transcriptase - 250 rxns
by New England Biolabs
Catalog #: [M0380L](#)

1 μ l



Isothermal Amplification Buffer - 6.0 ml
by New England Biolabs
Catalog #: [B0537S](#)

8 μ l



Deoxynucleotide Solution Mix - 8 μ mol of each
by New England Biolabs
Catalog #: [N0447S](#)

5.6 μ l [**M**] **0.7 Milimolar (mM)**



Magnesium Sulfate (MgSO₄) Solution - 6.0 ml
by New England Biolabs
Catalog #: [B1003S](#)

3.2 μ l [**M**] **4 Milimolar (mM)**



N gene primer mix

Catalog #: N/A



E gene primer mix

Catalog #: N/A



Saliva sample

Catalog #: N/A



Nuclease-free Water - 25 ml

by New England Biolabs

Catalog #: B1500S



Non-CRISPR nuclease

Catalog #: N/A



Manganese(II) chloride tetrahydrate

by Sigma Aldrich

Catalog #: M3634

 0.8 μ l [M] 0.5 Milimolar (mM)



Non-primer oligos

Catalog #: N/A



Reporter probe 1

Catalog #: N/A



Reporter Probe 2

Catalog #: N/A

- 3 Result ("Positive"/"Negative"/"Inconclusive") will be displayed on SPOT1 device LCD screen after completion of detection reaction and the 1-minute cooling period.