

Oct 03, 2024

## Measuring root system stiffness in maize or sorghum

DOI

[dx.doi.org/10.17504/protocols.io.yxmvmer9ng3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmer9ng3p/v1)

Ashley N. Hostetler<sup>1</sup>, Jonathan W. Reneau<sup>1</sup>, Erin E. Sparks<sup>1</sup>

<sup>1</sup>Department of Plant and Soil Sciences and the Delaware Biotechnology Institute, University of Delaware, Newark, DE, USA

 Ashley N Hostetler

University of Delaware

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.yxmvmer9ng3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmer9ng3p/v1)

**Protocol Citation:** Ashley N. Hostetler, Jonathan W. Reneau, Erin E. Sparks 2024. Measuring root system stiffness in maize or sorghum. [protocols.io https://dx.doi.org/10.17504/protocols.io.yxmvmer9ng3p/v1](https://dx.doi.org/10.17504/protocols.io.yxmvmer9ng3p/v1)

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 01, 2024

**Last Modified:** October 03, 2024

**Protocol Integer ID:** 109060

**Keywords:** Sorghum seed, DC motor, Flex Test, Sensors, Displacement rate

**Funders Acknowledgement:**

United States Department of Agriculture (USDA) National Institute of Food and Agriculture (NIFA)  
Postdoctoral Fellowship

Grant ID: 2022-67012-36840

## Abstract

Plant mechanical failure has been evaluated with proxy measures that quantify root failure strength. However, understanding the force-displacement rate to reach that failure strength is currently a gap in the field. The force-displacement rate is called the root system stiffness. Here, we outline the steps necessary measuring root system stiffness in maize or sorghum using a new tool developed by our lab called SMURF (Sorghum & Maize Under Rotational Force).

## Guidelines

- **Plant Material:**

Plant maize or sorghum seed under standard planting procedures. No specific field conditions are required. Maintain research fields.

### Note

It is important to remove weeds and debris around the base of plant to limit non-plant forces being recorded.

### Note

It is suggested that plants are not tested before the V8 stage.

- **Device:**

A device, called SMURF (**S**orghum and **M**aize **U**nder **R**otational **F**orce), was designed to non-destructively measure the root system stiffness of large grain root systems (Figure 1A). The SMURF applies rotational displacements and measures a force. Below we outline the components of the SMURF as a general overview of the device.

1. The SMURF consists of a physical body, electronics, powertrain, sensors, and software (**Figure 1**).

- The physical body consists of (**Figure 1A**):

1. A 3D-printed case housing the electronics and powertrain (**Figure 1A-1**)
2. A battery pack (**Figure 1A-2**)
3. A set of arms with straps to attach the device to the plant (**Figure 1A-3**)
4. A sensor housing case for the load cell (**Figure 1A-4**)
5. A rotational foot for device grounding (**Figure 1A-5**)


- The internal electronics are run via an Arduino microcontroller connecting a motor driver, load cell amplifier, and Bluetooth module (**Figure 1B**).

- The powertrain consists of:

1. A brushed DC motor
2. A lead screw
3. A battery

- A 10 kg load cell (sensor)

- A custom-built iPad application controls the device (**Video 1**)



**Figure 1. Overview of SMURF device.** (A) The physical body is composed of (1) a 3D-printed case housing the electronics and powertrain, (2) a battery pack, (3) a set of arms with straps to attach the device to the plant and limit stalk bending, (4) a sensor housing case for the 10 kg load cell, and (5) a rotational foot for device grounding. (B) The internal electronics, located within the 3D-printed case, are run via an Arduino microcontroller that connects a motor driver, load cell amplifier, and Bluetooth module. Dotted lines indicate a wireless connection, whereas solid lines indicate a hardwired connection between components.



**Video 1:** Custom-built iPad application controls the device

## Before Operating the SMURF:

- 1 Open the Device Controller iPad application (Video 1).
- 2 Hold the SMURF horizontally and power on.
- 3 Calibrate the load cell in the absence of a load.
- 4 Enter Research ID (e.g. researcher name/initials/identifier) in iPad application (Video 1).
- 5 Click button: "Save Researcher ID" (Video 1).
- 6 Click button: "Go to Testing Screen" (Video 1).

### Note

Once this is clicked, a pop-up box will prompt connection to the SMURF device via Bluetooth.

- 7 Click button: "Fully Retract Foot". This will re-home the device (Video 1).

## Operating the SMURF:

- 8 Attach the device with the two arms as low as possible on the plant while avoiding brace roots.

### Note

The straps should be tightly secured to prevent slippage.

- 9 Enter Plant ID into the text box (Video 1).

## 10 Select appropriate Test Type (Video 1).

**Note**

Test Type is an arbitrary, user-defined designation of A through F, which adds flexibility for the user to define different experiments within the same testing set.

## 11 Click button: "Get Initial Height" (Video 1)

**Note**

This instructs the device to lower the foot to the ground until the device weight is read on the load cell, which is the "tare" step of the measurement device.

## 12 Click button: "Small Flex Test" or "Large Flex Test" to start the testing cycle (Video 1)

- A Small Flex Test should be used for any small or mutant plants, whereas a Large Flex Test is the default for all healthy plants.

**Note**

The device should not be touched during the testing cycle to limit reading and recording external, non-plant forces. During the testing cycle, the lead screw extends 10.5 mm for a Small Flex Test or 15.8 mm for a Large Flex Test.

## 13 Click button: "Accept Results" or "Reject Results" (Video 1)

- "Accept Results" prompts the data directory to save the data.
- "Reject Results" prompts the device to go back to the original testing screen.

## 14 Click button: "Move to Next Plant," which will fully retract the foot and make the device ready for the next testing cycle (Video 1)

### Note

This process is repeated for each plant.