



Aug 23, 2021

Electrophysiological characterizations of pancreatic islet cells

Toshinori Hoshi Laboratory¹¹Department of Physiology, University of Pennsylvania, Philadelphia, PA 19104

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.bvyyn7xw

Human Islet Research Network

 Lili Liang

ABSTRACT

Hormone release from pancreatic cells, for example glucagon from α -cells and insulin from β -cells, is controlled by the cells' membrane potentials (V_m), which are in turn finely controlled by multiple classes of ion channels, transporters, and pumps. The perforated whole-cell patch-clamp method represents a direct and well-established way to monitor V_m in individual cells and also to investigate ionic currents (I_m) that contribute to determination of V_m . Here a step-by-step guide to obtain various patch-clamp measurements from pancreatic cells in isolated intact human islets is provided. The method requires a standard electrophysiological station and β escin is used as the perforating agent. Select exemplar current-clamp V_m results and voltage-clamp I_m results are provide. The electrophysiological characteristics may be used to distinguish pancreatic cell types, for example α -cells vs. β -cells in real time.

DOI

dx.doi.org/10.17504/protocols.io.bvyyn7xw

EXTERNAL LINK

<https://hpap.pmacs.upenn.edu/assets/workflow/Islet%20Physiology%20Studies/Intracellular%20Calcium%20Flux%20and%20Electrophysiology.pdf>

PROTOCOL CITATION

Toshinori Hoshi Laboratory 2021. Electrophysiological characterizations of pancreatic islet cells.

protocols.io

<https://dx.doi.org/10.17504/protocols.io.bvyyn7xw>

KEYWORDS

null, HPAP, HIRN, Human islet electrophysiology

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 21, 2021

LAST MODIFIED

Aug 23, 2021

PROTOCOL INTEGER ID

50936

Reagents/Buffers/Materials Needed for Experiments

1. A complete electrically quiet electrophysiological setup
 - a. A vibration isolate table
 - b. An inverted (or an upright) microscope
 - c. An inline-perfusion solution heater
 - d. A microscope stage heater
 - e. A low-speed peristaltic pump
 - f. A vacuum source (a house vacuum line or a vacuum pump) and a waste flask
 - g. A micromanipulator
 - h. A patch-clamp amplifier with its head stage assembly mounted on the aforementioned micromanipulator
 - i. An analog-to-digital/digital-to-analog converter (if not integrated into the patch-clamp amplifier)
 - j. A data acquisition program
 - k. A recording chamber with an Ag/AgCl pellet ground wire
2. pipette puller
3. Electrode tip polisher
4. Sylgard 184 elastomer or dental wax to coat electrodes to reduce capacitance
5. Patch-clamp glass
6. A small bath sonicator
7. A dissection microscope
8. A diamond glass cutting pen and a pair of tweezers
9. Recording solutions
 - a. No glucose external solution
130 mM NaCl, 4mM KCl, 2 CaCl₂, 2 MgCl₂, 30 mM mannitol, 10 mM HEPES, pH7.4 with *N-methyl-D*-glucamine at ♢ 35 °C
 - b. 5 mM glucose external solution
130 mM NaCl, 4mM KCl, 2 CaCl₂, 2 MgCl₂, 15 mM mannitol, 10 mM HEPES, pH7.4 with *N-methyl-D*-glucamine at ♢ 35 °C
 - c. 25 mM glucose external solution
130 mM NaCl, 4mM KCl, 2 CaCl₂, 2 MgCl₂, 10 mM HEPES, pH7.4 with *N-methyl-D*-glucamine at ♢ 35 °C
 - d. Sulfate internal solution
76 K₂SO₄ mM, 10 mM KCl, 10 mM NaCl, 6 MgCl₂, 30 mM mannitol, 10 mM HEPES, pH7.4 with *N-methyl-D*-glucamine at ♢ 35 °C

Note:

Based on the Ca²⁺ chelating abilities of sulfate, the free Ca²⁺ concentration is estimated to be in the low μM range and the free Mg²⁺ concentration should be about 2 mM

10. β-escin (8mM in water)

The final concentration is 8 μM; diluted with the internal solution before each recording session and sonicate


Procedure

- 1 Sonicate 8 μM escin diluted with the internal solution in the bath sonicator for > 10 min
4 μM escin is sometimes sufficient
- 2 Start the perfusion of the recording chamber so that the chamber solution is at ♢ 35 °C ; the perfusion is continuous


- 3 Take a culture dish with a coverslip with pancreatic islets (or pancreatic islet cells) using the dissection microscope
- 4 Cut out a small section, with an islet (or cells), using a diamond pen, and transfer the coverslip piece to the recording chamber filled with the desired recording solution (e.g., 5 mM glucose)
- 5 Equilibrate the islet (or cells) in the chamber for > 10 min
- 6 Fill the tip of a polished, sylgard- or wax-coated patch electrode with the internal recording without escin and back fill with the recording with escin
- 7 The input resistance of the electrode should be 3 to 6 Mohms, depending on the types of the data required; lower for voltage-clamp experiments and higher for current-clamp experiments.
- 8 Form a seal with a few Gohms in resistance
- 9 Change the holding voltage to -70 mV (in the whole-cell mode convention) and apply small short square pulses (-10 mV in size from the holding voltage and 20 ms in duration) to monitor the resistance and capacitance
- 10 Within 5 to 10 min, adequate perforated whole-cell access should be achieved
- 11 Apply short small square pulses in the voltage-clamp mode (-10 mV in size from the holding voltage and 20 ms in duration) to measure the input resistance and capacitance
- 12 Compensate the whole-cell capacitance and series resistance
- 13 Measure membrane potential (V_m) in the current-clamp mode; change the external solutions if desired
- 14 13 Measure membrane potential (I_m) in the voltage-clamp mode; change the external solutions if desired

Data analysis


- 15
 1. Export data to IgorPro (Wavemetrics) for data plotting and analysis
 2. Data visualization
 - a. Current-clamp V_m measurement example (from a cell in an intact islet)



b. Ramp voltage-clamp I_m measurement example (from a sorted β cell)



c. Step voltage-clamp I_m measurement example (from a sorted β cell)



3. Other custom analysis

- Igor scripts may be written for custom and automated analysis procedures as needed