




Jun 05, 2024 Version 1

Simplified ATAC-seq protocol on frozen brains of fruit flies V.1



Forked from [Isolation of nuclei from frozen tissue for ATAC-seq and other epigenomic assays](#)



This protocol is a draft, published without a DOI.

Ryan Corces¹, William J. Greenleaf², Howard Y. Chang², Siyuan Feng³, Jamie Freeman⁴

¹Gladstone Institutes; ²Stanford University; ³University of Wisconsin-Madison; ⁴University of Wisconsin Madison

Human Cell Atlas Method ...



Siyuan Feng

University of Wisconsin-Madison

OPEN  ACCESS



External link: <https://www.ncbi.nlm.nih.gov/pubmed/28846090>

Protocol Citation: Ryan Corces, William J. Greenleaf, Howard Y. Chang, Siyuan Feng, Jamie Freeman 2024. Simplified ATAC-seq protocol on frozen brains of fruit flies. **protocols.io** <https://protocols.io/view/simplified-atac-seq-protocol-on-frozen-brains-of-f-daw72fhn>

Manuscript citation:

Corces MR, Trevino AE, Hamilton EG, Greenside PG, Sinnott-Armstrong NA, Vesuna S, Satpathy AT, Rubin AJ, Montine KS, Wu B, Kathiria A, Cho SW, Mumbach MR, Carter AC, Kasowski M, Orloff LA, Risca VI, Kundaje A, Khavari PA, Montine TJ, Greenleaf WJ†, Chang HY†. An improved ATAC-seq protocol reduces background and enables interrogation of frozen tissues. 2017. Nature Methods. (PMID: 28846090).

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: March 19, 2024

Last Modified: June 05, 2024

Protocol Integer ID: 96959

Keywords: ATAC-seq, Nuclei Isolation, Frozen Tissue Homogenization



Abstract

This protocol enables the isolation of nuclei from frozen tissues. These nuclei are suitable for use in ATAC-seq, single-cell ATAC-seq, ChIP-seq, HiC/3C, and many other assays.

Modified for a small amount of fruit fly brains from the omni-ATACseq protocol from Ryan Cortes

Guidelines

The quality of the tissue at the time of freezing is a major factor in the quality of data downstream. If freezing tissue for use later, you may want to consider cryopreserving 50 mg tissue chunks in BAM Banker cryopreservative. If using previously frozen tissue, the ischemic time or post mortem interval should be kept as short as possible.

Nuclei can be cryopreserved in BAM Banker. In our experience, freeze-thaw cycles do not adversely affect data quality as long as thawing is carried out on ice and freezing is performed using a slow-freeze freezing container. We have had success isolating nuclei from every tissue we have tried including difficult to work with tissues such as heart.

Materials

STOCK BUFFER INGREDIENTS

A	B	C	D	E	F	G	H	I	J
Stock (planned Conc.)	Name	chemical info	Formula mass (g/mol)	Volume (ml) of solution (planned) 2024.03.20	Mass (g) of solube (planned) 2024.03.20	Volume (ml) of solution (actual) 2024.03.20	Mass (g) of solube (actual) 2024.03.20	Actual Conc. (M) 2024.03.20	autoclave bottle before usage
1	M Sucrose		342.30	100	34.23	100	34.34	1.003213 5553607946	1
2	M KCl	BDH, potassium chloride, purity 99 - 100.5%, lot No. 118446, opened 2024.03.20	74.47	100	14.894	100	15.37	2.063918 3563851216	1
1	M MgCl2 (hexahydrate)	dot scientific inc., magnesium Chloride hexahydrate, purity unknown, lot No. 61024- 91729, opened 2019.06.18	203.3	100	20.33	100	20.44	1.005410 7230693556	1
0.75	M Tricine- KOH pH 7.8		179.22	200	26.883	150	13.38	0.497712 3088940967	1
2.5	M NaCl	Fisher Scientific, Sodium Chloride, putiry 99.3%, log No. 145475, already opened 2024.03.20	58.44	200	29.22	250	29.22	2	1
60	% Iodixanol								
0.5	M Spermidin		254.63	2	0.25463	2	0.257	0.504653 81141263	

A	B	C	D	E	F	G	H	I	J
	e							8	
0.15	M Spermine		202.34	8	0.242808	8	0.2411	0.148944 84530987 448	
1	M DTT		154.25	3	0.46275	3	0.468	1.011345 21880064 82	
10	% wt./vol. Tween-20	opened 2024.04.0 2	density: 1.095 g/mL at 25 °C	50					
10	% wt./vol. IGEPAL® CA-630	opened 2024.04.0 4	density: 1.06 g/mL at 25 C	50					

STOCK BUFFER INGREDIENTS (dilute by density)

A	B	C	D	E	F	G	H	I	J	K	L	M
Stock (planned Conc.)	Name	chemical info	Density (g/mL at 25 C)	Volume (ml) of solution (planned) 2024.0 4.04	Mass (g) of soluble (planned) 2024.0 4.04	Volume (mL) of soluble (planned) 2024.0 4.04	Volume (mL) of water (planned) 2024.0 4.04	Mass (g) of soluble (actual) 2024.0 4.04	Volume (mL) of soluble (actual) 2024.0 4.04	Volume (mL) of water (actual) 2024.0 4.04	Sterile filter?	
10	% wt./vol. Tween- 20	opened 2024.0 4.02	1.095	50	5	4.5662 100456 62101	45.433 789954 3379	5.005	4.5707 762557 07763	45.429 223744 29224	yes	
10	% wt./vol. IGEPAL ® CA- 630	opened 2024.0 4.04	1.06	50	5	4.7169 811320 75471	45.283 018867 924525	5.005	4.7216 981132 07546	45.278 301886 792455	yes	

STOCK BUFFERS

All stock solutions should be filtered using a 0.22 um PVDF filter system. All solutions except for the 50% Iodixanol solution are stable at 4°C for at least 6 months.

1.034x Homogenization Buffer Stable Solution - 200 ml (2014 ul/sample)

A	B	C	D	E	F
Stock	Name	Final Conc.	Actual stock conc.(made on 2024.03.20)	actual Fold Dilution (x)	actual Total Vol. (ml)
1	M Sucrose	0.26	1.00	3.846153846 153846	52.00000000 000001
2	M KCl	0.03	2.06	68.66666666 666667	2.912621359 2233006
1	M MgCl2	0.01	1.01	101	1.980198019

A	B	C	D	E	F
					8019802
0.75	M Tricine-KOH pH 7.8	0.02	0.50	25	8
-	Water	-	-	-	135.1071806 209747

Diluent Buffer - 100 ml (16.4 ul/sample)

A	B	C	D	E	F
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Actual stock conc.(made on 2024.03.20)</i>	<i>actual Fold Dilution (x)</i>	<i>actual Total Vol. (ul)</i>
2	M KCl	0.15	2.06	13.73333333 3333334	7.281553398 058252
1	M MgCl ₂	0.03	1.01	33.66666666 666667	2.970297029 70297
0.75	M Tricine-KOH, pH 7.8	0.12	0.50	4.166666666 666667	24
-	Water	-	-	-	65.74814957 223879

50% Iodixanol Solution - 50 ml (Remake monthly for stability) (840 ul/sample)

A	B	C	D	E
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Fold Dilution (x)</i>	<i>Total Vol. (ml)</i>
-	Diluent Buffer	1	-	8.33333
60	% Iodixanol	50	1.20	41.66667

ATAC-RSB Buffer - 500 ml (2970 ul/sample)

A	B	C	D	E	F
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Actual stock conc.(made on 2024.03.20)</i>	<i>actual Fold Dilution (x)</i>	<i>actual Total Vol. (ml)</i>
1	M Tris-HCl pH 7.5	0.01	0.50	50	10
5	M NaCl	0.01	2	200	2.5
1	M MgCl ₂	0.003	1.01	336.6666666 666667	1.485148514 8514851
-	Water	-	-	-	486.0148514 8514854

SAME DAY BUFFERS

1x Homogenization Buffer Unstable Solution

****Note – cOmplete Protease Inhibitors come as tablets. It is difficult to use less than 1/2 tablet so we prepare the 1x Homogenization Buffer Unstable Solution in batches of 12 as outlined below.**

A	B	C	D	E	F	G
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Fold Dilution (x)</i>	<i>Vol per 12 samp. (ul)</i>	Vol (ul) made 2024.04.24 (14.4 sample)	Vol (ul) made 2024.05.01 (10 sample, assuming 200 ul per sample)
1.034 1	x HB Stable Solution	1	1.03	24175.00	29010	2417.5
1	M DTT	0.001	1000.00	25.00	30	2.5
500	mM Spermidine	0.5	1000.00	25.00	30	2.5
150	mM Spermine	0.15	1000.00	25.00	30	2.5
10	% NP40	0.3	33.33	750.00	900	75
-	cOmplete Protease Inhibitor mini tablets	-	-	2.50 Tablets	3 tablets	0.25 tablet
	TOTAL VOLUME			25000	30000	2500

30% Iodixanol Solution

A	B	C	D	E	F
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Fold Dilution (x)</i>	<i>Vol per sample (ul)</i>	Vol (ul) made 2024.04.04 (12 samp.)
-	1x Homog. Buffer Unstable	-	-	240.00	2880
50	% Iodixanol Solution	30	1.67	360.00	4320

40% Iodixanol Solution

A	B	C	D	E	F
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Fold Dilution (x)</i>	<i>Vol per sample (ul)</i>	Vol (ul) made 2024.04.04 (12 samp.)
-	1x Homog. Buffer Unstable	-	-	120.00	1440
50	% Iodixanol Solution	40	1.25	480.00	5760

ATAC-RSB-Tween Buffer

A	B	C	D	E	F	G
<i>Stock</i>	<i>Name</i>	<i>Final Conc.</i>	<i>Fold Dilution (x)</i>	<i>Vol per sample (ul)</i>	Vol (ul) made 2024.04.04 (12 samp.)	Vol (ul) made 2024.05.01 (11 samp., assuming

A	B	C	D	E	F	G
						1000 ul per sample)
-	ATAC-RSB	-	-	2970.00	35640	10890
10	% Tween-20	0.1	100.00	30.00	360	110

ATAC-seq Reaction Mix

A	B	C	D
<i>Reagent</i>	<i>Vol per sample (ul)</i>	Vol (ul) made 2024.04.24 (11 samp.)	Vol (ul) made 2024.05.01 (11 samp.)
H2O	5.25	57.75	57.75
PBS	16.5	181.5	181.5
2x TD Buffer	25	275	275
2% Digitonin	0.25	2.75	2.75
10% Tween-20	0.5	5.5	5.5
Tn5 Transposase	2.5	27.5	27.5

PCR master mix (MM)

A	B	C	D
<i>Reagent</i>	<i>Vol per sample 1x50 ul run (ul)</i>	Vol (ul) made 2024.04.25 (10 + 1 samp.)	Vol (ul) made 2024.05.05 (10 + 1 + 1 samp.)
H2O	25.5	280.5	306
2 X KAPA high fidelity buffer	10	110	120
dNTPs	1.5	16.5	18
Hifi enzyme	1	11	12
TOTAL MM VOLUME	38	418	456

Reagents used in this protocol



A	B	C
Item	Supplier	Cat Number
Eppendorf 2 ml Lo-Bind tubes	Sigma	Z666556-250EA
Eppendorf 1.5 ml Lo-Bind tubes	Sigma	Z666548-250EA
Nunc cryovials	Thermo	375418PK
Iodixanol (comes at 60%)	Sigma	D1556-250ML
Digitonin	Promega	G9441
Sucrose	Sigma	S7903-250G
IGEPAL® CA-630 (identical with NP40, which was discontinued)	Sigma	I8896-50ML
Tricine	Sigma	T0377-25G
Potassium Hydroxide (KOH)	Sigma	P5958-250G
cOmplete Protease Inhibitors mini tablets (1/5 size of the full tablets cat. 11697498001)	Roche	11836170001
MgCl ₂	Ambion (Thermo)	AM9530G
KCl	Ambion (Thermo)	AM9640G
DTT	Thermo	R0861
Spermidine	Sigma	S2501
Spermine	Sigma	S3256-1G
70 um Flowmi cell strainers	Fisher	03-421-228
70 um bucket-style cell strainers	BD Falcon	352350
Tris-HCl pH 7.5	Invitrogen	15567-027
NaCl	Ambion (Thermo)	AM9759
Tween 20	Sigma	P2287-100ML
H ₂ O	Invitrogen	10977-015
Dounce Tissue Grinder Set	Sigma	D8938-1SET
INCYTO Disposable hemocytometers	Fisher	22-600-100
BAM Banker	Wako Chemicals	302-14681
RiboLock	Thermo	E00384
0.22 um PVDF Filter Units (500 ml)	Millipore	SCGVU05RE
0.22 um PVDF Filter Units (50 ml)	Millipore	SE1M179M6



	A	B	C
	Tn5 Transposase (TDE1)	Illumina	15027865
	2x TD Tagment DNA Buffer	Illumina	15027866






Before you start the protocol:

- 1 All steps should be performed on ice or at 4 °C . Pre-chill a fixed angle centrifuge to 4°C.
- 2 Pre-chill all tubes. For each sample you are processing, you will need: (i) one 2 ml cryo tube with 10 1.6 mm beads for bead homogenization (ii) one 1.5 ml round-bottom LoBind tube for nuclei pelleting (iii) one 1.5 ml round-bottom LoBind tube for collecting purified DNA
- 3 Prepare all buffers. For faster dissolution, crush protease inhibitor tablets prior to addition to 1x Homogenization Buffer Unstable Solution. DTT, Spermidine, Spermine, and digitonin are stored at -20°C. All other detergents, ATAC-RSB, and other buffers are stored at 4 °C . Do not prepare transposition mix ahead of time.





Isolation of Nuclei via Bead Homogenization and Centrifuging

- 4 move samples from liquid nitrogen, or -80 °C freezer onto On ice (wet ice) wet ice 2 mins before homogenization. Add 200 µL HB to the tube wall, vortex and spin down to get the sample off the tube wall.
- 4.1 if samples were slow-frozen in Bambanker, put the sample on On ice (wet ice) until the Bambanker thaws, then spin down sample 00:05:00 at 4 °C at 750 RCF in a fixed angle centrifuge, and remove supernatant. Add 200 µL HB and 10 1.6 mm beads to the tube. 5m
- 4.2 if samples are fresh, put the sample on on On ice (wet ice) until homogenization (200 µL pre-chilled HB should have been added to the tube before dissection.
- 5 Homogenize in BeadRupter Elite bead mill homogenizer using 2 cycles of the following: 1m 20s
 - Speed 1.0 for 00:00:20
 - Dwell for 00:00:20
 - Speed 2.0 for 00:00:20
 - Dwell for 00:00:20
- 6 Pellet nuclei by spinning 00:05:00 at 4 °C at 350 RCF in a fixed angle centrifuge. 5m




- 7 Re-suspend the nuclei in  1000 μL of ATAC-RSB-Tween Buffer
- 8 Centrifuge nuclei for  00:10:00 at 500 RCF at  4 $^{\circ}\text{C}$ in a fixed angle centrifuge. At this point, the pellet should be clearly visible if 50,000 nuclei were used. Pellets of as few as 10,000 nuclei should be visible.

Transposition of Nuclei

- 9 Using a p1000 pipette, remove all but the last 100 μL of supernatant. Remove last 100 μL with p200 pipette set to 200 μL using a single fluid pipetting motion. Place the tip of your pipette on the opposite side of the tube to where the nuclei pellet is located during this final aspiration step.
- 10 Add 50 μL ATAC-seq Reaction Mix to each tube and pipette up and down 6 times to resuspend nuclei pellet.
- 10.1 Unlike the published ATAC-seq protocols, you do not need to do an individual lysis step in this protocol because the nuclei are exposed to NP40 throughout the Douncing portion of the protocol.
- 11 Incubate reactions at  37 $^{\circ}\text{C}$ for  00:30:00 in a heat block (ideally should use a thermoshaker with 1000 RPM constant shaking).
- 12 After incubation, follow instructions of Qiagen Minelute DNA clean-up Kit to purify the DNA
- 13 Store purified DNA under  -20 $^{\circ}\text{C}$, if not used for PCR for a more than a week. Otherwise the purified DNA can be stored at  4 $^{\circ}\text{C}$ temporarily

DNA Library Amplification

- 14 Remove PCR reagents from freezer to thaw at  Room temperature ; dNTPs, 2X KAPA High fidelity buffer (yellow lid), and primer plate. Once all traces of ice are gone, vortex and spin down. Turn on thermal cycler.
- 15 Remove Hifi enzyme from freezer and place on wet ice. Make PCR master mix (MM)

PCR master mix (MM)

A	B	C	D
1	<i>Reagent</i>	<i>Vol per sample 1x50 ul run (ul)</i>	Vol (ul) made 2024.04.25/04.29 (10 + 1 samp.)
2	H2O	25.5	280.5
3	2 X KAPA high fidelity buffer	10	110
4	dNTPs	1.5	16.5
5	Hifi enzyme	1	11
6	TOTAL MM VOLUME	38	418

16 Mix and spin down master mix, keep cold, and add 38 uL to each reaction tube.

17 Visually inspect primer plate to ensure evaporation has not occurred in the wells you intend to use. Spin down in plate spinner. Carefully open foil to needed wells.

17.1 primer plate in use:

A	B	C	D	E
	2024.04.29			
plate number	3 + 7	H2O	25.5	280.5
well number	1-10	2 X KAPA high fidelity buffer	10	110

18 Add 2 uL of primer mix to each reaction. Make sure to record which index primers are used!

19 Add 10 uL of cleaned tagmentation product. Mix and spin down. Make sure to close PCR tube lids tightly with a marker end.

20 Place tubes in thermal cycler, making sure lid is tightly against your tubes, and run program ML2.

(Optional) Cleanup and Freezing Down of Extra Nuclei



- 21 If you would like to save extra nuclei for other assays or to potentially use in additional ATAC-seq experiments downstream: