

Sep 11, 2024

Paramecium bursaria immunofluorescence protocol

DOI

dx.doi.org/10.17504/protocols.io.yxmvm2ozng3p/v1

David S Milner¹, Estelle Kiliass¹, Thomas A. Richards¹

¹University of Oxford

Richards Lab



David S Milner

University of Oxford

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.yxmvm2ozng3p/v1

Protocol Citation: David S Milner, Estelle Kiliass, Thomas A. Richards 2024. *Paramecium bursaria* immunofluorescence protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.yxmvm2ozng3p/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 01, 2022

Last Modified: September 11, 2024

Protocol Integer ID: 69436

Funders Acknowledgement:

European Research Council

Grant ID: CELL-in-CELL:

819507

Abstract

This protocol has been developed for immunofluorescent localisation of host proteins in *Paramecium bursaria* 186b.

Protocol materials

- ✕ Anti- α -Tubulin antibody Mouse monoclonal **Merck MilliporeSigma (Sigma-Aldrich) Catalog #T6199-25UL** Step 8
- ✕ CITIFLUOR Antifadent Mounting Medium AF2 **Electron Microscopy Sciences Catalog #17971-100** Step 13
- ✕ Phosphate buffered saline **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P4417-50TAB** Before starting
- ✕ Tween 20 **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P1379-500ml** Step 5
- ✕ Hoechst 33342 **Invitrogen - Thermo Fisher Catalog #H1399** Step 11
- ✕ Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) **Abcam Catalog #ab150077** Step 10
- ✕ Goat Anti-Mouse IgG H&L (Alexa Fluor® 488) preadsorbed **Abcam Catalog #ab150117** Step 10
- ✕ Pierce™ 16% Formaldehyde (w/v) Methanol-free **Thermo Fisher Scientific Catalog #28906** Step 3
- ✕ Bovine Serum Albumin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A2153** Step 7
- ✕ Poly-Prep Slides **Merck MilliporeSigma (Sigma-Aldrich) Catalog #P0425-72EA** Step 14

Before start


Prepare 200 mL of PBS

- ✕ Phosphate buffered saline **Sigma Aldrich Catalog #P4417-50TAB**




Fixation




52m

1 Take 6 mL of *P. bursaria* culture per treatment. Centrifuge  800 x g, 00:10:00 and remove 5 mL (i.e. concentrate sample 6-fold). Aliquot 1 mL of concentrated culture into a 1.5 mL micro-centrifuge tube per treatment for downstream steps.

10m

2 Centrifuge  800 x g, 00:10:00 to pellet the cells. Remove supernatant, avoiding the cell pellet.




10m

3 Re-suspend cells in  500 μ L PBS containing 4% formaldehyde [ 375 μ L PBS plus  125 μ L 16% formaldehyde]. Use a fresh ampule of formaldehyde for each experiment.

12m

Incubate for  00:12:00  Room temperature

 Pierce™ 16% Formaldehyde (w/v) Methanol-free **Sigma Aldrich Catalog #28906**

4 Centrifuge  800 x g, 00:10:00, remove supernatant and wash once in  500 μ L PBS. Repeat the centrifugation  800 x g, 00:10:00

20m

Permeabilisation




40m

5 Re-suspend cell pellet in  500 μ L PBS containing  0.05 % volume Tween 20.

10m

Incubate for  00:10:00  Room temperature



 Tween 20 **Sigma Aldrich Catalog #P1379-500ml**

6 Centrifuge  800 x g, 00:10:00, remove supernatant and wash **twice** in  500 μ L PBS. Repeat the centrifugation  800 x g, 00:10:00, x2

30m

Blocking

1h

7 Prepare PBS containing  3 % volume BSA. Filter through a 0.2 micron syringe filter prior to use. Re-suspend cells in  500 μ L of this solution.

1h



01:00:00




Room temperature

Bovine Serum Albumin **Sigma Aldrich Catalog #A2153**

Primary antibody

16h

- 8 Add  1 μL primary antibody directly to the cells in PBS+3% BSA.

16h

A good positive control here is  1 μL Anti- α -Tubulin antibody

Anti- α -Tubulin antibody Mouse monoclonal **Sigma Aldrich Catalog #T6199-25UL**

Overnight




4 °C

Note



You will need to titrate an appropriate concentration of the primary antibody here

Secondary antibody

3h

- 9 After the overnight incubation, wash **three times** in  500 μL PBS.

30m

Repeat the centrifugation  800 x g, 00:10:00 , x3 and re-suspend the cells in  500 μL PBS.

- 10 Add  1 μL secondary antibody.

2h

For rabbit primary antibodies:

Goat Anti-Rabbit IgG H&L (Alexa Fluor® 488) **Sigma Aldrich Catalog #ab150077**

For mouse primary antibodies (e.g. Anti- α -Tubulin antibody used above):

Goat Anti-Mouse IgG H&L (Alexa Fluor® 488) preadsorbed **Sigma Aldrich Catalog #ab150117**



02:00:00



Room temperature in the dark

- 11 If required, add  2 μL of [M] 10 mg/mL Hoechst 33342





5m

Hoechst 33342 **Sigma Aldrich Catalog #H1399**

00:05:00



Room temperature

- 12 Centrifuge  800 x g, 00:10:00 , remove supernatant and wash twice in  500 μL PBS. Repeat the centrifugation  800 x g, 00:10:00 , x2 and re-suspend in  20 μL of residual PBS

20m

Preparing slides

25m

- 13 Add a drop of CITIFLUOR AF2 to the residual PBS.

5m

CITIFLUOR Antifadent Mounting Medium AF2 **Electron Microscopy Sciences Catalog #17971-100**

- 14 Using a wide-bore tip (i.e. a 1 mL pipette tip with the end cut-off), spot approximately

5m

10 μL Poly-Prep Slides **Sigma Aldrich Catalog #P0425-72EA**

Note

Add only one sample to the centre of each slide to allow for sealing

- 15 Add a coverslip and seal with clear nail polish.

15m