




SEP 29, 2023

Manual extraction of High Molecular Weight DNA from single mosquitoes using the QIAGEN MagAttract HMW DNA kit

Fiona Teltscher¹, Mara Harriet Johnson¹, Lawniczak¹

¹Wellcome Sanger Institute



Fiona Teltscher
Wellcome Sanger Institute

OPEN ACCESS



DOI:
dx.doi.org/10.17504/protocols.io.n92ldp6ool5b/v1

Protocol Citation: Fiona Teltscher, Harriet Johnson, Mara Lawniczak 2023. Manual extraction of High Molecular Weight DNA from single mosquitoes using the QIAGEN MagAttract HMW DNA kit. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n92ldp6ool5b/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working
We use this protocol and it's working

Created: Mar 16, 2023

Last Modified: Sep 29, 2023

PROTOCOL integer ID: 78911

Keywords: magattract, HMW DNA, arthropods, insects, high molecular weight, reference genome, long read

ABSTRACT

This is a protocol for the manual extraction of high molecular weight (HMW) DNA from insects. It uses the Qiagen MagAttract kit and factors in modifications described in the Chromium™Genome Reagent Kits User Guide pages 6-8 (<https://support.10xgenomics.com/genome-exome/library-prep/doc/user-guide-chromium-genome-reagent-kit-v2-chemistry>). It includes further modifications that benefited two particular needs. First, most reagents are halved per extraction (apart from beads) due to small specimen size (2 mg is a typical weight of an *Anopheles* mosquito, which is far less than the 25 mg of tissue the MagAttract kit can support). Second, due to this small specimen size, we need to maximize DNA yield, so we also perform a second elution to release more DNA.

DNA resulting from this protocol can be further sheared and cleaned for successful PacBio HiFi sequencing. In our experience, a single fresh *Anopheles gambiae* mosquito weighing 2-3 mg yields about 600-800 ng of DNA (quantified using qubit HS DNA kit) using this protocol. Following shearing using G-tubes or the MegaRuptor and SPRI based clean up, we typically retain about 200 ng of sheared DNA, which is just about sufficient to reach 25x coverage PacBio HiFi of a 250 Mb genome. The quality and quantity of DNA is best when starting with a snap frozen from living specimen. However, we have also successfully extracted HMW DNA using this protocol from ethanol and DESS preserved specimens held at room temperature for as long as a week, as long as the specimen was punctured or gently squished to ensure rapid penetration of the preservative. See: *Squishing insects for preservation of HMW DNA in the field* <https://www.protocols.io/view/squishing-insects-for-preservation-of-hmw-dna-in-t-cyp3xvqn>

We normally perform up to 8 DNA extractions in parallel (this will also depend on which magnetic rack is being used). Please add a comment to let the wider community know if it has worked or not on your species.

NB – an automated version of this protocol that works on the Kingfisher Apex is also available here on protocols.io

GUIDELINES

This protocol is an adaptation of Chromium™Genome Reagent Kits User guide and Qiagen's MagAttract® HMW DNA Handbook. The Chromium protocol suggests the use of PBS at the beginning. Also, both SOPs recommend using a thermomixer for wash steps, but we find doing this by hand works fine. Finally, we use 1.5 mL microcentrifuge tubes because the plastic pestles that we use to grind insect tissue do not support effective grinding in the 2 mL microcentrifuge tubes recommended by both original protocols. There are, however, pestles available to use in 2 mL microcentrifuge tubes.

MATERIALS



MagAttract HMW DNA
kit Qiagen Catalog #67563



Ethanol absolute Merck Millipore (EMD
Millipore) Catalog #107017



1X PBS (Phosphate-buffered saline
)



Dry Ice Contributed by users

Equipment

ThermoMixer® C

NAME

Eppendorf

BRAND

5382000031

SKU

https://www.eppendorf.com/gb-en/eShop-Products/Temperature-Control-and-Mixing/Instruments/Eppendorf-ThermoMixerC-p-PF-19703?gclid=Cj0KCQjwn9CgBhDjARIsAD15h0CN0IZBUgxtlqb7IFhPATCy01lk4tpBsMrqDjU7eNcDDsJLBbTjArkaAsb7EALw_wcB&gclsrc=aw.ds

LINK

Equipment

Mini-Centrifuge 100-240V, 50/60Hz Universal Plug, Grey

NAME

minicentrifuge

TYPE

Fisherbrand™

BRAND

16617645

SKU

<https://www.fishersci.co.uk/shop/products/fisherbrand-standard-mini-centrifuge/16617645>

LINK

Equipment	
DNA LoBind® Tubes	NAME
microcentrifuge tubes	TYPE
Eppendorf	BRAND
0030108051	SKU
https://www.eppendorf.com/gb-en/eShop-Products/Laboratory-Consumables/Tubes/DNA-LoBind-Tubes-p-0030108051	LINK

















Equipment	
Pestle for 1.5 mL Microtube, 100/pk	NAME
pellet pestle	TYPE
Cole-Parmer Essentials	BRAND
WZ-44468-19	SKU
https://www.coleparmer.co.uk/i/cole-parmer-essentials-pestle-for-1-5-ml-microtube-100-pk/4446819	LINK

Equipment	
1,000 µl graduated TipOne® Filter Tip, Natural, Racks (sterile), Case	NAME
Pipette tips	TYPE
Starlab	BRAND
S1126-7810	SKU
https://www.starlabgroup.com/GB-en/product/1000-ul-graduated-tipone-filter-tip-natural-sterile-pf-sl-920414.html?childSku=S1126-7810	LINK


Equipment	
200 µl Filter Tip / Wide Orifice	NAME
Pipette tips	TYPE
Starlab	BRAND
E1011-8618	SKU
https://www.starlabgroup.com/GB-en/product/200-ul-filter-tip--wide-orifice-e1011-8618.html	LINK

Equipment	
DynaMag™-2 Magnet	NAME
Magnetic tube rack	TYPE
Invitrogen™	BRAND
12321D	SKU
https://www.thermofisher.com/order/catalog/product/12321D#:~:text=The%20DynaMag%E2%84%A2%2D2%20magnet%20combines%20a%20strong%20magnetic%20attraction,microcentrifuge%20tubes%20in%20numbered%20spaces	LINK

PROTOCOL MATERIALS

 Ethanol absolute Merck Millipore (EMD Millipore) Catalog #107017	Materials
 Dry Ice Contributed by users	Materials
 RNase A Qiagen Catalog #67563	Step 2
 Buffer AL Qiagen Catalog #67563	Step 2
 Proteinase K Contributed by users Catalog #67563	Step 2
 gDNA 165kb Analysis Kit 275 Samples Agilent Technologies Catalog #FP-1002-0275	
Step 31	
 Buffer PE Qiagen Catalog #67563	In <u>2 steps</u>
 Buffer MB Contributed by users Catalog #67563	Step 10
 Buffer AE Qiagen Catalog #19077	Step 29
 Qubit® dsDNA HS Assay Kit Thermo Fisher Scientific Catalog #Q32854	
Step 31	
 MagAttract Suspension G Qiagen Catalog #67563	Step 9
 Buffer AE Qiagen Catalog #67563	In <u>2 steps</u>
 Quant-iT™ PicoGreen™ dsDNA Assay Kit Invitrogen - Thermo Fisher Catalog #P11496	
Step 31	
 MagAttract HMW DNA kit Qiagen Catalog #67563	Materials
 Buffer MW1 Contributed by users Catalog #67563	Step 13
 Nuclease-Free Water Qiagen Catalog #67563	Step 23

SAFETY WARNINGS

-  Buffers AL, MB, and MW 1 contain guanidine hydrochloride/guanidine thiocyanate, which can form highly reactive compounds when combined with bleach. DO NOT add bleach or acidic solutions directly to the sample preparation waste. Waste needs to be collected in a suitable vessel and disposed of in accordance with local regulations.

BEFORE START INSTRUCTIONS

Ensure all surfaces have been cleaned with 70-80% Ethanol (and ideally bleach before that). Have cleaning wipes for forceps. All kit components, buffers, and RNase A stock solution can be stored at room temperature (15–25°C) for up to 1 year. The box should be labelled with received date. Mix Buffer AL thoroughly by shaking before use. Buffers MW1 and PE are supplied as a concentrate. Before using for the first time, be sure to add the appropriate amount of ethanol (96–100%) as indicated on the bottle. Many components of the kit are also available from Qiagen separately.

Procedure

2h 24m




- 1 Prepare an open insulated box of dry ice to store sample tubes on whilst working through steps 2-4.

- 2 Make mastermix of reagents for lysis.




Note


Make sure to choose the right size of tube for preparing the mastermix.


Calculate the mastermix volumes for the number of samples plus 1 for spare pipetting volume.

Volumes per sample are:  100 μ L  1X PBS (Phosphate-buffered saline) ;  10 μ L

 Proteinase K Contributed by  (Mix by inverting the tube 5 times);
users Catalog #67563

 2 μ L  RNase A Qiagen Catalog #67563 ;  75 μ L

 Buffer
AL Qiagen Catalog #67563 (mix by inversion)

- 3 For each sample, add  187 μ L of the mastermix from step 2 into a new 1.5 mL DNA LoBind tube.

Note

A 2 mL DNA LoBind tube can also be used but an appropriate size pestle will have to be used.


- 4 Carefully remove the **Sample** from the sample tube using clean forceps. If the **Sample** has been stored in a preservation liquid, lightly make contact on a clean piece of tissue to remove surface liquid from the sample. Submerge the **Sample** into the mastermix in a tube (see step 3) with clean forceps. Insert a sterile pestle in the tube and smash, smear, squash, twist, grind the tissue against the wall of the tube for **00:01:00**. There should be no recognisable body parts visible following pestle smashing, only flakes. Place the sample in a tube rack on the bench. Clean forceps with 100% ethanol.





Mosquito in lysis buffer.






Tissue disruption with an autoclavable pestle.








Mosquito debris after tissue disruption.

- 5 Repeat step 4 for the remaining samples.
- 6 Briefly spin all samples in a minicentrifuge or similar to collect solution at bottom before next step.
- 7 Incubate the sample at  25 °C for  02:00:00 . 2h
- 8 Briefly spin samples to collect solution at the bottom of the tube.

- 9 Vortex the  MagAttract Suspension
G Qiagen Catalog #67563 for  00:01:00 and add  15 µL to each sample. 1m

Note

If this is the first time using  MagAttract Suspension
G Qiagen Catalog #67563 ,
increase the vortexing time to
 00:03:00 . Briefly vortex  MagAttract Suspension
G Qiagen Catalog #67563 before
adding to each subsequent sample.

- 10 Add  140 µL  Buffer MB Contributed by users Catalog #67563 to each sample. Mix by gentle inversion, fully invert but don't shake. If you see the beads making flakey clumps, a little like gold leaf, this is a good sign. This is difficult to do simultaneously as you want to see the mixing. If you are doing 8 samples, give all samples another gentle inversion after the last one. 1m


Leave at least ⌚ 00:01:00 for the beads to bind. Doing multiple samples will often take more time than this.



Flakes of beads building in the lysis buffer after addition of buffer MB and careful inverting of tube.

- 11 Centrifuge the tube briefly and place on a DynaMagTM-2 Magnetic Rack for ⌚ 00:01:00 to

1m




Beads are collecting on the magnetic side.






Note






We use a DynaMag™-2 Magnetic Rack but other magnetic racks suitable for 1.5 or 2 mL microcentrifuge tubes will work as well.

- 12 Remove and discard the supernatant. Take care not to disturb the bead pellet.




Supernatant being removed without disturbing the pellet.

- 13 Remove the  Sample from the magnetic rack. Add  350 μL  Buffer MW1 Contributed by users Catalog #67563 directly to the bead pellet. Mix by inversion, ensuring that the beads have come away from the side of the tube. This often requires tapping the tube or swilling the contents. Try to be as gentle as possible. Again, this is difficult to do simultaneously as you need to check each sample.
- 14 Centrifuge the tube briefly and place on a DynaMagTM-2 Magnetic Rack for  00:01:00 to  allow bead capture.
- 15 Remove and discard the supernatant. Take care not to disturb the bead pellet.
- 16 Repeat steps 13-15 for a total of 2 washes.

- 17 Remove the sample from the magnetic rack. Add  350 µL  Buffer PE Qiagen Catalog #67563 directly to the bead pellet. Mix by inversion, ensuring that the beads have come away from the side of the tube. This often requires tapping the tube or swilling the contents. Try to be as gentle as possible. Again, this is difficult to do simultaneously as you need to check each sample.
- 18 Centrifuge the tube briefly and place on a DynaMagTM-2 Magnetic Rack for  00:01:00 to 
- 19 Remove and discard the supernatant. Take care not to disturb the bead pellet.
- 20 Repeat steps 17 and 18 for a total of 2 washes but do not remove supernatant immediately, proceed to step 21.
- 21 If you have more than four tubes, split them into groups of four or fewer. Perform steps 22-24 (water wash) on the first group of samples, whilst the remaining samples wait in  Buffer PE Qiagen Catalog #67563 on the magnet. [While the first group of samples are eluting in step 25-26 below, it is possible to perform the water wash on the second group of samples.]
- 22 Remove and discard the supernatant. Take care not to disturb the bead pellet. With a P20 pipette remove any remaining supernatant. Leave the sample on the magnetic rack for the next step. Do not pipette water directly onto the beads.








Note







The timing of the next step is extremely important. If a multichannel pipette is not available, ensure that each tube has the exact same incubation time. Do not exceed  00:01:00 .




Water being pipetted against the side opposite of the magnetic beads to avoid disturbing beads.

- 23** Carefully add 350 μ L Nuclease-Free Water Qiagen Catalog #67563 down the side of the tube opposite the magnetic pellet. Start a timer counting up from zero. After 00:00:15 add water to the second sample, after 00:00:30 add to the third sample, after 00:00:45 add to the fourth sample. At 00:01:00 remove and discard the water from the first sample, at 00:01:15 remove and discard water from the second sample, at 00:01:30 the third, and 00:01:45 the fourth. This will enable multiple samples to be incubated for exactly 00:01:00 . 8m
- 24** Repeat step 23 for a total of 2 washes.
- 25** Remove the samples from the magnetic rack. Add 100 μ L Buffer AE Qiagen Catalog #67563 directly to the bead pellet of each sample. Ensure that the pellet is submerged and has come away from the side of the tube. Incubate at 25 °C and 1400 rpm shaking for 00:03:00 . 3m

- 26 During this  00:03:00 incubation perform steps 22-25 on the second group of samples if present. 3m
- 27 Centrifuge each tube briefly and place them on a magnetic rack for  00:01:00 to allow bead capture. 1m
- 28 Using a wide-orifice pipette tip, carefully transfer the supernatant containing purified gDNA to a new labelled 1.5 mL LoBind microcentrifuge tube or barcoded (e.g. FluidX) tube.
- 29 Second elution: remove the samples from the magnetic rack. Add  100 μ L 3m
 Buffer AE Qiagen Catalog #67563 directly to the bead pellet of each sample. Ensure that the pellet is submerged and has come away from the side of the tube. Incubate at  25 °C and  1400 rpm shaking for  00:03:00 .
- Note**

Due to the second elution step, an additional bottle of  Buffer AE Qiagen Catalog #19077 will be necessary if you buy the kit.
- 30 Using a wide-orifice pipette tip, carefully transfer the supernatant containing purified gDNA to the same 1.5 ml LoBind microcentrifuge tube or Fluidx tube with a final volume of  200 μ L .
- 31 Store the extracted gDNA sample at  4 °C . Assess the quantity of DNA extracted using the  Qubit® dsDNA HS Assay Kit Thermo Fisher Scientific Catalog #Q32854 or  Quant-iT™ PicoGreen™ dsDNA Assay Kit Invitrogen - Thermo Fisher Catalog #P11496 and assess the quality of the DNA using the Femto Pulse  gDNA 165kb Analysis Kit 275 Samples Agilent Technologies Catalog #FP-1002-0275 .



Example of a Femtopulse profile of DNA extracted from a single snap frozen *Anopheles* mosquito with this protocol.