

Jun 28, 2024

Yale Murine TMC - H&E staining following Phenocycler-Fusion Imaging

DOI

dx.doi.org/10.17504/protocols.io.ewov19zjklr2/v1

Jungmin Nam¹, Yale Pathology Tissue Services¹, Rong Fan¹, Archibald Enniful¹

¹Yale University

Cellular Senescence Net...



Jungmin Nam

Yale University

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.ewov19zjklr2/v1

Protocol Citation: Jungmin Nam, Yale Pathology Tissue Services, Rong Fan, Archibald Enniful 2024. Yale Murine TMC - H&E staining following Phenocycler-Fusion Imaging. [protocols.io <https://dx.doi.org/10.17504/protocols.io.ewov19zjklr2/v1>](https://dx.doi.org/10.17504/protocols.io.ewov19zjklr2/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 27, 2024

Last Modified: June 28, 2024

Protocol Integer ID: 102545

Keywords: phenocycler-fusion, codex, h&e

Funders Acknowledgement:

Vishwa Deep Dixit

Grant ID: U54 AG079759



Abstract

Flow cell removal after Phenocycler-Fusion imaging and following H&E staining steps from Yale Pathology Tissue Services (YPTS).



Removing flow cell following CODEX imaging

- 1 Remove slides from storage buffer.
- 2 Immerse slides to xylene for 10 mins.
- 3 Separate flow cell from the slide using a razorblades.
- 4 Rinse slide with 1x PBS.
- 5 Store in storage buffer until H&E staining.

10m

Hematoxylin and eosin staining

11m 30s

- 6 Stain with hematoxylin for 5 min.
- 7 Rinse with water.
- 8 Clarifier for 30 seconds.
- 9 Rinse with water.
- 10 Bluing for 1 min.
- 11 Rinse with water.

5m

30s

1m



12 Stain with eosin 5 min.

5m

13 Immerse slides in graded ethanol to xylene.

14 Add coverslip with resin mounting medium.