



Aug 17, 2021

# Protocol for Image processing and analysis of VPS13D recruitment to mitochondria

Marianna Leonzino<sup>1,2</sup>, Andrés Guillén-Samander<sup>1,2</sup>, Ni Tang<sup>1,2</sup>, Pietro De Camilli<sup>1,2</sup>

<sup>1</sup>Departments of Neuroscience and of Cell Biology, Howard Hughes Medical Institute, Program in Cellular Neuroscience, Neurodegeneration and Repair, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

<sup>2</sup>Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.bvgnn3ve](https://doi.org/10.17504/protocols.io.bvgnn3ve)

Andrés Guillén-Samander

## ABSTRACT

This protocol details the image processing and analysis of VPS13D recruitment to mitochondria as it was performed in <https://doi.org/10.1083/jcb.202010004>. The first part details the analysis of acute optogenetic recruitment, and the second part of basal recruitment under different conditions.

## ATTACHMENTS

[dn3wbgtzx.pdf](#)

## DOI

[dx.doi.org/10.17504/protocols.io.bvgnn3ve](https://doi.org/10.17504/protocols.io.bvgnn3ve)

## EXTERNAL LINK

<https://doi.org/10.1083/jcb.202010004>

## PROTOCOL CITATION

Marianna Leonzino, Andrés Guillén-Samander, Ni Tang, Pietro De Camilli 2021. Protocol for Image processing and analysis of VPS13D recruitment to mitochondria. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bvgnn3ve>

## KEYWORDS

Image processing, VPS13D recruitment, Mitochondria, ASAPCRN

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Jun 02, 2021

## LAST MODIFIED

Jul 02, 2021

## OWNERSHIP HISTORY

Jun 02, 2021 Urmilas

Jun 08, 2021 Andrés Guillén-Samander

## PROTOCOL INTEGER ID

50414

Fluorescence images were processed using FIJI (ImageJ; NIH) software. Gaussian blur filters were applied on some of the images presented.

#### For the analysis of optogenetic experiments:

- 1 Build kymographs tracing a line ROI across the mitochondria that was illuminated with Blue Light.
- 2 Measure an intensity profile by tracing a line ROI across the center of the mitochondrial signal on the kymograph. Take a background profile using the same line ROI outside the illumination area.
- 3 Normalize the fluorescence intensities to the min and max values using the following formula:

$$\text{Normalized } f[t] = \frac{f[t] - f_{\min}[t]}{f_{\max}[t] - f_{\min}[t]}$$

where  $f[t] = F_{\text{mito}}[t] - F_{\text{bkgd}}[t]$ .

#### For the analysis of VPS13D mitochondrial recruitment:

- 4 Analyze the cells coexpressing VPS13D<sup>EGFP</sup> and mitoBFP individually. Store the each cell outline as an ROI in the ROI Manager of Fiji.
- 5 Apply the "Gaussian Blur" tool and "Enhance Local Contrast" tool to the MitoBFP channel to average the brightness of the whole image and to clarify the shape and margin of mitochondria.
- 6 Use an automated macro designed in Fiji to measure VPS13D<sup>EGFP</sup> signal overlapping with mitochondria and that localized in the cytosol surrounding them.
- 7 Calculate the ratio between the intensity of VPS13D<sup>EGFP</sup> on mitochondria and the intensity of VPS13D<sup>EGFP</sup> in the cytosol surrounding mitochondria for each cell and plot a graph indicating the enrichment of VPS13D<sup>EGFP</sup> on mitochondria.

The macro (which can be adapted to analyse recruitment of any given protein to mitochondria) is included below:

```
macro "ROI in batches auto measure VPS13 binding with mitochondrial " {
```

```
Dialog.create("reminder");//create an interface
```

```
Dialog.addMessage("choose saving path");//an interactive interface for user to set the saving directory
```

```
Dialog.show();//an interactive interface for user to set the saving directory
```

```

    dir_saving = getDirectory("Choose a Directory to save");//an interactive interface for user to set
the saving directory

count = roiManager("Count");//save the information from ROI manager into the term "count"

for (i = 0; i < count; i++){

    name = getTitle;

    selectWindow(name);

    roiManager("Select", i);//save the selected ROIs in the ROI manager in batch

    Stack.getPosition(channel, slice, frame);//save the selected ROIs in the ROI manager in batch

    run("Duplicate...", "duplicate slices=" + slice);//save the selected ROIs in the ROI manager in
batch

    saveAs("Tiff", dir_saving + "ROIslice"+i+name); //save the selected ROIs in the ROI manager in
batch

    name = replace(name, ".tif", ""); //motify the name of images

    C1 = "C1-ROIslice"+i+name+".tif";//motify the name of images

    C2 = "C2-ROIslice"+i+name+".tif";//motify the name of images

    run("Split Channels");//split the integrated image into single channel

    Dialog.create("Content of the channel");//create an interface

        Dialog.addMessage("Make sure the content of the channels. If same as described, click
OK. If different, click Cancel and make changes.\n");

        Dialog.addMessage("Channel 1 is VPS13D channel. Channel 2 is Mitochondrial
channel.");
        Dialog.show();//make sure the content of each channels

    selectWindow(C1);//choose the cytosol slice

    run("Duplicate...", "");//duplicate a seperate image to process

    run("Brightness/Contrast...");//preprocess the VPS13D signal in the cytosol area

//measure the intensity of VPS13D signal in mitochondria area

    selectWindow(C2);//choose the mitochondria slice

    run("Duplicate...", "");//duplicate a seperate image to process

    run("Gaussian Blur...", "sigma=2 slice");//balance the brightness of image

    run("Enhance Local Contrast (CLAHE)", "blocksize=127 histogram=256 maximum=3
mask=*None* fast_(less_accurate)");//clarify the outline of mitochondria

    run("Enhance Local Contrast (CLAHE)", "blocksize=127 histogram=256 maximum=3

```

```

mask=*None* fast_(less_accurate)");//clarify the outline of mitochondria

run("Enhance Local Contrast (CLAHE)", "blocksize=127 histogram=256 maximum=3
mask=*None* fast_(less_accurate)");//clarify the outline of mitochondria

run("Threshold...");//choose the mitochondria area

run("Convert to Mask");//generate an accurate mitomask

saveAs("Tiff", dir_saving + "ROI"+i+" mitomask "+name);//generate an accurate mitomask

imageCalculator("AND create", mitomask, C1);//map C1 onto the mitomask

setThreshold(1, 255);//choose the calculated area

saveAs("Tiff", dir_saving + "ROI"+i+" result of mitochondria intensity "+name);

mitointensity=getTitle;

run("Measure");//measure the intensity of VPS13D signal in mitochondria area

//measure the VPS13D intensity in the cytosol area

run("Dilate");

run("Dilate");

run("Dilate");//duplicate three times to eliminate the interference from mitochondria signal

run("Convert to Mask");

run("Invert");//invert the enlarged mitomask

saveAs("Tiff", dir_saving + "ROI"+i+" mitomask x3 "+name);//generate the enlarged inverted mitomask

mitomask3 = getTitle;

imageCalculator("AND create", mitomask, mitomask3);//map the enlarged inverted mitomask
onto the original mitomask to get the enlarged outline of mitochondria in the cytosol area

saveAs("Tiff", dir_saving + "ROI"+i+" cytomask "+name);//generate outline profile of mitomask,
which refers to mask of cytosol area

cytomask = getTitle;

imageCalculator("AND create", cytomask, C1);//map the cytomask onto cytosol

setThreshold(1,255);

saveAs("Tiff", dir_saving + "ROI"+i+" result of cytosolic intensity "+name);

cytointensity=getTitle;

run("Measure");//measure the VPS13D intensity in the cytosol area

saveAs("Results", dir_saving + "ROI"+i+" measurement results of VPS13 binding to
mitochondrial.csv");

```

```
        close()

    }

    Dialog.create("reminder");

        Dialog.addMessage("Mission accomplished.Congratulations!");

        Dialog.show();

    }
```