

SEP 30, 2023

# OPEN ACCESS



Protocol Citation: Annan SI Cook 2023. Production of GTPase Deficient RAB1A(Q70L) Protein. protocols.io

https://protocols.io/view/production-of-gtpase-deficient-rab1a-q70l-protein-c2suyeew

License: This is an open access protocol distributed under the terms of the Creative Commons Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working We use this protocol and it's working

Created: Sep 18, 2023

Last Modified: Sep 30,

2023

**PROTOCOL** integer ID:

88628

## Production of GTPase Deficient RAB1A(Q70L) Protein

Annan SI Cook<sup>1,2</sup>

<sup>1</sup>University of California, Berkeley; <sup>2</sup>ASAP CRN



Annan SI Cook

**ABSTRACT** 

Bacterial expression of the GTPase deficient Q70L mutant of human RAB1A.

**ATTACHMENTS** 

844-2181.pdf

### **Materials and Reagents**

- BL21 strain of *E. coli*
- pET vector containing 6x-His-RAB1A(Q70L)
- Ampicillin-resistant plates
- LB broth
- Isopropyl β-D-1-thiogalactopyranoside (IPTG)
- Phosphate Buffered Saline (PBS)
- Lysis buffer
- EDTA-free protease inhibitor tablet (Thermo Fisher)
- Phenylmethylsulfonyl fluoride (PMSF)
- Ni-NTA column
- Wash buffer (containing 300 mM imidazole)
- S75 10/300 column
- Liquid nitrogen

### **Lysis buffer**

A	В
HEPES pH 7.5	25 mM
NaCl	300 mM
MgCl <sub>2</sub>	2 mM
TCEP	1 mM

## **Transformation and Colony Selection**

- 1 Transform *E. coli* BL21 cells with the pET vector containing 6x-His-RAB1A(Q70L) by heat shock.
- 2 Plate the transformed cells on ampicillin-resistant plates.

3 Incubate the plates Overnight at \$\mathbb{I}\$ 37 °C



## **Pre-culture and Expansion**

- 4 Inoculate a single colony into 🔼 10 mL of LB broth.
- 5 Incubate the culture Overnight at 37 °C



## **Main Culture Preparation**

- **6** Transfer the entire pre-culture into a 1 L culture flask.
- Grow the A 1 L culture at 37 °C with shaking at 5 220 rpm until the optical density (OD) reaches approximately 0.6.

## **Cooling and Induction**

15m

Place the flask in an ice bath for 00:15:00 to cool rapidly.

10111

9 Induce protein expression by adding [M] 0.1 millimolar (mM) IPTG

10 Continue incubating the culture Overnight (approximately 16 hours) at \$\mathbb{E}\$ 18 °C with shaking at (5 180 rpm **Cell Harvest and Pelleting** 11 Harvest the cells by centrifugation at 4000 rpm in a BECKMAN Coulter centrifuge. **Cell Washing and Freezing** 12 Wash the cell pellet once with PBS. 13 Pellet the cells again at 4000 x g 14 Flash freeze the pellet in liquid nitrogen. Store the frozen pellet at [ -80 °C until purification. 15 6m **Cell Lysis** 

Thaw the frozen cell pellet to Room temperature

16

17 Resuspend the cells in lysis buffer.

### **Lysis buffer**

A	В
HEPES pH 7.5	25 mM
NaCl	300 mM
MgCl <sub>2</sub>	2 mM
TCEP	1 mM

- Add an EDTA-free protease inhibitor tablet (Thermo Fisher) and PMSF to a concentration of IMI 1 millimolar (mM).
- Lyse the cells by sonication using a 5-second on, 5-second off cycle at 50% power for a cumulative time of 00:06:00.

6m

## **Clarification by Centrifugation**

45m

20 Centrifuge the lysate at 17000 rpm, 00:45:00 to clarify.

15m



### **Protein Purification on Ni-NTA Column**

- 21 Apply the supernatant to a Ni-NTA column (3x).
- Wash the column until the wash buffer is free of protein, as measured by a Bradford assay.



Elute the protein using a wash buffer containing [M] 300 millimolar (mM) imidazole.

## **Size Exclusion Chromatography**

- 24 Concentrate the eluted protein.
- Apply the concentrated protein to an S75 10/300 column in buffer containing [M] 25 millimolar (mM) HEPES (PH 7.5), [M] 150 millimolar (mM) NaCl, [M] 1 millimolar (mM) MgCl<sub>2</sub>, and [M] 1 millimolar (mM) TCEP.

### **Final Protein Concentration**

- 26 Concentrate the purified protein to a final concentration of [M] 50 micromolar (µM)
- 27 Divide the protein into A 50 µL aliquots.

## Cryopreservation

- 28 Flash freeze the protein aliquots in liquid nitrogen.
- Store the frozen protein aliquots at \$\circ\$ -80 °C for future use.

