

Jun 02, 2021

Protocol to process Gastruloids for FACS (reporters only)

Stefano Vianello¹, Tania Hübscher¹, Giuliana Rossi¹, Matthias Lutolf²

¹Laboratory of Stem Cell Bioengineering, Institute of Bioengineering, School of Life Sciences (SV) and School of Engineering (STI), École Polytechnique Fédérale de Lausanne (EPFL), Lausanne, Switzerland;

²Laboratory of Stem Cell Bioengineering, Institute of Bioengineering, School of Life Sciences (SV) and School of Engineering (STI), École Polytechnique Fédérale de Lausanne (EPFL), Lausanne, Switzerland and Institute of Chemical Sciences and Engineering, School of Basic Science (SB), EPFL, Lausanne, Switzerland

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.bvgrn3v6



Stefano Vianello
EPFL - EPF Lausanne

ABSTRACT

Protocol to process reporter Gastruloids for FACS analysis of their endogenous reporters.

DOI

dx.doi.org/10.17504/protocols.io.bvgrn3v6

PROTOCOL CITATION

Stefano Vianello, Tania Hübscher, Giuliana Rossi, Matthias Lutolf 2021. Protocol to process Gastruloids for FACS (reporters only) . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bvgrn3v6>

KEYWORDS

Gastruloids, FACS

LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jun 02, 2021

LAST MODIFIED

Jun 02, 2021

PROTOCOL INTEGER ID

50417

GUIDELINES

This protocol assumes standard aseptic technique, and work in a laminar flow hood.

MATERIALS TEXT

Recipes:

Digestion solution (2mL per sample)

final concentrations: Collagenase IV [3mg/mL], Dispase [4mg/mL], DNase I [100ug/mL] in PBS

- [3mg/mL] **Collagenase IV** [CAT#17104-019; Gibco™/Thermo Fisher Scientific] 1:4 from 12mg/mL stock
- [4mg/mL] **Dispase in H2O** [CAT#17105-041; Gibco™/Thermo Fisher Scientific] 1:2.5 from 10mg/mL stock
- [100ug/mL] **DNase I** [CAT#11284932001; Roche] 1:500 from 50mg/mL stock
- top up with **PBS**-/-

Staining buffer

Final concentrations: 10% ES-FBS, [100U/mL] Pen-Strep, [1mM] EDTA in PBS

- 1:10 final volume: **ES-FBS** [CAT#16141079; Gibco™/Thermo Fisher Scientific]
- 1:100 final volume: **Pen-Strep** [CAT#15140122; Gibco™/Thermo Fisher Scientific]
- 1:500 final volume: 500mM **EDTA** [CAT#15575020; Gibco™/Thermo Fisher Scientific]
- top up with **PBS**-/- [CAT#10010056; Gibco™/Thermo Fisher Scientific]

2% PFA in PBS-/-

(e.g. diluted 1:2 in PBS-/- from 4% PFA solution, CAT#15434389, Alfa-Aesar/Fisher Scientific)

BEFORE STARTING


Prepare enough FACS tubes/samples/controls to have a complete setup with both positive, negative, FMO (Fluorescence Minus One) references.

For a simple FACS setup of Gastruloids with a double reporter, in addition to

- your experimental samples (i.e. dissociated Gastruloids for each timepoint, with DAPI as live-dead discriminant)

you will thus also need

- double negative cells (ideally, the parental cells of the reporter, with no fluorescent insert; but use alternatively 2D cultured reporter cells if they do not express the reporters in the pluripotent state)
- cells that are RFP+ and have no GFP reporter (or RFP+ beads)
- cells that are GFP+ and have no RFP reporter (or GFP+ beads)
- Make sure to have a sample that not only has no reporter expression at all, but also no DAPI staining (to help set the gate for the DAPI signal).



Gastruloid collection

- 1 Harvest Gastruloid from each well of a 96 well plate and collect them in a 15mL Falcon tube (pooled by timepoint or condition)
- 2 Once all Gastruloids have settled to the bottom of the tube, aspirate out the supernatant (N2B27 medium carried over with each Gastruloids)
- 3 Resuspend the Gastruloids in 5mL PBS^{-/-}, to wash away traces of N2B27
- 4 Once all Gastruloids have re-settled to the bottom of the tube, aspirate out the PBS^{-/-}. Gastruloids are ready to be digested.

Gastruloid digestion

8m

- 5 Digest the Gastruloids by adding 1mL **Digestion Solution** (3mg/mL Collagenase IV, 4mg/mL Dispase, 100ug/mL^{4m} DNaseI, in PBS^{-/-}), **37 °C (waterbath)**, **00:04:00** min
- 6 After the 4min incubation has elapsed, use a P1000 coated in Staining Buffer to try to mechanically disrupt the Gastruloids
- 7 Place the tube back again at **37 °C (waterbath)**, **00:04:00** min; to complete the digestion 4m

- 8 After the last 4min of incubation have elapsed, use a P1000 coated in Staining Buffer (see recipe in "Materials" section) to mechanically disrupt the Gastruloids. Gastruloids should break easily and give rise to a single cell suspension.

Cell filtering 4m

- 9 Prewet the blue cap of a filter-FACS tube (e.g. CAT#352235, Falcon/Corning), by pipetting 1mL of Staining Buffer through it.



This step is very important. Not prewetting the filter can lead to loss of a high number of cells at the filtering step.

- 10 *Working* ⚠ **On ice**

Filter your cell suspension (digested Gastruloids) through the cap of the FACS tube

- 11 ⚠ **On ice**

Wash the filter by passing 1mL more of Digestion Solution through it.

- 12 ⚠ **On ice**

Transfer the contents of the FACS tube to a clean 15mL Falcon tube.

- 13 ⚠ **On ice**

Collect any leftover cells by flushing the old FACS tube with 2mL Staining Buffer, and transferring in the 15mL Falcon tube with the rest of the cells

- 14 ⚠ **On ice**

Add 5mL more Staining buffer into the tube Falcon tube with your cells suspension, to completely stop the digestion reaction. This tube now contains a filtered, single-cell suspension of your sample for FACS.

- 15 Spin down the cell suspension, 🌀 **200 x g, 4°C** , ⌚ **00:04:00** min

4m

Cell collection 4m

- 16 For samples that come from 2D culture (e.g. 2D-grown reporter cells as negative controls, or 2D-grown single-color ^{9m} cells as positive controls):

- aspirate out the medium from the culture dish
- add 3mL PBS/- to wash away any trace of medium, and aspirate it out again
- Add 500uL Accutase, to detach the cells, ⌚ **00:05:00** min, ⚠ **Room temperature**
- Transfer the detached cells +Accutase to a clean Falcon tube
- Add ~5mL 10% Serum to stop the reaction
- Spin down the cell suspension, 🌀 **200 x g, 4°C** , ⌚ **00:04:00** min

The cells are ready to be processed just as any Gastruloid sample (proceed to the next step)

DAPI staining and fixation



23m

- 17 For each sample, prepare 500uL of 0.2ug/mL DAPI in Staining Buffer

From a stock solution of DAPI 1mg/mL, use 1:5000.
E.g. add 0.1uL DAPI in 500uL Staining Buffer

- 18  **On ice**

15m

Resuspend the cell pellet in the DAPI solution, and incubate  **00:15:00** min,  **On ice**

- 19 Spin down the cell suspension,  **200 x g, 4°C** ,  **00:04:00** min

4m

- 20 **Fixation:**

 **On ice**

Resuspend the DAPI-stained cell pellet in 2mL 1%PFA,


- 21 Spin down the cell suspension,  **200 x g, 4°C** ,  **00:04:00** min

4m

- 22 Resuspend in 300uL Staining Buffer and transfer to a clean FACS tube pre-coated with Staining Buffer

Processing/Storage

23m

- 23 Cells can be kept in Staining Buffer,  **4 °C** , in the dark (wrapped in aluminium foil) for at least a week (e.g. while waiting to collect samples from later on timepoints)

- 24 When ready for FACS,  **On ice** , distribute your cells in 300uL per FACS tube and proceed for sorting/analysis.


Beads samples

23m

- 25 If using beads as positive references, prepare bead tubes just before FACS processing.

10s

For GFP BrightComp eBeads™ Compensation Bead Kit (Invitrogen/Thermo Fisher Scientific, CAT#A10514):

- vortex vial for  **00:00:10** s
- Dispense 1 drop into a clean FACS tube
- Resuspend in 1mL Staining Buffer