



Oct 09, 2024

## Permanent specimen preparation by protargol staining

DOI

[dx.doi.org/10.17504/protocols.io.q26g71or9gwz/v1](https://doi.org/10.17504/protocols.io.q26g71or9gwz/v1)

Marek Valt<sup>1</sup>, Michael Kotyk<sup>1</sup>

<sup>1</sup>Charles University, Faculty of Science, Department of Zoology, Prague, Czechia

SUM-K



Marek Valt

Charles University, Faculty of Science, Department of Zoolog...

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.q26g71or9gwz/v1](https://doi.org/10.17504/protocols.io.q26g71or9gwz/v1)

**Protocol Citation:** Marek Valt, Michael Kotyk 2024. Permanent specimen preparation by protargol staining. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.q26g71or9gwz/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 09, 2024

**Last Modified:** October 09, 2024


**Protocol Integer ID:** 109479

### Abstract

Protargol staining according to the protocol by Bodian (Bodian 1936) and later modified by Nie (Nie 1950), further optimized for staining of small flagellates.

## Protargol staining optimized for a small flagellates, used for staining of *Solarion arienae*.

2d 1h 15m

- 1 Mix 1  $\mu\text{m}$  of culture sediment with a drop of sterile egg albumin, on a coverslip.  
For *Solarion arienae*, egg albumin was not diluted to achieve better adherence to the cells to coverslip.
- 2 Fix in Bouin-Hollande's solution overnight. (0,175 mol/l picric acid, 0,138 mol/l  $\text{Cu}^{+2}$  acetate, 4% formaldehyde, and 5% acetic acid in water solution) 
- 3 Next day, transfer coverslips throught a graded ethanol series (70%, 50%).
- 4 Wash in distilled water.
- 5 Bleach in 0,5%  $\text{KMnO}_4$  solution for 5 minutes. 5m
- 6 Wash in distilled water (5x).
- 7 Treat coverslips with 5% oxalic acid for 5 minutes. 5m
- 8 Wash again in distilled water (5x).
- 9 Stain in 1% protargol solution, for 48 hours at 37°C in a beaker with copper wire pieces placed between the coverslips. 2d  
  
For *Solarion*, protargol produced by Bayer, I. G. Farbenindustrie Actinengesellschaft (out of business since 1952) was used.
- 10 Wash in distilled water (2x).
- 11 Treat with freshly prepared reducing solution (1% hydroquinone and 5%  $\text{Na}_2\text{SO}_3$ ), for 10 minutes. 10m



- 12 Wash in distilled water (5x).
- 13 Tone with 1%  $\text{AuCl}_3$  for 5 minutes. 5m
- 14 Wash in distilled water (2x).
- 15 Treat with 2% oxalic acid, for 5 minutes. 5m
- 16 Wash in distilled water (5x).
- 17 Treat with 5%  $\text{Na}_2\text{S}_2\text{O}_3$  for 10 minutes 10m
- 18 The final wash is done under a constant stream of tap water for 20 minutes. 20m
- 19 Dehydrate the coverslips in an ethanol series (50%, 70%, 80%, 96%, 100%) and 3 times in xylene (5 minutes each).
- 20 Dehydrate in xylene, 3 times (5 minutes each). 15m
- 21 Finally, the stained, dehydrated coverslips are mounted on glass slides with DPX mounting medium (Sigma–Aldrich).