



# RNA Integrity Check V.1

Elizabeth Fozo<sup>1</sup>

<sup>1</sup>In-house protocol

Version 1

Nov 23, 2020

1

Works for me

This protocol is published without a DOI.

Eadewunm

## ABSTRACT

RNA Integrity Check: Modification of Protocol from Qiagen

## PROTOCOL CITATION

Elizabeth Fozo 2020. RNA Integrity Check. **protocols.io**  
<https://protocols.io/view/rna-integrity-check-bpwjmpcn>

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Nov 19, 2020

## LAST MODIFIED

Nov 23, 2020

## PROTOCOL INTEGER ID

44715

## SAFETY WARNINGS

Formamide and formaldehyde are toxic and dangerous chemicals. Handle with care in the hood!

## DISCLAIMER:

DISCLAIMER: This is a work in progress.

## ABSTRACT

RNA Integrity Check: Modification of Protocol from Qiagen

## BEFORE STARTING

Pretreat gel containers and gel flask with RNase Zap

## RNA Integrity Check

- 1 Determine the concentration of RNA; ideally want to run approximately **1 µg** of RNA on the gel (should be about **1 µl** of your sample)
- 2 Spray gel tray, comb, dams, and a small flask with RNase Zap.
- 3 For small gels, melt 0.3 grams of RNA ONLY agarose in 30 ml 1X TBE (or 1X TAE) to the cleaned flask (for a 1% gel using Horizon gel boxes. If using something else, adjust accordingly).

4 Let the solution cool slightly, add **1 µl EtBr** . Pour gel and let solidify.

5 For samples:

1. RNA sample in a total of **8 µl** of water (if RNA concentration is about **1 µg /ml** , add **1 µl** of RNA to **7 µl** of RNase-free water)
2. **4 µl** of sample buffer
3. Mix samples

Sample Buffer (make FRESH!)

- **80 µl of concentrated dye** (80% glycerol, 0.2% bromophenol blue, 0.2% xylene cyanol; this does NOT need to be made fresh!)
- **500 µl of formamide**
- **120 µl of 37% formaldehyde**

**(1/10 of the recipe above is normally enough)**

6 Load gel and include a 1 kb DNA ladder for size approximation; run at 100 volts about 45-60 minutes

7 Gel: should see 2 distinct bands corresponding to the 16s and 23s rRNA species. If RNA is intact, will not see smearing. **Should also see lots of tRNA at the bottom.**



Formamide and formaldehyde are toxic and dangerous chemicals.  
Handle with care in the hood!

On occasion, I have run a mini-1X MOPS gel to validate RNA integrity especially if I do not have much sample left and really needed to use the RNA for an important experiment.