



Dec 23, 2021

🌐 Isopropanol DNA Precipitation - Best for DNA Concentration

Lynn Doran¹¹Realizing Increased Photosynthetic Efficiency (RIPE)

1

dx.doi.org/10.17504/protocols.io.bxg7pjzn**Burgess Lab UIUC**

Lynn Doran

Realizing Increased Photosynthetic Efficiency (RIPE)

Protocol to precipitate extracted DNA from an aqueous solution to increase concentration or resolubilize in a different storage buffer.

*Use isopropanol DNA precipitation if your DNA is suspended in a very large volume, if your DNA concentration is low, or you are trying to concentrate large molecular weight DNA fragments and remove smaller fragments.

*Use ethanol DNA precipitation if you are trying to remove salt contamination or precipitate small DNA fragments.

Additional Resources:

- [New England Biolabs, "DNA Precipitation: Ethanol vs. Isopropanol". June 23, 2015](#)
- [Green, Michael R. and Joseph Sambrook, "Precipitation of DNA with Isopropanol". doi:10.1101/pdb.prot093385 Cold Spring Harb Protoc 2017.](#)
- [Qiagen, "How can I precipitate genomic DNA using isopropanol?".](#)

DOI

dx.doi.org/10.17504/protocols.io.bxg7pjzn

Lynn Doran 2021. Isopropanol DNA Precipitation - Best for DNA Concentration.

protocols.io<https://dx.doi.org/10.17504/protocols.io.bxg7pjzn>

Realizing Increased Photosynthetic Efficiency (RIPE) that is funded by the Bill & Melinda Gates Foundation, Foundation for Food and Agriculture Research, and the U.K. Foreign, Commonwealth & Development Office

Grant ID: OPP1172157

DNA Precipitation, DNA, Isopropanol, DNA Concentration

_____ protocol ,

Aug 18, 2021

Dec 23, 2021

52479

DNA precipitation performed at room temperature minimizes co-precipitation of salt that interferes with downstream applications.

Reagents:

- Isopropanol, room temperature
- 70% Ethanol,

*Measure 70 ml of ethanol and 30 ml of distilled water separately and then mix. Do not measure 70 ml of ethanol and then bring to volume.

NOTE: 70% ethanol is hygroscopic. When opened the ethanol will both evaporate and absorb water over time. Re-use eventually will be at a lower concentration. There is also miscibility of ethanol and water. For example, measuring out 70 mL of ethanol and topping off to 100 mL with water will generate ~65% ethanol. Measuring 70 mL ethanol and 30 mL water separately, then combining them will generate ~95 mL of 70% ethanol.

- Sodium acetate (3 M, pH 5.2): 0.246g/mL (0.738g/3mL dH₂O)
- DNA resuspension solution.

*Molecular Grade Water: **Safest option for downstream applications.**

*TE buffer (pH 8.0), 100 mL: **Allows easier resuspension of DNA pellet.**

Inhibits DNases. Good for long term storage but may interfere with enzymes in downstream applications.

A	B	C
Reagent	Volume	Final concentration
1M Tris-Cl (pH 8.0)	1 mL	10 mM
0.5M EDTA (pH 8.0)	0.02 mL (20 ul)	0.1 mM
Distilled H ₂ O	99 mL	

1) Mix all components. Autoclave to sterilize.

A) 1 M Tris-Cl (pH 8.0), 10 mL

A	B	C
Reagent	Volume	Final concentration
Tris base	1.21 g	10 M
Distilled H ₂ O	10 mL	

- i) Weigh out 1.21 g Tris and add to a 15 mL centrifuge tube.
- ii) Measure out 8 mL of distilled water and add to the tube.
- iii) Stir to mix the solution.
- iv) Add a pH meter into the solution to observe the pH.
- v) Slowly add concentrated hydrochloric acid (HCl) solution using a Pasteur pipette to reduce the pH to 8.0, or another desired pH. Be careful not to add too much at a time, since the pH will change rapidly.
- vi) Once the desired pH has been reached, top up the solution to 10 mL using distilled water.

B) 0.5M EDTA pH 8.0, 5 mL

A	B	C
Reagent	Volume	Final concentration
EDTA disodium salt, dihydrate	0.93 g	0.5 M
Distilled H ₂ O	5 mL	

- i) Weigh out 0.93 g EDTA disodium salt, dihydrate and add to a centrifuge tube.
- ii) Measure out 3.5 mL distilled water and add to the tube.
- iii) Shake to mix. The EDTA salt will not go into solution until the pH reaches 8.0.
- iv) Add a pH meter into the solution to observe the pH.
- v) To dissolve the salt, add sodium hydroxide (NaOH) pellets to the solution. Add a few pellets at a time and wait until the pellets have fully dissolved before adding more. It may take around 0.1 g of NaOH pellets before the pH is at 8.0.
- vi) Once fully dissolved (this will take some time so be patient), top up the solution to 1 mL using distilled water, if necessary.

Materials:

- Analytical balance
- pH meter
- Graduated cylinder
- Pipette, single channel, 1-100 µl
- Pipette tips, 1-100 µl
- 1.5 mL microcentrifuge tube
- Microcentrifuge
- Paper towel or kim wipe

- 1 Transfer desired amount of aqueous DNA solution to a fresh microtube. (70 ul)
- 2 Add 3.0 M pH 5.2 sodium acetate to the DNA solution to result in a final concentration of 0.3 M sodium acetate. (7 ul)
- 3 Add 0.6–0.7 volume of isopropanol at room temperature and mix well. (50 ul)
- 4 Point caps of tubes out to help identify pellet location below cap. Centrifuge the sample immediately at 10,000–15,000 x g for 20–30 min in a microcentrifuge tube at 4°C.
- 5 Identify the pelleted DNA. Carefully decant the supernatant fluid into a fresh labeled tube without losing the pellet.

If concerned about losing DNA pellet during decanting, the supernatant can be stored in a fresh tube until the recovery of the precipitated DNA is confirmed.

- 6 Wash the DNA pellet by adding 500 ul or more (depending on the size of the preparation) of room-temperature 70% ethanol. Flick tube and invert several times until pellet is removed from the side of the tube and all sides can be washed.

Thorough 70% ethanol washing is critical to remove residual isopropanol which can hinder the redissolution of the DNA.

- 7 Centrifuge the sample at 10,000–15,000 x g for 20–30 min at 4°C.

- 8 Visually identify the pellet. Pour out ethanol. Tap tube gently on paper towel, being sure to keep track of the pellet.





Optional: For high quality DNA applications requiring very clean DNA extractions, repeat steps 6-8. More ethanol washes will increase quality but will negatively affect recovery.

- 10 After last ethanol wash, use a pipette to remove any ethanol that remains after decanting.

- 11 Allow samples to air dry at room temperature for 10-20 min. When pellet turns from white to clear the pellet is dry. Do not over dry.

Overdrying of the pellet will make resuspension of the DNA more difficult.

- 12 Resuspend the DNA pellet in Tris, TE (pH 8.0), or molecular grade water. Incubate in  4 °C overnight for more complete resuspension. For long term storage, store at  -20 °C .