



Version 2 ▼

Jul 19, 2022

Preparation of M9 worm buffer V.2

Alfonso Pérez Escudero¹, Alid Al-Asmar¹, Gabriel Madirolas¹¹CNRS, Université Paul Sabatier (Toulouse III)

1 Works for me



Share

dx.doi.org/10.17504/protocols.io.x54v9yo6pg3e/v2

Perez-Escudero Lab



Alfonso Pérez Escudero

ABSTRACT

Preparation of M9 worm buffer, used to maintain *Caenorhabditis elegans* in liquid for short periods of time. Protocol taken from the WormBook Methods

http://www.wormbook.org/chapters/www_strainmaintain/strainmaintain.html

DOI

dx.doi.org/10.17504/protocols.io.x54v9yo6pg3e/v2

PROTOCOL CITATION

Alfonso Pérez Escudero, Alid Al-Asmar, Gabriel Madirolas 2022. Preparation of M9 worm buffer. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.x54v9yo6pg3e/v2>

Version created by [theapelab](#)



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Jul 04, 2022

LAST MODIFIED

Jul 19, 2022

OWNERSHIP HISTORY

Jul 19, 2022



Alfonso Pérez Escudero

MATERIALS TEXT

Our references (also check our magnesium sulfate solution protocol):

Monopotassium phosphate (KH_2PO_4):

[Potassium phosphate monobasic](#) **Sigma**

Aldrich Catalog #795488-500G






Disodium phosphate dihydrate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$):

[Sodium phosphate dibasic](#)

dihydrate Sigma Catalog #71643-250G

[Sodium](#)

Sodium chloride (NaCl): [chloride](#) **Sigma Catalog #S5886-1KG**

- 1 Add  **3 g ±0.03** of monopotassium phosphate (KH_2PO_4) in a 1L clean bottle.
- 2 Add  **7.52 g ±0.05** of disodium phosphate dihydrate ($\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$).
- 3 Add  **5 g ±0.05** of sodium chloride (NaCl).
- 4 Add  **1 L** of milliQ water.
- 5 Autoclave.
- 6 Add  **1 mL** of magnesium sulfate solution. Our protocol for making this solution can be found here:



Preparation of 1M magnesium sulfate solution (MgSO₄)
by theapelab

PREVIEW

RUN

