



2 ▼

Dec 17, 2021

# Isolation of nuclei from frozen human skeletal and cardiac muscle for single nucleus RNA and chromatin assays V.2

James R Hocker<sup>1</sup>, Sebastian Preissl<sup>1</sup>, Dinh H Diep<sup>1</sup><sup>1</sup>University of California, San Diego

1

[dx.doi.org/10.17504/protocols.io.b2z4qf8w](https://dx.doi.org/10.17504/protocols.io.b2z4qf8w)

Dinh Diep  
University of California, San Diego

This protocol was adapted for the isolation of single nuclei from frozen skeletal and cardiac muscle tissues for molecular characterization with the SNARE-seq2, sci-ATAC-seq, and snRNA-seq assays.

## References:

- (1) Preissl et al (2015). Circulation Research. Doi: 10.1161/CIRCRESAHA.115.306337
- (2) Hocker et al (2021). Science Advances. Doi: 10.1126/sciadv.abf1444

DOI

[dx.doi.org/10.17504/protocols.io.b2z4qf8w](https://dx.doi.org/10.17504/protocols.io.b2z4qf8w)

James R Hocker, Sebastian Preissl, Dinh H Diep 2021. Isolation of nuclei from frozen human skeletal and cardiac muscle for single nucleus RNA and chromatin assays . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.b2z4qf8w>

Dinh Diep



protocol

Hocker, J. D., Poirion, O. B., Zhu, F., Buchanan, J., Zhang, K., Chiou, J., Wang, T.-M., Zhang, Q., Hou, X., Li, Y. E., Zhang, Y., Farah, E. N., Wang, A., McCulloch, A. D., Gaulton, K. J., Ren, B., Chi, N. C., & Preissl, S. (2021). Cardiac cell type-specific gene regulatory programs and disease risk association. Science Advances, 7(20), eabf1444. <https://doi.org/10.1126/sciadv.abf1444>

protocol ,

Dec 17, 2021

Dec 17, 2021

56092

Required consumables per sample:

1. gentleMACS M Tube (Miltenyi Biotec, 130-096-335)
2. CellTrics filter, 30 uM (CellTrics, 04-0042-2316)
3. Eppendorf Tubes, 5 mL

Required instrument:

1. gentleMACS Tissue Dissociator (Miltenyi Biotec, 130-095-937)

Buffers:

A	B	C	D	E	F
Reagent	Stock conc.	Final conc.	1 mL	5 mL	10 mL
PBS	1X	1X	0.925 mL	4.625 mL	9.25 mL
BSA (Sigma Aldrich)	n/a	5%	50 mg	250 mg	500 mg
IGEPAL CA-630 / NP 40	10%	0.20%	20 uL	100 uL	200 uL
DTT	100 mM	1 mM	10 uL	50 uL	100 uL
Protease Inhibitor	25X	1X	40 uL	200 uL	400 uL
Enzymatics RNase In	40 U/uL	0.2 U/uL	5 uL	25 uL	50 uL

**Nuclear Permeabilization Buffer (NPB). Make fresh for each use.**

A	B	C	D	E	F
Reagent	Stock conc.	Final conc.	1 mL	4 mL	12 mL
Nuclease free water	n/a	n/a	0.927 mL	3.708 mL	11.124 mL
CaCl <sub>2</sub>	1 M	5 mM	5 uL	20 uL	60 uL
MgOAc	1 M	3 mM	3 uL	12 uL	36 uL
Tris-HCl, pH 8.0	1 M	10 mM	10 uL	40 uL	120 uL
EDTA	500 mM	2 mM	4 uL	16 uL	48 uL
DTT	100 mM	0.6 mM	6 uL	24 uL	72 uL
Protease Inhibitor	25X	1X	40 uL	160 uL	480 uL
Enzymatics RNase In	40 U/uL	0.2 U/uL	5 uL	20 uL	60 uL

**MACS Buffer. Make a stock without Protease In and add when using.**

A	B	C	D	E	F
Reagent	Stock conc.	Final conc.	1 mL	5 mL	10 mL
PBS	1X	1X	0.93 mL	4.65 mL	9.3 mL
Suprase In	20 U/uL	0.05 U/uL	2.5 uL	12.5 uL	25 uL
Enzymatics RNase In	40 U/uL	0.05 U/uL	1.25 uL	6.75 uL	12.5 uL

**PBS + Rnase In. Make fresh for each use.**

- 1 Section flash-frozen skeletal muscle into aliquots according to desired nuclei yield and store on dry ice or at **-80 °C**.

Expect yield of 2000-4000 nuclei per mg of skeletal muscle tissue.

- 2 Prepare buffers and chill components:

A	B	C
Tissue mass	MACS buffer	NPB buffer
10-50mg	1 mL	2 mL
50-100mg	2 mL	4 mL
>100mg	3 mL	1 mL per 25 mg

Note: May require further optimization.

- Prepare buffer fresh on day of nuclei isolation
- Precool centrifuge to **4 °C**
- Precool all buffers at **4 °C**
- Place gentleMACS dissociator in cold room / chiller
- Chill gentleMACS M tubes on ice
- Chill 5.0 mL Eppendorf tubes on ice

3 Transfer sectioned frozen tissue to gentleMACS M tubes.

4 Add the recommended volume of MACS buffer to gentleMACS M tube.

Allow tissue to thaw in MACS buffer for approximately 60 seconds on ice.

5 Homogenize with Miltenyi tissue dissociation protocol: "Protein\_01\_01" on the gentleMACS instrument in the chiller/cold room.

6 Briefly centrifuge the gentleMACS M tube to pull all the homogenate to the bottom. 15s  
**160 x g, 4°C, 00:00:15**

7 Filter homogenate through 30 µm CellTrics filter into 5 mL Eppendorf Tubes.

8 Wash gentleMACS M tube with another 1 mL of MACS buffer and filter the wash.

9 **900 x g, 4°C, 00:10:00** 10m

Use ramp rate: 3/9 acceleration and 3/9 deceleration.

10 Decant and discard the supernatant.

11 Resuspend the pellet with the recommended NPB buffer.


12 Gently rotate the sample in cold room/chiller for  **00:10:00** . 10m

Users should optimize lysis timing for different samples. (i.e., 5 minutes for cardiac tissues).

13 Centrifuge permeabilized nuclei at  **900 x g, 4°C, 00:10:00** . 10m

14 Decant and discard the supernatant.


15 Resuspend the pellet in 500 uL PBS + RNase In.


For SNARE-seq2: immediately fix nuclei by adding (1:1) 500 uL PBS + Rnase In + 1% paraformaldehyde. (Make 1 mL : 937.5 uL PBS+RI + 62.5 uL 16% methanol-free PFA). Incubate on ice for  **00:10:00**

16 **QA/QC:** Count the nuclei. Mix 1:1 volume of cellular suspension with a staining solution (ie DAPI). Load 10 uL of mixture onto a Biorad Cell Counter slide then count with a BioRad TC20 Cell Counter. Gate 4 uM-6 uM for nuclei sizes.

**QA/QC:** Check nuclei integrity under fluorescent microscope using DAPI channel. Nuclei

17 should appear distinct, have rounded borders and the majority occurring as singlets.

If high clumping is observed: bring total volume to ~1 mL with PBS+0.1%BSA and filter the sample through a 30 uM CellTrics filter. Pellet  **900 x g, 4°C, 00:10:00** . Then resuspend again in 100 uL PBS + Rnase In.

If high debris is observed (low DAPI+): bring volume to ~ 1 mL with PBS+0.1%BSA and pellet:  **900 x g, 4°C, 00:10:00** . Then resuspend again in 100 uL PBS + RNase In.