

Jun 18, 2024

🌐 Gene Expression Dual Index Library Construction

📁 In 1 collection

DOI

dx.doi.org/10.17504/protocols.io.n92ld8r7xv5b/v1



Heidi Monroe¹, Nayra Cardenes², Melanie Königshoff², koenigshoffm², Robert Lafyatis¹

¹University of Pittsburgh;

²Division of Pulmonary, Allergy, Critical Care and Sleep Medicine, School of Medicine, University of Pittsburgh

TriState SenNet

Cellular Senescence Net...

1 more workspace



Nayra Cardenes

Division of Pulmonary, Allergy, Critical Care and Sleep Med...

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.n92ld8r7xv5b/v1

Protocol Citation: Heidi Monroe, Nayra Cardenes, Melanie Königshoff, koenigshoffm, Robert Lafyatis 2024. Gene Expression Dual Index Library Construction. [protocols.io https://dx.doi.org/10.17504/protocols.io.n92ld8r7xv5b/v1](https://dx.doi.org/10.17504/protocols.io.n92ld8r7xv5b/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

Created: June 12, 2024

Last Modified: June 18, 2024

Protocol Integer ID: 101659

Keywords: Fragmentation, End Repair, A-tailing, Library construction, Adaptor ligation, Sample index PCR, SenNet, TriState, Lung, PCLS, SenNet, TriState

Funders Acknowledgement:
TriState SenNET (Lung and Heart) Tissue Map and Atlas consortium - NIA
Grant ID: U54AG075931

Abstract

The Chromium Single Cell Gene Expression Solution upgrades short read sequencers to deliver a scalable microfluidic platform for 3' digital gene expression by profiling 500-10,000 individual cells per sample.

Once cDNA is amplified, enzymatic fragmentation and size selection are used to optimize the cDNA amplicon size. P5, P7, i7 and i5 sample indexes, and TruSeq Read 2 (read 2 primer sequence) are added via End Repair, A-tailing, Adaptor Ligation, and PCR. The final libraries contain the P5 and P7 primers used in Illumina amplification.

This protocol details the gene expression dual index library construction and sequencing.

Attachments



snRNAseq_ProtocolsIO..

66KB

Image Attribution

Nayra Cardenes, PhD

Materials

- Fragmentation Buffer (2000091)
- Adaptor Oligos (2000094)
- Ligation Buffer (2000092)
- Dual Index Plate TT Set A (3000431)
- Fragmentation Enzyme (2000090/2000104)
- DNA Ligase (220110/220131)
- Amp Mix (2000047/2000103)
- Cleanup Buffer (2000088)

Fragmentation Mix:

A	B	C	D
Reagents	1X (μl)	4X+10% (μl)	8X+10% (μl)
Fragmentation Buffer	5	22	44
Fragmentation Enzyme	10	44	88
Total	15	66	132

Ligation Mix:

A	B	C	D
Reagents	1X (μl)	4X+10% (μl)	8X+10% (μl)
Ligation Buffer	20	88	176
DNA Ligase	10	44	88
Adaptor Oligos	20	88	179
Total	50	220	440

Safety warnings

- !** The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to **protocols.io** is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with **protocols.io**, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.



Before start

- **Equilibrate to Room temperature (RT)** – Fragmentation Buffer (2000091), Adaptor Oligos (2000094), Ligation Buffer (2000092) and Dual Index Plate TT Set A (3000431)
- **Place on ice**– Fragmentation Enzyme (2000090/2000104), DNA Ligase (220110/220131) and Amp Mix (2000047/2000103)
- **Thaw at 65°C**- Cleanup Buffer (2000088)

Gene Expression Dual Index Library Construction

1w 2d 0h 34m 30s

1 Fragmentation, End Repair & A-tailing

1.1 Prepare a thermal cycler with the following incubation protocol:

A	B	C
Lid Temperature	Reaction Volume	Run Time
65°C	50 µl	~35 min
Step	Temperature	Time
Pre-cool block	4°C	Hold
Fragmentation	32°C	5 min
End Repair & A-tailing	65°C	30 min
Hold	72°C	Hold

Thermocycler protocol.

1.2 Vortex fragmentation buffer. Verify there is no precipitate.






1.3 Prepare fragmentation mix  On ice . Pipette mix and centrifuge briefly.



Fragmentation Mix


A	B	C	D
Reagents	1X (µl)	4X+10% (µl)	8X+10% (µl)
Fragmentation Buffer	5	22	44
Fragmentation Enzyme	10	44	88
Total	15	66	132

Calculation for the Fragmentation Mix preparation.

1.4 Transfer ONLY  10 µL purified cDNA sample from pellet cleanup to a tube strip. The remaining  30 µL (75%) cDNA sample can be stored at  4 °C for up to  72:00:00 or at  -20 °C for up to 4 weeks for generating additional 3' Gene Expression libraries.

3d




1.5 Add  25 μL buffer EB to each sample.



1.6 Add  15 μL fragmentation mix to each sample.




1.7 Pipette mix 15 \times (pipette set to 35 μl)  On ice . Centrifuge briefly.





1.8 Transfer into the pre-cooled thermal cycler ( 4 $^{\circ}\text{C}$) and press "SKIP" to initiate the protocol.

2 Post Fragmentation, End Repair & A-tailing Double Sided Size Selection – SPRIselect:

2.1 Vortex to resuspend SPRIselect reagent. Add  30 μL SPRIselect (0.6X) reagent to each sample. Pipette mix 15 \times (pipette set to 75 μl).




2.2 Incubate  00:05:00 at  Room temperature .



5m



2.3 Place on the magnet. High until the solution clears. DO NOT discard supernatant.

2.4 Transfer  75 μL supernatant to a new tube strip.


2.5 Vortex to resuspend SPRIselect reagent. Add  10 μL SPRIselect reagent (0.8X) to each transferred supernatant. Pipette mix 15 \times (pipette set to 80 μl).



2.6 Incubate  00:05:00 at  Room temperature .

5m



2.7 Place on the magnet. High until the solution clears.

2.8 Remove  80 μL supernatant. DO NOT discard any beads.

2.9 Add  125 µL 80% ethanol to the pellet. Wait  00:00:30 .

30s




2.10 Remove the ethanol.



2.11 **Repeat sub-steps 2.9 and 2.10 for a total of 2 washes.**

2.12 Centrifuge briefly. Place on the magnet. Low until the solution clears. Remove remaining ethanol.



2.13 Remove from the magnet. Add  50.5 µL buffer EB to each sample. Pipette mix 15× (pipette set to 45 µl).




2.14 Incubate  00:02:00 at  Room temperature .

2m



2.15 Place on the magnet. High until the solution clears.

2.16 Transfer  50 µL sample to a new tube strip.



3 Adaptor Ligation

3.1 Prepare adaptor ligation Mix. Pipette mix and centrifuge briefly.



Adaptor Ligation Mix:

A	B	C	D
Reagents	1X (µl)	4X+10% (µl)	8X+10% (µl)
Ligation Buffer	20	88	176
DNA Ligase	10	44	88
Adaptor Oligos	20	88	179
Total	50	220	440

3.2 Add  50 µL adaptor ligation mix to  50 µL sample. Pipette mix 15× (pipette set to 90 µl). Centrifuge briefly.




3.3 Incubate in a thermal cycler with the following protocol:





A	B	C
Lid Temperature	Reaction Volume	Run Time
30°C	100 µl	15 min
Step	Temperature	Time
1	20°C	15 min
2	4°C	Hold

4 Post Ligation Cleanup – SPRIselect:

4.1 Vortex to resuspend SPRIselect reagent. Add  80 µL SPRIselect (0.8X) reagent to each sample. Pipette mix 15× (pipette set to 75 µl).





4.2 Incubate  00:05:00 at  Room temperature .

5m



4.3 Place on the magnet. High until the solution clears.

4.4 Remove the supernatant.

4.5 Add  200 µL 80% ethanol to the pellet. Wait  00:00:30 .

30s




4.6 Remove the ethanol.


4.7 Repeat sub-steps 4.5 and 4.6 for a total of 2 washes.



4.8 Centrifuge briefly. Place on the magnet. Low until the solution clears.



4.9 Remove remaining ethanol. Air dry for  00:02:00 .

2m

4.10 Remove from the magnet. Add  30.5 μ L buffer EB to each sample. Pipette mix 15x.

4.11 Incubate  00:02:00 at  Room temperature .

2m





4.12 Place on the magnet. Low until the solution clears.


4.13 Transfer  30 μ L sample to a new tube strip.

5 Sample Index PCR

5.1 Choose the appropriate sample index sets to ensure that no sample indices overlap in a multiplexed sequencing run. Record the 10 \times Sample Index name (PN-3000431 Dual Index Plate TT Set A well ID) used.

5.2 Add  50 μ L Amp mix to  30 μ L sample.



5.3 Add  20 μ L of an individual Dual Index TT Set A to each sample and record the well ID used. Pipette mix 5 \times (pipette set to 90 μ L). Centrifuge briefly.



5.4 Incubate in a thermal cycler with the following protocol:



A	B	C
Lid Temperature	Reaction Volume	Run Time
105°C	100 μ L	~25-40 min
Step	Temperature	Time
1	98°C	45 sec



A	B	C
2	98°C	20 sec
3	54°C	30 sec
4	72°C	20 sec
5	Go to Step 2 – * # cycles calculated below	
6	72°C	1 min
7	4°C	Hold

Note

The total cycles should be optimized based on 25% carry forward cDNA yield/input calculated during cDNA QC & Quantification

A	B
cDNA Input	Total cycles
0.25-25 ng	14-16
25-150 ng	Dec-14
150-500 ng	10-Dec
500-1,000 ng	08-Oct
1,000-1,500 ng	06-Aug
>1,500 ng	5

5.5 Store at 4 °C for up to 72:00:00 or proceed to the next step.

3d

6 Post Sample Index PCR Double Sided Size Selection – SPRIselect

6.1 Vortex to resuspend SPRIselect reagent. Add 60 µL SPRIselect (0.6X) reagent to each sample. Pipette mix 15× (pipette set to 150 µl).

6.2 Incubate 00:05:00 at Room temperature .


5m







6.3 Place on the magnet. High until the solution clears. DO NOT discard supernatant.

6.4 Transfer  150 μL supernatant to a new tube strip.


6.5 Vortex to resuspend SPRIselect reagent. Add  20 μL SPRIselect reagent (0.8X) to each transferred supernatant. Pipette mix 15 \times (pipette set to 150 μL).



6.6 Incubate  00:05:00 at  Room temperature .

5m



6.7 Place on the magnet. High until the solution clears.

6.8 Remove  165 μL supernatant. DO NOT discard any beads.

6.9 Add  200 μL 80% ethanol to the pellet. Wait  00:00:30 .

30s




6.10 Remove the ethanol.



6.11 **Repeat sub-steps 6.9 and 6.10 for a total of 2 washes.**

6.12 Centrifuge briefly. Place on the magnet. Low until the solution clears. Remove remaining ethanol.



6.13 Remove from the magnet. Add  35.5 μL buffer EB to each sample. Pipette mix 15 \times (pipette set to 35 μL).



6.14 Incubate  00:02:00 at  Room temperature .




2m






6.15 Place on the magnet. Low until the solution clears.

6.16 Transfer  35 µL sample to a new tube strip.

6.17 Store at  4 °C for up to  72:00:00 or at  -20 °C for long-term storage.

3d

7 Post Library Construction QC

7.1 Run  1 µL sample at 1:10 dilution on an Agilent Bioanalyzer High Sensitivity chip.

7.2 Determine the average fragment size from the Bioanalyzer trace. This will be used as the insert size for library quantification.

Sequencing

8 3' Gene Expression Library Sequencing Depth & Run Parameters:

A	B
Sequencing Depth	Minimum 20,000 read pairs per cell
Sequencing Type	Pair-end, dual indexing
Sequencing Read	Recommended Number of cycles
Read 1	28 cycles
i7 Index	10 cycles
i5 Index	10 cycles
Read 2	90 cycles

9 Once quantified and normalized, the 3' Gene Expression libraries should be denatured and diluted as recommended for Illumina sequencing platforms.



Protocol references

https://cdn.10xgenomics.com/image/upload/v1660261285/support-documents/CG000505_Chromium_Nuclei_Isolation_Kit_UG_RevA.pdf