



1 ▼

Jan 23, 2022

# 🌐 General proteomics FASP (Filter-Aided Sample Preparation) V.1

Jacob Waldbauer<sup>1</sup>, Lichun Zhang<sup>1</sup><sup>1</sup>University of Chicago

1

[dx.doi.org/10.17504/protocols.io.if7cbrn](https://dx.doi.org/10.17504/protocols.io.if7cbrn) Jacob Waldbauer

Spin-filter based protein digestion and purification for bottom-up proteomics

DOI

[dx.doi.org/10.17504/protocols.io.if7cbrn](https://dx.doi.org/10.17504/protocols.io.if7cbrn)

Jacob Waldbauer, Lichun Zhang 2022. General proteomics FASP (Filter-Aided Sample Preparation). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.if7cbrn>



protocol ,

Jun 14, 2017

Jan 23, 2022

6367

It is important that nearly all liquid passes through spin filters after each centrifugation step. If a substantial amount of liquid (more than 75  $\mu$ L) remains on filter after 10 minutes of centrifugation, spin for longer.

## MATERIALS

☒ Pierce trypsin Protease Thermo Scientific

☒ Tween 20 Bio-rad

Laboratories Catalog #170-6606-MSDS

☒ Optima LC/MS grade Water Fisher

Scientific Catalog #W6-4

☒ Urea Fisher

Scientific Catalog #U15-500

☒ Ammonium bicarbonate Fisher

Scientific Catalog #A643-500

☒ Deoxycholic

Acid bioworld Catalog #40430004-1

☒ Optima LC/MS grade Acetonitrile Fisher

Scientific Catalog #A955-212

☒ Formic Acid (Optima LC/MS) Fisher

Scientific Catalog #A117-50

☒ Ethyl Acetate Fisher

Scientific Catalog #E145-4

Passivate Vivacon filter unit (Satorius Vivicon 500, 30,000 MWCO) and collection tube overnight in 5% (v/v) Tween-20 in MilliQ water. The next day, rinse with nanopure water until no suds, then rinse 3x with LC-MS water.

## Day 1

- 1 Prepare fresh buffers in LC-MS water:  
Exchange buffer: 8M urea, 0.2% (w/v) deoxycholate, 1M ammonium bicarbonate (pH 8)  
Digestion buffer: 0.2% (w/v) deoxycholate, 50 mM ammonium bicarbonate (pH 8)  
Trypsin dissolution/peptide recovery buffer: 50 mM ammonium bicarbonate (pH 8)
- 2 Dispense 25-50µL (depending on concentration) reduced, alkylated protein extract to filter unit.
- 3 Fill filter units with 200 µL exchange buffer, mix with pipet.
- 4 Spin at 14,000g for 10 min, discard filtrate.

- 5 Repeat steps 3 and 4 three more times.
- 6 Wash filter unit with 200  $\mu$ L digestion buffer, and spin at 14,000g for 10 min. Do this a total of 2 times.
- 7 Transfer filter unit to passivated collection tube.
- 8 Add 110  $\mu$ L digestion buffer and desired amount of trypsin (~1:50 trypsin/protein), dissolved in 50 mM ammonium bicarbonate, keeping trypsin on ice or at 4° C. Incubate at 37° C overnight. Parafilm tubes is recommended to prevent evaporation.

#### Day 2

- 9 Remove parafilm if applicable and centrifuge tubes at 14,000g for 10 min. **Do not discard filtrate!**
- 10 Add 50  $\mu$ L of peptide recovery buffer to filters, spin at 14,000g for 10 min. Do this step a total of 2 times.
- 11 Transfer the filtrate to a LoBind tube.
- 12 Add 900 $\mu$ L ethyl acetate and 2.5  $\mu$ L TFA to filtrate, vortex. Precipitate may be visible.
- 13 Sonicate in a bath sonicator for 10s, and centrifuge at 16,000g for 10 min.
- 14 Carefully remove upper organic layer with a pipet without disturbing phase boundary.

- 15 Repeat steps 12, 13 and 14 two more times (Note: Only add 900µL ethyl acetate, no more TFA)
- 16 Place sample uncovered in Thermomixer and dry at 60° C for 5-10 minutes until ethyl acetate is gone.
- 17 Freeze sample solid at -80° C, Centrивap to dryness.