



plasma preparation_exRNAQC V.7

Anneleen Decock¹

¹Center for Medical Genetics Ghent (CMGG), Cancer Research Institute Ghent (CRIG)

Version 7

Aug 10, 2020

1

Works for me

dx.doi.org/10.17504/protocols.io.n99dh96

CMGG



Anneleen Decock

Center for Medical Genetics Ghent (CMGG), Cancer Research In...

ABSTRACT

This protocol describes how to prepare plasma from venous blood draw.

ATTACHMENTS

[plasma preparation.pdf](#)

DOI

dx.doi.org/10.17504/protocols.io.n99dh96

PROTOCOL CITATION

Anneleen Decock 2020. plasma preparation_exRNAQC. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.n99dh96>

KEYWORDS

plasma

LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Apr 06, 2018

LAST MODIFIED

Aug 10, 2020

PROTOCOL INTEGER ID

11297

GUIDELINES

- In the Center for Medical Genetics Gent (CMGG), plasma preparation is performed in the blood lab. Guidelines on how to handle blood samples and how to work in this lab are described in H5.3-OP2-B1.
- Take pictures (optional) on a white background, upright position

MATERIALS

NAME	CATALOG #	VENDOR
Liquid Nitrogen		
Pipettes		
Safe-Lock cup DNA LoBind 1,5ml PCR clean	A08970	Eppendorf
Cryotube	10674511	Thermo Fisher Scientific

NAME	CATALOG #	VENDOR
Centrifuge with swinging bucket rotor (centrifugation speed up to 2500 g (rcf)) and buckets for blood collection tubes and conical 15 ml tubes		
Disposable lab coat		
Nitrile powder-free gloves		
Nalgene® 62080-00 VERSI-DRY® Lab Soaker Bench Protector Mat, White Standard Absorbency	NAL-62080-00	Thermo Fisher Scientific
Filter tips		
Racks for the different types of tubes		
Conical tube, 15 ml	188271	greiner bio-one

DISCLAIMER:

DISCLAIMER – FOR INFORMATIONAL PURPOSES ONLY; USE AT YOUR OWN RISK

The protocol content here is for informational purposes only and does not constitute legal, medical, clinical, or safety advice, or otherwise; content added to protocols.io is not peer reviewed and may not have undergone a formal approval of any kind. Information presented in this protocol should not substitute for independent professional judgment, advice, diagnosis, or treatment. Any action you take or refrain from taking using or relying upon the information presented here is strictly at your own risk. You agree that neither the Company nor any of the authors, contributors, administrators, or anyone else associated with protocols.io, can be held responsible for your use of the information contained in or linked to this protocol or any of our Sites/Apps and Services.

Centrifugation step 1: 20 min at 400 g (rcf)

- 1
 - Invert tubes 5 times before centrifugation
 - Spin tubes for 20 min at 400 g (rcf) (without brake), at room temperature. Note time point of start of centrifugation
 - Pipette platelet-rich plasma (PRP) carefully into a new collection tube, leave ± 0.5 cm above the buffy coat (do not disturb the buffy coat)
 - Invert the PRP tube before aliquoting. Aliquot the PRP into cryovials or LoBind tubes, snap freeze in liquid nitrogen and store at -80 °C (note time point of snap freeze), and/or continue to prepare platelet-poor plasma (PPP)

Centrifugation step 2: 10 min at 800 g (rcf)

- 2
 - Spin the PRP for 10 min at 800 g (rcf) (without brake) to obtain platelet-poor plasma (PPP), at room temperature
 - Pipette PPP carefully into a new collection tube, leave ± 0.5 cm above pellet (do not disturb pellet)
 - *Optional: If platelets need to be collected, resuspend pellet in remaining volume above pellet, transfer to 1.5 ml tube, and spin again for 10 min at 800 g (rcf) (without brake), remove volume and snap freeze pellet in 1.5 ml tube. Note time point of snap freeze*
 - Invert the PPP tube before aliquoting. Aliquot the PPP into cryovials or LoBind tubes, snap freeze in liquid nitrogen and store at -80°C (note time point of snap freeze), and/or continue to prepare platelet-free plasma (PFP)

Centrifugation step 3: 15 min at 2500 g (rcf)

- 3
 - Spin the PPP for 15 min at 2500 g (rcf) (without brake) to obtain platelet-free plasma (PFP), at room temperature
 - Pipette PFP carefully into a new collection tube, leave ± 0.5 cm above pellet (do not disturb pellet)
 - Aliquot the PFP into cryovials or LoBind tubes, snap freeze in liquid nitrogen and store at -80°C (note time point of snap freeze)