



Upload image

Version 2

Aug 08, 2020

# Cytotoxicity assay using LLC-MK2 cell line V.2

luanaborba<sup>1</sup><sup>1</sup>Universidade Federal do Rio de Janeiro

In Development

dx.doi.org/10.17504/protocols.io.bjgbkjsn

L. P. Borba-Santos

luanaborba

## ABSTRACT

[Protocol that described the cytotoxicity assay](#) using a mammalian epithelial cell line (LLC-MK2; ATCC CCL-7).

## DOI

dx.doi.org/10.17504/protocols.io.bjgbkjsn

## DOCUMENT CITATION

luanaborba 2020. Cytotoxicity assay using LLC-MK2 cell line. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.bjgbkjsn>

## LICENSE

This is an open access document distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Aug 07, 2020

## LAST MODIFIED

Aug 08, 2020

## DOCUMENT INTEGER ID

40163

## Cytotoxicity assay using LLC-MK<sub>2</sub> cell line

1. The mammalian epithelial cell line LLC-MK<sub>2</sub> (ATCC CCL-7) was cultivated in RPMI 1640 medium<sup>1</sup> supplemented with 2 mM L-glutamine and heat-inactivated 10% fetal bovine serum and buffered with sodium bicarbonate at pH 7.5;
2. Aliquots of 100 µl containing 2.5 x 10<sup>5</sup> cells/ml were added into flat-bottom 96-well microplates (TPP™) and incubated for 48 h at 37 °C, in a 5% CO<sub>2</sub> atmosphere, to obtain confluent monolayers of LLC-MK<sub>2</sub> cells;
3. The supernatant was gently discarded, samples were washed in sterile PBS, and 100 µl of compounds previously diluted<sup>2</sup> in RPMI 1640 medium were added on confluent monolayers;
4. Cells were treated for 48 h at 37 °C, in a 5% CO<sub>2</sub> atmosphere;
5. The supernatant was gently discarded and samples were washed in sterile PBS;
6. Cell viability was analyzed by the tetrazolium (XTT) reduction assay and 150 µl of XTT solution (1 mg/ml XTT<sup>1</sup> and 1mM menadione<sup>1</sup> in PBS) were added in each well;
6. Microplates were incubated for 2 h at 37°C, in a 5% CO<sub>2</sub> atmosphere, in the dark;

7. Microplates were centrifugated at 4000 rpm for 5 min and the supernatant was added into a new microplate;
8. Spectrophotometric readings at 492 nm were performed using a microtiter plate reader<sup>3</sup>;
9. The absorbance value for each well was subtracted from the value for the negative controls<sup>4</sup> and inhibition of cell growth (I) relative to positive controls<sup>5</sup> was calculated according to the following equation:  $I = 100 - (A \times 100/C)$ , where A is the absorbance of treated wells, and C is the absorbance of untreated control wells;
10. The concentration of compounds that elicited 50% cytotoxicity (CC<sub>50</sub>) was estimated by linear regression;
11. The diluent control containing 1% DMSO was included in experiments;
12. Experiments were performed in triplicate in two independent moments.

<sup>1</sup> Sigma Chemical Co., USA.

<sup>2</sup> Stock solutions of compounds in dimethyl sulfoxide (DMSO) at 1 mM were diluted in RPMI 1640 medium supplemented with 2 mM L-glutamine to obtain concentrations of 0.1, 1, 2, 4, 5, 7, 8, 9, and 10 µM.

<sup>3</sup> EMax Plus, Molecular Devices, USA.

<sup>4</sup> Wells without cells containing only RPMI media that were incubated in the same conditions of wells with cells.

<sup>5</sup> Wells with untreated cells.