



Jan 21, 2022

diDO-IPTL protocol

diDO-IPTL protocol

Jacob Waldbauer¹, Lichun Zhang¹¹University of Chicago

1

dx.doi.org/10.17504/protocols.io.b36bqran

Jacob Waldbauer

Isobaric peptide termini labeling (IPTL) is a quantification method which permits relative quantification using quantification points distributed throughout the whole tandem mass spectrometry (MS/MS) spectrum. It is based on the complementary derivatization of peptide termini with different isotopes resulting in isobaric peptides.

DOI

dx.doi.org/10.17504/protocols.io.b36bqran

Jacob Waldbauer, Lichun Zhang 2022. diDO-IPTL protocol. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.b36bqran>



diDO-IPTL protocol, Lichun Zhang

protocol ,

Jan 21, 2022

Jan 21, 2022

57251

DAY 1

- 1 Prepare dried MS-grade trypsin, 2ug per sample to be labeled. Note that separate aliquots of dried trypsin must be used for ¹⁶O- and ¹⁸O-labeling.
- 2 Prepare two separate N-methylmorpholine (NMM) and acetic acid O-isotope exchange buffers with ¹⁶O and ¹⁸O (>98%) enriched water. Volume ratio of buffer components is 48.5 H₂O : 1 NMM : 0.5 acetic acid (pH 7.4). The total volume of buffer to prepare depends on how many

peptide samples are to be labeled. The protocol below is for 40µL of O-isotope exchange buffer per sample (~10-20 µg of peptide). For samples with small amounts of peptide, the concentration can be raised by adding only 20µL buffer per sample and halving the volumes of all subsequent reagent additions.

- 3 Redissolve dried trypsin aliquots in their respective O-isotope exchange buffers and transfer 40µl to each tube of dried peptide sample. KEEP TRACK of which tubes received which buffer.
- 4 Parafilm the tops of the tubes and incubate at 37°C overnight.

DAY 2

- 5 After overnight incubation, add 2µl of 11.3M monochloroacetic acid (prepared in LC-MS grade water) to each peptide tube to lower pH down to 2.6.
- 6 Prepare 16% CH₂O and CD₂O formaldehyde solutions, at least 2µl per sample, diluting stocks with LC-MS grade water if necessary.
- 7 Prepare 4.8M NaBH₃CN in LC-MS water. Tare 1.5mL tube, add a small amount of NaBH₃CN powder in fume hood, weigh, dissolve.
- 8 Add 2µL of CH₂O to each ¹⁸O tube and 2µL of CD₂O to each ¹⁶O tube.
- 9 Add 2µL 4.8M NaBH₃CN to each tube.
- 10 Mix well and incubate at 45°C for 1 hr. While waiting, prepare 5M ammonium formate in LC-MS water.
- 11 Add 2µL 5M ammonium formate to each tube and mix well.

- 12 Add 8 μ L of formic acid to each tube and mix well. Total sample volume is \sim 56 μ L. ^{16}O - and ^{18}O -labeled samples are ready to be mixed and analyzed by LC-MS.