



# Plasmid DNA extraction V.2

Shuning Guo<sup>1</sup>

<sup>1</sup>2021 iDEC NEFU\_China

Version 2 ▾

Oct 03, 2021

1 Works for me

Share

[dx.doi.org/10.17504/protocols.io.byqvpvw6](https://dx.doi.org/10.17504/protocols.io.byqvpvw6)

2021 iDEC NEFU\_China

Shuning Guo

## ABSTRACT

This protocol is used to extract plasmid DNA from E. coli.

## DOI

[dx.doi.org/10.17504/protocols.io.byqvpvw6](https://dx.doi.org/10.17504/protocols.io.byqvpvw6)

## PROTOCOL CITATION

Shuning Guo 2021. Plasmid DNA extraction. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.byqvpvw6>  
Version created by Shuning Guo

## KEYWORDS

Pasmid

## LICENSE

This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Oct 03, 2021

## LAST MODIFIED

Oct 03, 2021

## PROTOCOL INTEGER ID

53749

## GUIDELINES

This protocol is based on E.Z.N.A.® Plasmid DNA Mini Kit I Protocol - Spin Protocol.

## MATERIALS TEXT

- 100% ethanol
- Isopropanol
- Microcentrifuge capable of at least 13,000 x g
- Nuclease-free 1.5 mL or 2 mL microcentrifuge tubes
- Culture tubes
- Optional: sterile deionized water
- Optional: water bath or incubator capable of 70°C
- Optional: 3M NaOH solution

## SAFETY WARNINGS


Please wear gloves before start.

#### BEFORE STARTING

Prepare DNA Wash Buffer, HBC Buffer, and Solution I .

1. Add the vial of RNase A to the bottle of Solution I and store at 2-8°C;
2. Dilute DNA Wash Buffer with 100% ethanol 120ml and store at room temperature;
3. Dilute HBC Buffer with isopropanol and store at room temperature;
4. Check Solution II and Solution III for precipitation before use. Redissolve any precipitation by warming to 37°C.

- 1 Solate a single colony from a freshly streaked selective plate, and inoculate a culture of 1- 5 mL LB medium containing the appropriate selective antibiotic. Incubate for ~12-16 hours at 37°C with vigorous shaking (~ 300 rpm). Use a 10-20 mL culture tube or a flask with a volume of at least 4 times the volume of the culture. It is strongly recommended that an endA negative strain of E. coli be used for routine plasmid isolation. Examples of such strains include DH5a® and JM109®.

- 2 Centrifuge at 10,000 x g for 1 minute at room temperature. 1m  
 **10000 x g, Room temperature , 00:01:00**

- 3 Decant or aspirate and discard the culture media.

- 4 Add 250 µL Solution I/RNase A. Vortex or pipet up and down to mix thoroughly. Complete resuspension of cell pellet is vital for obtaining good yields.

RNase A must be added to Solution I before use.

- 5 Transfer suspension into a new 1.5 mL microcentrifuge tube.

- 6 Add 250 µL Solution II. Invert and gently rotate the tube several times to obtain a clear lysate. A 2-3 minute incubation may be necessary.

Avoid vigorous mixing as this will shear chromosomal DNA and lower plasmid purity. Do not allow the lysis reaction to proceed more than 5 minutes. Store Solution II tightly capped when not in use to avoid acidification from CO<sub>2</sub> in the air.

- 7 Add 350 µL Solution III. Immediately invert several times until a flocculent white precipitate forms.

It is vital that the solution is mixed thoroughly and immediately after the addition of Solution III to avoid localized precipitation.

- 8 Centrifuge at maximum speed ( $\geq 13,000 \times g$ ) for 10 minutes. A compact white pellet 10m

will form. Promptly proceed to the next step.

🌀 **15000 x g, Room temperature , 00:10:00**

- 9 Insert a HiBind® DNA Mini Column into a 2 mL Collection Tube.
- 10 Transfer the cleared supernatant from Step 8 by CAREFULLY aspirating it into the HiBind® DNA Mini Column. Be careful not to disturb the pellet and that no cellular debris is transferred to the HiBind® DNA Mini Column.
- 11 Centrifuge at maximum speed for 1 minute. 1m  
🌀 **15000 x g, Room temperature , 00:01:00**
- 12 Discard the filtrate and reuse the collection tube.
- 13 Add 500 µL HBC Buffer

HBC Buffer must be diluted with isopropanol before use. Please see Page 6 for instructions.
- 14 Centrifuge at maximum speed for 1 minute. 1m  
🌀 **15000 x g, Room temperature , 00:01:00**
- 15 Discard the filtrate and reuse collection tube.
- 16 Add 700 µL DNA Wash Buffer.
- 17 Centrifuge at maximum speed for 1 minute. 1m  
🌀 **15000 x g, Room temperature , 00:01:00**
- 18 Discard the filtrate and reuse the collection tube.
- 19 Centrifuge the empty HiBind® DNA Mini Column for 2 minutes at maximum speed to dry the column matrix.

- 20 Transfer the HiBind® DNA Mini Column to a clean 1.5 mL microcentrifuge tube.
- 21 Add 30-100 µL Elution Buffer or sterile deionized water directly to the center of the column membrane.

The efficiency of eluting DNA from the HiBind® DNA Mini Column is dependent on pH. If using sterile deionized water, make sure that the pH is around 8.5.

- 22 Let sit at room temperature for 1 minute.

- 23 Centrifuge at maximum speed for 1 minute.

1m

 **15000 x g, Room temperature , 00:01:00**

This represents approximately 70% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

- 24 Test the concentration and purity of DNA using NanoDrop.

- 25 Store DNA at -20°C.