



Apr 09, 2021

# Untargeted Top-down Proteomics by CZE-MS/MS on Eclipse

Forked from [Untargeted Top-down Proteomics by LC-MS/MS on Eclipse](#)

Kevin Jooß<sup>1</sup>, Bryon Drown<sup>1</sup>, Rafael Melani<sup>1</sup>, Kelleher KRG Research Group<sup>1</sup>

<sup>1</sup>Northwestern University

**1** Works for me [dx.doi.org/10.17504/protocols.io.btxdnpi6](https://dx.doi.org/10.17504/protocols.io.btxdnpi6)

[Human BioMolecular Atlas Program \(HuBMAP\) Method Development Community](#) [Kelleher Research Group](#)



Kelleher Research Group  
Northwestern University, National Resource for Translational...

## ABSTRACT

Describes the CZE-MS/MS data acquisition procedure for top-down proteomics samples using the Thermo Scientific Orbitrap Eclipse Tribrid mass spectrometer

## DOI

[dx.doi.org/10.17504/protocols.io.btxdnpi6](https://dx.doi.org/10.17504/protocols.io.btxdnpi6)

## PROTOCOL CITATION

Kevin Jooß, Bryon Drown, Rafael Melani, Kelleher KRG Research Group 2021. Untargeted Top-down Proteomics by CZE-MS/MS on Eclipse. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.btxdnpi6>

## FORK NOTE

## FORK FROM

Forked from [Untargeted Top-down Proteomics by LC-MS/MS on Eclipse](#), Kelleher Research Group

## LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

## CREATED

Apr 03, 2021

## LAST MODIFIED

Apr 09, 2021

## PROTOCOL INTEGER ID

48837

## MATERIALS TEXT

CE vendor: Sciex

CE model: CESI 8000 Plus

CE capillary: Neutral OptiMS™ Capillary Cartridges (30 µm ID, L = 90 cm), neutral-coated

Background Electrolyte: 3% Acetic Acid

Conductive Liquid: 3% Acetic Acid

Optima Acetic Acid (Fisher Scientific, cat. no. A11310X1AMP)

Optima Water (Fisher Scientific, cat. no. W6-4)

**1** Samples were analyzed on a Thermo Scientific Orbitrap Eclipse Tribrid mass spectrometer in line with a Sciex CESI 8000

Plus

2 Samples were prepared in (  10 µl ) 0.3% acetic acid and transferred to sample vials.

3 Setup buffers in tray:

A	B
Location	Buffer composition
BI:A1	Background electrolyte
BI:B1	Background electrolyte
BI:C1	0.1 M HCl
BI:D1	Water
BO:A1	Conductive liquid
BO:B1	Water

Background electrolyte: 3.0% acetic acid; Conductive liquid: 3.0% acetic acid

4 Prior to first use, the CESI cartridge were conditioned. This process can be divided into two methods: 1) washing method, 2) electrical conditioning method. Each morning, a shortened version of the electrical conditioning method was performed (+15 kV for 30 min). For long term storage, capillaries were rinsed with water and kept at 4 °C.

A	B	C	D	E	F	G
Time [min]	Event	Value	Duration	Inlet	Outlet	Summary
	Rinse - Pressure	100 psi	5.00 min	BI: C1	BO: B1	forward
	Rinse - Pressure	100 psi	10.00 min	BI: A1	BO: B1	forward
	Rinse - Pressure	100 psi	5.00 min	BI: D1	BO: B1	reverse
	Rinse - Pressure	100 psi	30.0 min	BI: D1	BO: B1	forward

Washing method

A	B	C	D	E	F	G
Time [min]	Event	Value	Duration	Inlet	Outlet	Summary
	Rinse - Pressure	100 psi	3.00 min	BI: A1	BO: A1	reverse
	Rinse - Pressure	100 psi	5.00 min	BI: A1	BO: A1	forward
0.00	Separate - Voltage	15.0 kV	60.00 min	BI: B1	BO: A1	1.00 min ramp, normal polarity, both
60.00	Separate - Voltage	1.0 kV	5.00 min	BI: B1	BO: A1	5.00 min ramp, normal polarity, both
65.00	End					

Electrical conditioning method

- 5 Prior each sample injection, the capillary is flushed with 0.1 M HCl, filled with new background electrolyte, and conductive liquid was replaced. Sample is injected hydrodynamically at 2.5 psi for 60 sec (estimated 20 nL injection volume). Separation is performed at 15kV with 0.5 psi supplemental pressure for 60 minutes.

Method for CESI 8000 Plus is summarized below:

A	B	C	D	E	F	G
Time [min]	Event	Value	Duration	Inlet	Outlet	Summary
	Rinse - Pressure	100 psi	5.00 min	BI: C1	BO: A1	forward
	Rinse - Pressure	100 psi	3.00 min	BI: A1	BO: A1	reverse
	Rinse - Pressure	100 psi	5.00 min	BI: A1	BO: A1	forward
	Inject - Pressure	2.5 psi	60 sec	sample vial	BO: A1	forward
	Wait		0.00 min	BI: D1	BO: A1	dipping
	Inject - Pressure	2.5 psi	10 sec	BI: B1	BO: A1	forward
0.00	Separation - Voltage	15 kV 0.5 psi	60.00 min	BI: B1	BO: A1	1.0 min ramp, normal polarity, both
1.00	Relay On					
60.00	Separation - Voltage	1.0 kV 5 psi	5.00 min	BI: B1	BO: A1	5.0 min ramp, normal polarity, both
65.00	End					

CESI 8000 method

- 6 Eluted proteins were ionized in positive ion mode nanoelectrospray ionization (nESI) using a sheathless emitter.

A	B
	<b>High-High</b>
<b>Spray voltage</b>	1400 - 1800
<b>Sweep gas</b>	0
<b>Ion transfer tube temp</b>	320
<b>Application mode</b>	Intact Protein
<b>Pressure mode</b>	Low Pressure
<b>Advanced Peak Determination</b>	True
<b>Default charge state</b>	15
<b>S-lens RF</b>	30
<b>Source fragmentation</b>	15 eV

Global MS parameters

- 7 Precursor (intact protein) spectra were acquired at 120k FTRP.

A	B
	<b>High-High</b>
<b>Detector type</b>	Orbitrap
<b>Resolving power</b>	120000
<b>m/z RP measured</b>	200 m/z
<b>Scan range</b>	450-2000
<b>Mass range</b>	Normal
<b>AGC target</b>	2000000
<b>Normalized AGC target</b>	500%
<b>Max Injection Time</b>	50 ms
<b>Microscans</b>	1
<b>Data type</b>	Profile
<b>Polarity</b>	Positive
<b>Use wide quad isolation</b>	True

Parameters for MS1 acquisition

- 8 The mass spectrometer was operated using a TopN 3 sec data-dependent acquisition mode  
Precursor ions were filtered by intensity, charge state, and dynamic exclusion:

A	B
<b>Intensity minimum</b>	5000
<b>Intensity maximum</b>	1E20
<b>Included charge states</b>	4-60
<b>Include undetermined charge states</b>	False
<b>Dynamic exclusion after n times</b>	1
<b>Dynamic exclusion duration</b>	45 s
<b>Mass tolerance</b>	0.5 m/z
<b>Exclude isotopes</b>	True

Precursor selection filters for DDA

- 9 Ions for fragmentation were isolated and fragmented via higher energy dissociation (HCD):

A	B
	<b>High-High</b>
<b>Detector type</b>	Orbitrap
<b>Isolation mode</b>	Quadrupole
<b>Resolving power</b>	60000
<b>m/z RP measured</b>	200 m/z
<b>Scan range</b>	350-2000
<b>AGC target</b>	400000
<b>Normalized AGC target</b>	800%
<b>Max injection time</b>	600 ms
<b>Microscans</b>	1
<b>Isolation window</b>	3 m/z
<b>Activation type</b>	HCD
<b>Collision energy</b>	32
<b>Collision energy mode</b>	Fixed
<b>Polarity</b>	Positive

Parameters for MS2 acquisition