

DEC 08, 2023

OPEN ACCESS



Protocol Citation: Andrea R
Daniel 2023. TFome
CRISPRko screens.
protocols.io

https://protocols.io/view/tfom e-crisprko-screens-c549y8z6

MANUSCRIPT CITATION:

McCutcheon, S.R., Swartz, A.M., Brown, M.C. *et al.* Transcriptional and epigenetic regulators of human CD8⁺ T cell function identified through orthogonal CRISPR screens. *Nat Genet* (2023). https://doi.org/10.1038/s4158 8-023-01554-0

License: This is an open access protocol distributed under the terms of the Creative Commons
Attribution License, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

TFome CRISPRko screens

Andrea R Daniel¹

¹Duke University

Andrea R Daniel: This protocol was adapted from work of Sean McCutcheon and colleagues in the Gersbach lab at Duke University.



Andrea R Daniel
Duke University

ABSTRACT

This protocols describes methods for transcription factor CRISPR knockout screens to reveal cofactors of BATF3 in T cells.

Created: Dec 08, 2023

Last Modified: Dec 08,

2023

PROTOCOL integer ID:

92033

Keywords: CRISPR, screen, knockout, BATF3, T cells

Funders

Acknowledgement:

NIH

gRNA library construction

- 1 The Brunello genome-wide KO^{83} library (four gRNAs per gene) was subset for 1,612 TFs⁴⁵ and *IL7R*.
- A total of 550 NT gRNAs were included in the library for a total of 7,000 gRNAs (available in Supplementary Table 6 of McCutcheon et al. Nature Genetics, 2023. https://doi.org/10.1038/s41588-023-01554-0).
- 3 This gRNA library was cloned into SpCas9 gRNA lentiviral plasmids with either mCherry or BATF3.

gRNA library cloning

- 4 Oligonucleotide pools containing variable gRNA sequences and constant regions for polymerase chain reaction (PCR) amplification were synthesized by Twist Bioscience.
- gRNA amplicons were gel extracted, PCR purified and input into 20 µl Gibson reactions (5:1 molar ratio of insert to backbone) with 200 ng of Esp3l digested and 1 × solid-phase reversible immobilization (SPRI)-selected (Beckman Coulter) plasmid backbone.
- **6** Gibson reactions were purified using ethanol precipitation and transformed into Lucigen's Endura ElectroCompetent Cells.

7 Transformed cells were cultured overnight and plasmids were isolated using Qiagen Midi Kits.

TFome CRISPRko screens

- A total of 20×10^6 CD8⁺ T cells from two donors were activated with CD3/CD28 dynabeads at a 1:1 ratio.
- At 24 h post-activation, CD8⁺ T cells were split evenly and transduced in parallel with TFome CRISPRko gRNA libraries with mCherry or BATF3.
- At 48 h post-activation, cells were electroporated with Cas9 protein. Briefly, the cells were collected, spun down at 90 g for 10 min, resuspended in 100 μl of Lonza P3 Primary Cell buffer with 3.2 μg Cas9 (Thermo) per 10⁶ cells, and electroporated with the pulse code EH115.
- After electroporation, warm medium was immediately added to each cuvette and cells were recovered at 37 °C for 20 min before being transferred into a six-well plate.
- On day 3 post transduction, cells were selected with 2 μg ml⁻¹ of puromycin for 3 days. On day 9 post transduction, cells were stained for CD8, IL7R and a viability dye.

12.1

	Antibod y Target	Fluorop hore/Se quence	Clone	Isotype	Dilution	Applicati on	Manufa cturer	Catalog #
	IL7RA	PE-Cy5	eBioRDR 5	Mouse / lgG1, kappa	1:100	Flow cytomet ry	Thermo	15- 1278-42
	CD8	bv-421	HIT8a	Mouse lgG1, κ	1:50	Flow cytomet ry	BD Bioscien ces	740078

