



Aug 26, 2022

PSI FluorCAM protocol for Rapid Light Curve with far red pre-illumination.

Andrei Herdean¹¹University of Technology Sydney

1 Works for me

 Sharedx.doi.org/10.17504/protocols.io.x54v9deypg3e/v1 Andrei Herdean
University of Technology Sydney

ABSTRACT

PSI FluorCAM protocol for Rapid Light Curve with far red pre-illumination.

DOI

dx.doi.org/10.17504/protocols.io.x54v9deypg3e/v1

PROTOCOL CITATION

Andrei Herdean 2022. PSI FluorCAM protocol for Rapid Light Curve with far red pre-illumination.. **protocols.io**
<https://protocols.io/view/psi-fluorcam-protocol-for-rapid-light-curve-with-f-cfuetnte>



LICENSE

————— This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](https://creativecommons.org/licenses/by/4.0/), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

CREATED

Aug 26, 2022

LAST MODIFIED

Aug 26, 2022

PROTOCOL INTEGER ID

69222

PSI FluorCAM protocol for Rapid Light Curve with far red pre-illumination.

- 1 Load the following script in to the FluorCam software

2 TS=50ms
 ;Rapid light curve in Actinic light 2 for TOMI-2
 ;version November 11, 2020
 ;with FilterWheel
 ;high-resolution CCD - TOMI-2
 ;estimated time: app. 372 s
 ;for successful measurement check Light curve settings manual
 ADD1=0
 Act2=0
 Act1=0
 ;
 include default.inc ;Includes standard options, do not remove it !
 include light.inc ;Includes standard options, do not remove it !
 include FW.inc ;Includes standard options, do not remove it !

 Shutter=2
 Sensitivity=30.6
 Super=60
 ADD2=10

 __LightA=24.122
 __LightB=-61.463

 ;3 = 10.903
 ;4 = 35.025
 ;5 = 59.147
 ;6 = 83.269
 ;7 = 107.391
 ;8 = 131.513
 ;9 = 155.635
 ;10 = 179.757
 ;15 = 300.367
 ;20 = 420.977
 ;25 = 541.587
 ;30 = 662.197
 ;35 = 782.807
 ;40 = 903.417

 ;These are the light intensities used for the Rapid light curve
 I1 = 3
 I2 = 3
 I3 = 3
 I4 = 3
 I5 = 3
 I6 = 3

I7 = 3
 I8 = 3
 I9 = 3
 I10 = 3
 I11 = 4
 I12 = 5
 I13 = 6
 I14 = 7
 I15 = 8
 I16 = 9
 I17 = 10
 I18 = 12
 I19 = 17
 I20 = 20
 I21 = 25
 I22 = 30
 I23 = 35
 I24 = 40
 I25 = 45

__LightIntensity=

;

;There is a pre-illumination phase with only far red light, followed by determination of Fo and Fm

<0s>=>SET_FILTER(CHL)

Preillumination=60s

<0s>=>add2(Preillumination)

start = Preillumination;

;

*** Fo Measurement

;

F0duration = 5s

F0measure = 1s

start + <0,F0measure..F0duration>=>mfmsub

start + <0s>=>checkPoint,"startFo"

start + =>checkPoint,"endFo"

;

;

*** Saturating Pulse & Fm Measurement

;

PulseDuration = 800ms;

Fmstart = start + F0duration+mfmsub_length

;

=>SI_Act2(I25)

```

=>SatPulse(PulseDuration)
=>act2(PulseDuration)
Fmstart+<200ms, 200ms + mfmsub_length .. PulseDuration-mfmsub_length=>mfmsub
;
=>checkPoint,"startFm"
=>checkPoint,"endFm"
=>checkPoint,"timeVisual"
;
;*** Actinic light Exposure *****
Pause = 2*TS
ALduration = 30s
ALmeasure = ALduration/6
;Protocol contains 5 minutes of illumination with low light in order to keep photosynthesis
active while samples are brogth to target temperatures by the thermocycler
ALstart1 = Fmstart + PulseDuration + Pause
ALstart2 = ALstart1 + ALduration + Pause
ALstart3 = ALstart2 + ALduration + Pause
ALstart4 = ALstart3 + ALduration + Pause
ALstart5 = ALstart4 + ALduration + Pause
ALstart6 = ALstart5 + ALduration + Pause
ALstart7 = ALstart6 + ALduration + Pause
ALstart8 = ALstart7 + ALduration + Pause
ALstart9 = ALstart8 + ALduration + Pause
ALstart10 = ALstart9 + ALduration + Pause
;this is where the Light Curve starts
ALstart11 = ALstart10 + ALduration + Pause
ALstart12 = ALstart11 + ALduration + Pause
ALstart13 = ALstart12 + ALduration + Pause
ALstart14 = ALstart13 + ALduration + Pause
ALstart15 = ALstart14 + ALduration + Pause
ALstart16 = ALstart15 + ALduration + Pause
ALstart17 = ALstart16 + ALduration + Pause
ALstart18 = ALstart17 + ALduration + Pause
ALstart19 = ALstart18 + ALduration + Pause
ALstart20 = ALstart19 + ALduration + Pause
ALstart21 = ALstart20 + ALduration + Pause
ALstart22 = ALstart21 + ALduration + Pause
ALstart23 = ALstart22 + ALduration + Pause
ALstart24 = ALstart23 + ALduration + Pause
ALstart25 = ALstart24 + ALduration + Pause
;
=>SI_Act2(I1)
=>SI_Act2(I2)
=>SI_Act2(I3)
=>SI_Act2(I4)
=>SI_Act2(I5)
=>SI_Act2(I6)

```

```

=>SI_Act2(I7)
=>SI_Act2(I8)
=>SI_Act2(I9)
=>SI_Act2(I10)
=>SI_Act2(I11)
=>SI_Act2(I12)
=>SI_Act2(I13)
=>SI_Act2(I14)
=>SI_Act2(I15)
=>SI_Act2(I16)
=>SI_Act2(I17)
=>SI_Act2(I18)
=>SI_Act2(I19)
=>SI_Act2(I20)
=>SI_Act2(I21)
=>SI_Act2(I22)
=>SI_Act2(I23)
=>SI_Act2(I24)
=>SI_Act2(I25)
;
ALstart =
ALstart => act2(ALduration)
;
;***** Kautsky Effect Measurement *****
;
;ALstart + TS => mfmsub
ALstart#=>mfmsub
;
;***** Ft' & Fm' definition *****
PulseLStart = ALduration - PulseDuration - mfmsub_length
;
FtL = ALstart + PulseLStart - PulseDuration
FtL#=>mfmsub
;
PulseL = ALstart + PulseLStart
PulseL - TS=> SI_Act2(I25)
PulseL=>SatPulse(PulseDuration)
PulseL#<200ms, 200ms + mfmsub_length .. PulseDuration-mfmsub_length>=>mfmsub
;
Last = ALstart25+ALduration+2*mfmsub_length
Last + <0s, 1s .. 5s>=>mfmsub
;
;
=>checkPoint,"startFm_Lss1"
=>checkPoint,"endFm_Lss1"
;
=>checkPoint,"startFm_Lss2"

```

```

=>checkPoint,"endFm_Lss2"
;
=>checkPoint,"startFm_Lss3"
=>checkPoint,"endFm_Lss3"
;
=>checkPoint,"startFm_Lss4"
=>checkPoint,"endFm_Lss4"
;
=>checkPoint,"startFm_Lss5"
=>checkPoint,"endFm_Lss5"
;
=>checkPoint,"startFm_Lss6"
=>checkPoint,"endFm_Lss6"
;
=>checkPoint,"startFm_Lss7"
=>checkPoint,"endFm_Lss7"
;
=>checkPoint,"startFm_Lss8"
=>checkPoint,"endFm_Lss8"
;
=>checkPoint,"startFm_Lss9"
=>checkPoint,"endFm_Lss9"
;
=>checkPoint,"startFm_Lss10"
=>checkPoint,"endFm_Lss10"
;
=>checkPoint,"startFm_Lss11"
=>checkPoint,"endFm_Lss11"
;
=>checkPoint,"startFm_Lss12"
=>checkPoint,"endFm_Lss12"
;
=>checkPoint,"startFm_Lss13"
=>checkPoint,"endFm_Lss13"
;
=>checkPoint,"startFm_Lss14"
=>checkPoint,"endFm_Lss14"
;
=>checkPoint,"startFm_Lss15"
=>checkPoint,"endFm_Lss15"
;
=>checkPoint,"startFm_Lss16"
=>checkPoint,"endFm_Lss16"
;
=>checkPoint,"startFm_Lss17"
=>checkPoint,"endFm_Lss17"
;

```

```

=>checkPoint,"startFm_Lss18"
=>checkPoint,"endFm_Lss18"
;
=>checkPoint,"startFm_Lss19"
=>checkPoint,"endFm_Lss19"
;
=>checkPoint,"startFm_Lss20"
=>checkPoint,"endFm_Lss20"
;
=>checkPoint,"startFm_Lss21"
=>checkPoint,"endFm_Lss21"
;
=>checkPoint,"startFm_Lss22"
=>checkPoint,"endFm_Lss22"
;
=>checkPoint,"startFm_Lss23"
=>checkPoint,"endFm_Lss23"
;
=>checkPoint,"startFm_Lss24"
=>checkPoint,"endFm_Lss24"
;
=>checkPoint,"startFm_Lss25"
=>checkPoint,"endFm_Lss25"
;
;***** Ft' definition *****
;
=>checkPoint,"startFt_Lss1"
=>checkPoint,"endFt_Lss1"
;
=>checkPoint,"startFt_Lss2"
=>checkPoint,"endFt_Lss2"
;
=>checkPoint,"startFt_Lss3"
=>checkPoint,"endFt_Lss3"
;
=>checkPoint,"startFt_Lss4"
=>checkPoint,"endFt_Lss4"
;
=>checkPoint,"startFt_Lss5"
=>checkPoint,"endFt_Lss5"
;
=>checkPoint,"startFt_Lss6"
=>checkPoint,"endFt_Lss6"
;
=>checkPoint,"startFt_Lss7"
=>checkPoint,"endFt_Lss7"
;

```

```

=>checkPoint,"startFt_Lss8"
=>checkPoint,"endFt_Lss8"
;
=>checkPoint,"startFt_Lss9"
=>checkPoint,"endFt_Lss9"
;
=>checkPoint,"startFt_Lss10"
=>checkPoint,"endFt_Lss10"
;
=>checkPoint,"startFt_Lss11"
=>checkPoint,"endFt_Lss11"
;
=>checkPoint,"startFt_Lss12"
=>checkPoint,"endFt_Lss12"
;
=>checkPoint,"startFt_Lss13"
=>checkPoint,"endFt_Lss13"
;
=>checkPoint,"startFt_Lss14"
=>checkPoint,"endFt_Lss14"
;
=>checkPoint,"startFt_Lss15"
=>checkPoint,"endFt_Lss15"
;
=>checkPoint,"startFt_Lss16"
=>checkPoint,"endFt_Lss16"
;
=>checkPoint,"startFt_Lss17"
=>checkPoint,"endFt_Lss17"
;
=>checkPoint,"startFt_Lss18"
=>checkPoint,"endFt_Lss18"
;
=>checkPoint,"startFt_Lss19"
=>checkPoint,"endFt_Lss19"
;
=>checkPoint,"startFt_Lss20"
=>checkPoint,"endFt_Lss20"
;
=>checkPoint,"startFt_Lss21"
=>checkPoint,"endFt_Lss21"
;
=>checkPoint,"startFt_Lss22"
=>checkPoint,"endFt_Lss22"
;
=>checkPoint,"startFt_Lss23"
=>checkPoint,"endFt_Lss23"

```



```
;
=>checkPoint,"startFt_Lss24"
=>checkPoint,"endFt_Lss24"
;
=>checkPoint,"startFt_Lss25"
=>checkPoint,"endFt_Lss25"

;
;
;END *****
```